

Technical Data Sheet

Doc No – Ab/TDS/LH-80

Rev No : 03 Date 11/2019

TECHNICAL DATA SHEET

Abdos Labtech Pvt. Ltd.

Rave Pipette Controller, PP

Valid for Catalog No. E11113

1. Description/Specification

1.1 Description:

Rave Pipette controller ideal for Repetitive Pipetting performance

Maximizing Pipetting comfort

1.2 Design:

Robust, Light weight & simple to use operation.

1.3 Performance Specifications:

1.3.1 Can use volumetric or serological pipettes ranging from 0.1 ml, to 100 ml.

1.3.2 It contains integral & replacement 0.3 μ m Filter.

1.3.3 Autoclavable silicon pipette holder

1.4 Material : PP (polypropylene), 21 CFR

1.5 Color : Grey

1.6 Basic Fetures

1.6.1 Made up of highly Virgin thermoplastic polypropylene

1.6.2 Autoclavability : Yes (Only Silicon Pipette Holder – Don't autoclave whole instrument)

1.6.3 CE certified



Shiv Ganga Industrial Estate,
Plot no.1 Lakeshwari, Bhagwanpur,
Roorkee – 247667 Uttrakhand, India
An ISO 9001 & 13485 Certified
A Division of Abdos


Regd. Office:
13A, Government Place (East)
Kolkata – 700069, India
CIN : U85100WB2008PTC123903
GST : 05AAGCA8605A1Z7


P : +91 (0) 1332 235277
E : labtech@abdosindia.com
W: www.abdoslifesciences.com

Technical Data Sheet

1.7 Dimension:

Code	Dimension
E11113	: Length- 210 mm, ϕ -60mm
E11119	: (Spare) 3.0 μ m Filter


Quality Assurance Incharge
Abdos Labtech Pvt. Ltd.




ABDOS

Shiv Ganga Industrial Estate,
Plot no.1 Lakeshwari, Bhagwanpur,
Roorkee – 247667 Uttrakhand, India
An ISO 9001 & 13485 Certified
A Division of Abdos

Regd. Office:
13A, Government Place (East)
Kolkata – 700069, India
CIN : U85100WB2008PTC123903
GST : 05AAGCA8605A1Z7

P : +91 (0) 1332 235277
E : labtech@abdosindia.com
W: www.abdoslifesciences.com