

Technical Data Sheet

Doc No – Ab/TDS/GLP-181

Rev No : 04 Date 11/2019

TECHNICAL DATA SHEET

Abdos Labtech Pvt. Ltd.

Triangular Magnetic Stirrer Bar, PTFE/Alnico v Magnet

Valid for Catalog No.

P80130, P80131, P80132, P80182

1. Description/Specification

1.1 Description:

Triangular Magnetic Stirrer Bar

1.2 Design :

Triangular

1.3 Performance Specification :

1.3.1 Provides low speed mixing of high volume of media and buffer solution.

1.3.2 Reduces marring of plastic containers.

1.3.3 Multifaceted surface adds turbulence

1.4 Material :

1.4.1 Pure Virgin Grade PTFE (Teflone) – Durable and inert, ANICO V Magnet

1.5 Color: White

1.6 Quality Control :

Testing of various characteristic in accordance with the valid Specification

1.7 Basic Features:

1.7.1 Made up of Highly Pure Virgin PTFE Material.

1.7.2 Micro Magnetic Stirrer bars used for stirring of liquids on magnetic stirrers & hot plate stirrers

1.7.3 Resistance to virtual all chemical attacks.



Shiv Ganga Industrial Estate,
Plot no.1 Lakeshwari, Bhagwanpur,
Roorkee – 247667 Uttrakhand, India
An ISO 9001 & 13485 Certified
A Division of Abdos


Regd. Office:
13A, Government Place (East)
Kolkata – 700069, India
CIN : U85100WB2008PTC123903
GST : 05AAGCA8605A1Z7


P : +91 (0) 1332 235277
E : labtech@abdosindia.com
W: www.abdoslifesciences.com

Technical Data Sheet

1.8. With Stand Temperature up to 260°C.

Code	Type	H x L (mm)	No. pre Case	Basic Features
P80130	Triangular Double ended Magnetic Stirrer Bar	14x25	2	-Ideal for mixing difficult solutions e.g Dissolving solvents or mixing sediments -Spins well giving the best possible mixing action -Resistance with chemicals & with stand upto 260°C
P80131		14x40	2	
P80132		14x55	2	
P80182		14x80	2	


Quality Assurance Incharge
Abdos Labtech Pvt. Ltd.




ABDOS

Shiv Ganga Industrial Estate,
Plot no.1 Lakeshwari, Bhagwanpur,
Roorkee – 247667 Uttrakhand, India
An ISO 9001 & 13485 Certified
A Division of Abdos

Regd. Office:
13A, Government Place (East)
Kolkata – 700069, India
CIN : U85100WB2008PTC123903
GST : 05AAGCA8605A1Z7

P : +91 (0) 1332 235277
E : labtech@abdosindia.com
W: www.abdoslifesciences.com