

Technical Data Sheet

Doc No–Ab/TDS/UTL-225

Rev No: 04 Date 11/2019

TECHNICAL DATA SHEET

Abdos Labtech Pvt. Ltd.

Abdos Brand – Derma Gaurd Blue Nitrile Gloves – 9.5 inches length

Valid for Catalog No. U20341, U20342, U20343, U20447, U20448

1. Description/Specification

1.1 Description: Abdos Brand – Derma Gaurd Blue Nitrile Gloves – 9.5 inches length

1.1.1 These gloves are ideal product to protect you exposure to bacteria, viruses and chemicals. It also enhances your ability to offer maximum security against cross contamination.

1.1.2 Its unique blue color provides quick, visual differentiation from natural rubber latex and other synthetic gloves.

1.1.3 Ambidextrous, textured fingertips and beaded cuff enhances wet & dry grip.

1.1.4 Custom design enhances comfort & fit.

1.1.5 These gloves are powder free, minimizing the potential for powder related complications.

1.1.6 Conforms to ASTM D6319 and EN 455 Standards.

1.1.7 Derma Gaurd Blue Nitrile Gloves Powder free- sample meets the overall migration requirement under Commission Regulation (EU) No 10/2011- Plastic materials and articles shall not transfer their constituents to foodstuff in quantities exceeding 10mg of total constituents released per dm² of food contact surface (mg/dm²) (overall migration limit)

1.2 Design: Hand gloves design

1.3 Performance Specifications:

Code	Size	No. Per Case
U20341	Small	100
U20342	Medium	100
U20343	Large	100
U20447	X Small	100
U20448	X Large	100

ABDOS

Shiv Ganga Industrial Estate,
Plot no.1 Lakeshwari, Bhagwanpur,
Roorkee – 247667 Uttrakhand, India
An ISO 9001 & 13485 Certified
A Division of Abdos

Regd. Office:
13A, Government Place (East)
Kolkata – 700069, India
CIN : U85100WB2008PTC123903
GST : 05AAGCA8605A1Z7

P : +91 (0) 1332 235277
E : labtech@abdosindia.com
W: www.abdoslifesciences.com

Technical Data Sheet

1.4 Material: Nitrile gloves

1.5 Color: Blue

1.6 Quality control:

Testing of various characteristics in accordance with the valid Specification


1.7 Other information

1.7.1 These gloves are used for working on lab, machine and any other work.

1.7.2 These gloves are ideal product to protect you any of microbial contamination.

1.7.3 These Flexible gloves are fit in hand & easy to use on working.

1.7.4 Powder free Chlorinated gloves



ABDOS

Shiv Ganga Industrial Estate,
Plot no.1 Lakeshwari, Bhagwanpur,
Roorkee – 247667 Uttrakhand, India
An ISO 9001 & 13485 Certified
A Division of Abdos

Regd. Office:
13A, Government Place (East)
Kolkata – 700069, India
CIN : U85100WB2008PTC123903
GST : 05AAGCA8605A1Z7

P : +91 (0) 1332 235277
E : labtech@abdosindia.com
W: www.abdoslifesciences.com