

# Technical Data Sheet

Doc No–Ab/TDS/GLP-173

Rev No: 01 Date 11/2019

## TECHNICAL DATA SHEET

**Abdos Labtech Pvt. Ltd.**

**Animal Cage PP/Stainless Steel**

**Valid for Catalog No. P20501, P20502, P20503, P20504**

### **1. Description/Specification**

#### **1.1 Description:** Animal Cage PP/Stainless Steel

1.1.1 These provide excellent shelter for animal species in the laboratory.

1.1.2 Metallic top provides good ventilation & provides for feeding the animals species.

#### **1.2 Design:**

Animal Cage PP/Stainless Steel

#### **1.3 Performance Specifications:**

<b>Code</b>	<b>Places</b>	<b>L x W x H (mm)</b>	<b>No. Per Case</b>
<b>P20501</b>	Single Grill	290 x 220 x 140	1
<b>P20502</b>	Single Grill	430x270x150	1
<b>P20503</b>	Twin Grill	290x220x140	1
<b>P20504</b>	Twin Grill	430x270x150	1

#### **1.4 Material:**

1.4.1 Polypropylene/Stainless steel

1.4.2 Comply with requirement of Directive 2002/72/EC and requirement of FDA Regulation 21 CFR olefin polymers.

1.4.3 Complies with the requirement of Indian Standard IS: 10910 on “ Specification for Polypropylene and its copolymer for its safe use in contact with food stuffs, Pharmaceutical and Drinking water”



**Shiv Ganga Industrial Estate,**  
Plot no.1 Lakeshwari, Bhagwanpur,  
Roorkee – 247667 Uttarakhand, India  
**An ISO 9001 & 13485 Certified**  
A Division of Abdos

**Regd. Office:**  
13A, Government Place (East)  
Kolkata – 700069, India  
**CIN :** U85100WB2008PTC123903  
**GST :** 05AAGCA8605A1Z7

**P :** +91 (0) 1332 235277  
**E :** labtech@abdosindia.com  
**W:** www.abdoslifesciences.com

# Technical Data Sheet

1.4.4 Free from heavy metals according to ELV directive 2000/53/EC and its following amendments.

**1.5 Color:** White

**1.6 Quality control:**

Testing of various characteristics in accordance with the valid Specification

**1.7 Other information**

1.7.1 Strong and lightweight material.


1.7.2 Good impact strength

1.7.3 Better chemical resistance.

**1.8. Basic Features**

**1.8.1 Autoclavability:** autoclavable.

**1.8.2 Chemical resistance:** see website.

  
Quality Assurance Incharge  
Abdos Labtech Pvt. Ltd.



  
**ABDOS**

**Shiv Ganga Industrial Estate,**  
Plot no.1 Lakeshwari, Bhagwanpur,  
Roorkee – 247667 Uttrakhand, India  
**An ISO 9001 & 13485 Certified**  
A Division of Abdos

**Regd. Office:**  
13A, Government Place (East)  
Kolkata – 700069, India  
**CIN :** U85100WB2008PTC123903  
**GST :** 05AAGCA8605A1Z7

**P :** +91 (0) 1332 235277  
**E :** labtech@abdosindia.com  
**W:** www.abdoslifesciences.com