

# Technical Data Sheet

Doc No– Ab/TDS/LH-81

Rev No: 02 Date 11/2019

## TECHNICAL DATA SHEET

**Abdos Labtech Pvt. Ltd.**

**Fast Release Pipette Pump**

**Valid for Catalog No. E11219, E11220, E11221**

### **1. Description/Specification**

#### **1.1 Description:** Fast Release Pipette Pump

1.1.1 Fast Release Pipette Pump ensures quick & accurate operation.

1.1.2 These light weight pipette pump have a fast release trigger that dispenses the entire pipette content.

1.1.3 The thumb wheel easily rotates for precision aspirating or dispensing, down to a volume of one drop at a time.

1.1.4 The thumb wheel can be used to blow out the contents of the pipette.

1.1.5 Resistant to acids, alkalis & easily disassembly for cleaning.

#### **1.2 Design:**

Fast Release Pipette Pump

#### **1.3 Performance Specifications:**

<b>Code</b>	<b>Capacity (ml)</b>	<b>No. per Case</b>
<b>E11219</b>	2	4
<b>E11220</b>	10	4
<b>E11221</b>	25	4



**Shiv Ganga Industrial Estate,**  
Plot no.1 Lakeshwari, Bhagwanpur,  
Roorkee – 247667 Uttrakhand, India  
**An ISO 9001 & 13485 Certified**  
A Division of Abdos

**Regd. Office:**  
13A, Government Place (East)  
Kolkata – 700069, India  
**CIN :** U85100WB2008PTC123903  
**GST :** 05AAGCA8605A1Z7

**P :** +91 (0) 1332 235277  
**E :** labtech@abdosindia.com  
**W:** www.abdoslifesciences.com

# Technical Data Sheet

## 1.4 Quality control:


Testing of various characteristic in accordance with the valid Specification

## 1.5. Basic Features:

**1.5.1 Autoclavability:** NO, Non autoclavable

**1.5.2 Chemical resistance:** see website

**1.5.3 CE certified**



**ABDOS**

**Shiv Ganga Industrial Estate,**  
Plot no.1 Lakeshwari, Bhagwanpur,  
Roorkee – 247667 Uttrakhand, India  
**An ISO 9001 & 13485 Certified**  
A Division of Abdos

**Regd. Office:**  
13A, Government Place (East)  
Kolkata – 700069, India  
**CIN :** U85100WB2008PTC123903  
**GST :** 05AAGCA8605A1Z7

**P :** +91 (0) 1332 235277  
**E :** labtech@abdosindia.com  
**W:** www.abdoslifesciences.com