

Technical Data Sheet

Doc No– Ab/TDS/LH-82

Rev No: 01 Date 11/2019

TECHNICAL DATA SHEET

Abdos Labtech Pvt. Ltd.

Handypette Pipette Aid.

Valid for Catalog No. E11218, E11216, E11222

1. Description/Specification

1.1 Description: Handypette Pipette Aid.

1.1.1 Handypette Pipette Aid ensures quick & accurate operation.

1.1.2 Colour Coded by volume with green for 10 ml, red for 25 ml & blue for 2ml.

1.1.3 Users can easily insert glass & plastic pipette into the collar with stabilizing fingers for secure hold.

1.1.4 Resistant to acids, alkalis & easily disassembly for cleaning.

1.1.5 Rotating the knurled thumbwheel draws liquid to desired volumes. To release liquid apply light pressure to the fast release lever. For graduation dispensing rotate the thumb wheel.

1.2 Design:

Handypette Pipette

1.3 Performance Specifications:

Code	Capacity (ml)	No. per Case
E11218	10	4
E11216	25	4
E11222	2	4



Shiv Ganga Industrial Estate,
Plot no.1 Lakeshwari, Bhagwanpur,
Roorkee – 247667 Uttrakhand, India
An ISO 9001 & 13485 Certified
A Division of Abdos

Regd. Office:
13A, Government Place (East)
Kolkata – 700069, India
CIN : U85100WB2008PTC123903
GST : 05AAGCA8605A1Z7

P : +91 (0) 1332 235277
E : labtech@abdosindia.com
W: www.abdoslifesciences.com

Technical Data Sheet

1.4 Quality control:


Testing of various characteristic in accordance with the valid Specification

1.5. Basic Features:

1.5.1 Autoclavability: NO, Non autoclavable

1.5.2 Chemical resistance: see website

1.5.3 CE certified



ABDOS

Shiv Ganga Industrial Estate,
Plot no.1 Lakeshwari, Bhagwanpur,
Roorkee – 247667 Uttrakhand, India
An ISO 9001 & 13485 Certified
A Division of Abdos

Regd. Office:
13A, Government Place (East)
Kolkata – 700069, India
CIN : U85100WB2008PTC123903
GST : 05AAGCA8605A1Z7

P : +91 (0) 1332 235277
E : labtech@abdosindia.com
W: www.abdoslifesciences.com