

Technical Data Sheet

Doc No– Ab/TDS/LH-55

Rev No: 02 Date 11/2019

TECHNICAL DATA SHEET

Abdos Labtech Pvt. Ltd.

Pipette Rack Stand, Acrylic

Valid for Catalog No. P70201

1. Description/Specification

1.1 Description: Pipettes Rack Stand, Acrylic

1.1.1 Keep your pipettes within easy inside this convenient rack.

1.1.2 Holds five standard size pipettes.

1.1.3 Racks are mounted on their own stand.

1.1.4 Keep working area or desk free from clutter

1.2 Design:

Pipette Rack Stand

1.3 Performance Specifications:

Code	Places	L x W x H (mm)	No. per Case
P70201	5	190x135x265	1

1.4 Material:

Acrylic, Poly (methyl methacrylate) PMMA, Confirms 21 CFR

1.5 Color:

Transparent



Shiv Ganga Industrial Estate,
Plot no.1 Lakeshwari, Bhagwanpur,
Roorkee – 247667 Uttrakhand, India
An ISO 9001 & 13485 Certified
A Division of Abdos

Regd. Office:
13A, Government Place (East)
Kolkata – 700069, India
CIN : U85100WB2008PTC123903
GST : 05AAGCA8605A1Z7

P : +91 (0) 1332 235277
E : labtech@abdosindia.com
W: www.abdoslifesciences.com

Technical Data Sheet

1.6 Quality control:

Testing of various characteristic in accordance with the valid Specification

1.7. Basic Features

1.7.1 Made up of Acrylic, Poly (methyl methacrylate) PMMA


1.7.2 Strong & Lightweight material.

1.7.3 Good impact strength, higher than both glass and polystyrene.

1.7.4 Fitment with ABDOS Pipettes only

1.7.5 **Autoclavability:** NO, Non autoclavable

1.7.6 **Chemical resistance:** see website.



Shiv Ganga Industrial Estate,
Plot no.1 Lakeshwari, Bhagwanpur,
Roorkee – 247667 Uttrakhand, India
An ISO 9001 & 13485 Certified
A Division of Abdos

Regd. Office:
13A, Government Place (East)
Kolkata – 700069, India
CIN : U85100WB2008PTC123903
GST : 05AAGCA8605A1Z7

P : +91 (0) 1332 235277
E : labtech@abdosindia.com
W: www.abdoslifesciences.com