

# Technical Data Sheet

Doc No– Ab/TDS/LH-56

Rev No: 02 Date 11/2019

## TECHNICAL DATA SHEET

**Abdos Labtech Pvt. Ltd.**

**Pipette Rack Stand (Carosol Type), PP**

**Valid for Catalog No. P70203**

### **1. Description/Specification**

**1.1 Description:** Pipettes Rack Stand, PP

1.1.1 Keep your pipettes within easy inside this convenient rack.

1.1.2 Holds six standard size pipettes.

1.1.3 Racks are mounted on their own stand.

1.1.4 Keep working area or desk free from clutter

### **1.2 Design:**

Pipette Rack Stand (Caeosol Type)

### **1.3 Performance Specifications:**

<b>Code</b>	<b>Places</b>	<b>No. per Case</b>
<b>P70203</b>	6	1

### **1.4 Material:**

Polypropylene (PP), Confirms 21 CFR

### **1.5 Color:**

Transparent



**Shiv Ganga Industrial Estate,**  
Plot no.1 Lakeshwari, Bhagwanpur,  
Roorkee – 247667 Uttrakhand, India  
**An ISO 9001 & 13485 Certified**  
A Division of Abdos

**Regd. Office:**  
13A, Government Place (East)  
Kolkata – 700069, India  
**CIN :** U85100WB2008PTC123903  
**GST :** 05AAGCA8605A1Z7

**P :** +91 (0) 1332 235277  
**E :** labtech@abdosindia.com  
**W:** www.abdoslifesciences.com

# Technical Data Sheet

## 1.6 Quality control:

Testing of various characteristic in accordance with the valid Specification

## 1.7. Basic Features

1.7.1 Made up of Polypropylene (PP)

1.7.2 Strong & Lightweight material.

1.7.3 Good impact strength, higher than both glass and polystyrene.


1.7.4 Fitment with ABDOS Pipettes only


1.7.5 **Autoclavability:** Yes, Autoclavable at 15 psi and 121<sup>0</sup>C for 20 minutes

1.7.6 **Chemical resistance:** see website.

## Autoclaving Guideline for Plastics

All plastic consumables are autoclavable at 121<sup>0</sup>C for 20 minutes, 15 psi/1atm. Do not use a dry cycle as this may cause deformation or weakening of the plastic. Allow temperature within the autoclave to return to at least 50<sup>0</sup>C before removing product.

  
Quality Assurance Incharge  
Abdos Labtech Pvt. Ltd.



  
**ABDOS**

**Shiv Ganga Industrial Estate,**  
Plot no.1 Lakeshwari, Bhagwanpur,  
Roorkee – 247667 Uttrakhand, India  
**An ISO 9001 & 13485 Certified**  
A Division of Abdos

**Regd. Office:**  
13A, Government Place (East)  
Kolkata – 700069, India  
**CIN :** U85100WB2008PTC123903  
**GST :** 05AAGCA8605A1Z7

**P :** +91 (0) 1332 235277  
**E :** labtech@abdosindia.com  
**W:** www.abdoslifesciences.com