

Technical Data Sheet

Doc No– Ab/TDS/LH-86

Rev No: 01 Date 11/2019

TECHNICAL DATA SHEET

Abdos Labtech Pvt. Ltd.

Pipette Storage Rack with Magnet, ABS

Valid for Catalog No. P70202

1. Description

1.1.1. Pipette Storage Rack with Magnet, ABS

1.1.2 These racks help to protect & organize glass and plastic pipettes.

1.1.3 Angled four shelved compartments keep 0.1 to 50 ml pipettes or cans within easy reach while preventing damage and clutter.

1.1.4 These have cut-out viewing window and come with two removable attachments that hold pipette bulbs, fillers and pumps.

1.1.5 These unique innovative magnetized racks have strong magnets on both sides and mounts securely to most metal surfaces.

1.2 Design: Pipette Storage Rack with Magnet, ABS

1.3 Performance Specifications:

Code	Description	L x W x H (mm)	No. per Case
P70202	Storage 1-50ml pipette or cans	291x88x401	1

1.4 Quality control:

Testing of various characteristic in accordance with the valid Specification

1.5 Material:

1.5.1 Acrylonitrile Butadiene Styrene (ABS), confirm US FDA 21CFR



Shiv Ganga Industrial Estate,
Plot no.1 Lakeshwari, Bhagwanpur,
Roorkee – 247667 Uttrakhand, India
An ISO 9001 & 13485 Certified
A Division of Abdos

Regd. Office:
13A, Government Place (East)
Kolkata – 700069, India
CIN : U85100WB2008PTC123903
GST : 05AAGCA8605A1Z7

P : +91 (0) 1332 235277
E : labtech@abdosindia.com
W: www.abdoslifesciences.com

Technical Data Sheet


1.5.2 Materials has successful passed the biological tests according to ISO 10993 external communication devices for indirect blood contact for Prolonged period.

1.5.3 Free from heavy metals according to ELV Direction 2000/53/EC, EU CONEG & it's following amendments. It is safe use in contact with food stuff, Pharmaceuticals & drinking water.

1.6. Basic Features

1.6.1 Autoclavability: NO, Non autoclavable

1.6.2 Chemical resistance: see website



ABDOS

Shiv Ganga Industrial Estate,
Plot no.1 Lakeshwari, Bhagwanpur,
Roorkee – 247667 Uttrakhand, India
An ISO 9001 & 13485 Certified
A Division of Abdos

Regd. Office:
13A, Government Place (East)
Kolkata – 700069, India
CIN : U85100WB2008PTC123903
GST : 05AAGCA8605A1Z7

P : +91 (0) 1332 235277
E : labtech@abdosindia.com
W: www.abdoslifesciences.com