

Technical Data Sheet

Doc No – Ab/TDS/LH-50

Rev No: 01 Date 11/2019

TECHNICAL DATA SHEET

Abdos Labtech Pvt. Ltd.

Premium Pipettes, PP

Valid for Catalog No. E11100, E11101, E11102, E11106, E11103, E11104, E11112, E11105, E11131, E11132

1. Description/Specification

1.1 Description: Premium Pipettes, PP

1.1.1 Abdos premium pipettes are semi autoclavable.

1.1.2 Abdos premium pipettes are lightweight, ergonomic, with low force design.

1.1.3 Volume settings can be clearly read through digital display.

1.1.4 Abdos premium pipettes are easy to calibrate & maintain with the tool supplied.

1.1.5 Abdos premium pipettes are designed to avoid repetitive strain injuries.

1.1.6 Abdos premium pipettes are supplied with operator manual that provides complete information on calibration & maintenance of the pipettes.

1.1.7 Abdos premium pipettes are suitable for transfer liquid volumes.

1.1.8 Matching universal style pipette tips, calibrated in accordance with ISO 8655/DIN 12650. The quality control involves gravimetric testing of each pipette with distilled water at 22°C.

1.2 Design:

Light weight & simple to use operation.

1.3 Performance Specifications:

Push Button, Cap, Push Button Rod, Tip ejector, Grip cover, Handle, Display cover, Tip ejector collar, Tip cone, Volume setting, Sealing & ejecting tips

1.4 Material: PP (polypropylene)



Shiv Ganga Industrial Estate,
Plot no.1 Lakeshwari, Bhagwanpur,
Roorkee – 247667 Uttrakhand, India
An ISO 9001 & 13485 Certified
A Division of Abdos

Regd. Office:
13A, Government Place (East)
Kolkata – 700069, India
CIN : U85100WB2008PTC123903
GST : 05AAGCA8605A1Z7

P : +91 (0) 1332 235277
E : labtech@abdosindia.com
W: www.abdoslifesciences.com

Technical Data Sheet

1.5 Color: White & Black

1.6 Basic Features

1.6.1 Made up of highly Virgin thermoplastic polypropylene

1.6.2 Abdos premium pipettes are semi autoclavable.


1.6.3 Autoclavability: Yes


1.6.4 CE certified

1.6.5 Chemical resistance: see website

1.7. Pipettes

Code	Capacity	No. Per Case
E11100	0.1-2.5 µl	1
E11101	0.5-10 µl	1
E11102	2-20 µl	1
E11106	5-50 µl	1
E11103	10-100 µl	1
E11104	20-200 µl	1
E11112	50-200 µl	1
E11105	100-1000 µl	1
E11131	1-5ml	1
E11132	2-10ml	1


Quality Assurance Incharge
Abdos Labtech Pvt. Ltd.



ABDOS

Shiv Ganga Industrial Estate,
Plot no.1 Lakeshwari, Bhagwanpur,
Roorkee – 247667 Uttrakhand, India
An ISO 9001 & 13485 Certified
A Division of Abdos

Regd. Office:
13A, Government Place (East)
Kolkata – 700069, India
CIN : U85100WB2008PTC123903
GST : 05AAGCA8605A1Z7

P : +91 (0) 1332 235277
E : labtech@abdosindia.com
W: www.abdoslifesciences.com