

Technical Data Sheet

Doc No – Ab/TDS/LH-79

Rev No: 03 Date 11/2019

TECHNICAL DATA SHEET

Abdos Labtech Pvt. Ltd.

Rave Plus Pipette Filler, PP

Valid for Catalog No. E11114, E11120

1. Description/Specification

1.1 Description: Rave Plus Pipette Filler, PP

1.1.1 Rave Plus Pipette Filler is powerful & can fill a 50 ml pipette in less than 10 seconds motor and pump operate quietly with very low vibration.

1.1.2 Users can use volumetric and serological Pipette ranging from 1 ml to 100 ml.

1.1.3 Rave Pipette Controller is supplied with integral and replaceable 0.45 µm filler.

1.2 Design:

Robust, Lightweight & simple to use operation.

1.3 Performance Specifications:

Battery	Lithium-Ion
Battery Service Life	More than 8 hours of intermittence use.
Charging Time	Approx. 2 hour
Pipette Types	Glass & Plastic pipette (0.1-100ml)
Aspirate Speeds	8
Dispense Speed	8 Gravity Dispense
Filter	0.45 µm Hydrophobic

1.4 Material: PP (polypropylene)



Shiv Ganga Industrial Estate,
Plot no.1 Lakeshwari, Bhagwanpur,
Roorkee – 247667 Uttrakhand, India
An ISO 9001 & 13485 Certified
A Division of Abdos

Regd. Office:
13A, Government Place (East)
Kolkata – 700069, India
CIN : U85100WB2008PTC123903
GST : 05AAGCA8605A1Z7

P : +91 (0) 1332 235277
E : labtech@abdosindia.com
W: www.abdoslifesciences.com

Technical Data Sheet

1.5 Color: Grey & White

1.6 Basic Features

1.6.1 Made up of highly Virgin thermoplastic polypropylene


1.6.2 Autoclavability : Yes (Only Silicon Pipette Holder – Don't autoclave whole instrument)


1.6.3 Chemical resistance: see website

1.6.4 CE certified

1.7. Spare Filters

Code	Product	No. Per Case
E11116	Spare 0.45 µm Filter	1
E11117	Spare Lithium-ion	1
E11118	Spare Adaptor	1


Quality Assurance Incharge
Abdos Labtech Pvt. Ltd.



ABDOS

Shiv Ganga Industrial Estate,
Plot no.1 Lakeshwari, Bhagwanpur,
Roorkee – 247667 Uttrakhand, India
An ISO 9001 & 13485 Certified
A Division of Abdos

Regd. Office:
13A, Government Place (East)
Kolkata – 700069, India
CIN : U85100WB2008PTC123903
GST : 05AAGCA8605A1Z7

P : +91 (0) 1332 235277
E : labtech@abdosindia.com
W: www.abdoslifesciences.com