

Technical Data Sheet

Doc No– Ab/TDS/LH-83

Rev No: 02 Date 11/2019

TECHNICAL DATA SHEET

Abdos Labtech Pvt. Ltd.

Rubber Pipette Pump

Valid for Catalog No. E11217

1. Description/Specification

1.1 Description: Rubber Pipette Pump

1.1.1 Rubber Pipette Pump ensures quick & accurate operation.

1.1.2 Pump opening is standard to fit almost all glass & plastic pipette.

1.1.3 One size can be used for liquid up to 100 ml.

1.1.4 The pump has a heavy wall and easy to grip body.

1.2 Design: Rubber Pipette Pump

1.3 Performance Specifications:

Code	Capacity (ml)	No. per Case
E11217	Up to 100ml	4

1.4 Quality control:


Testing of various characteristic in accordance with the valid Specification


1.5. Basic Features:

1.8.1 Autoclavability: NO, Non autoclavable

1.8.2 Chemical resistance: see website

1.8.3 Colour- Available in various colour


Quality Assurance Incharge
Abdos Labtech Pvt. Ltd.



ABDOS

Shiv Ganga Industrial Estate,
Plot no.1 Lakeshwari, Bhagwanpur,
Roorkee – 247667 Uttrakhand, India
An ISO 9001 & 13485 Certified
A Division of Abdos

Regd. Office:
13A, Government Place (East)
Kolkata – 700069, India
CIN : U85100WB2008PTC123903
GST : 05AAGCA8605A1Z7

P : +91 (0) 1332 235277
E : labtech@abdosindia.com
W: www.abdoslifesciences.com