

Technical Data Sheet

Doc No–Ab/TDS/GLP-150

Rev No: 02 Date 11/2019

TECHNICAL DATA SHEET

Abdos Labtech Pvt. Ltd.

Scoop, HDPE

Valid for Catalog No. P11901, P11902, P11903, P11904, P11905

1. Description/Specification

1.1 Description: Scoop, HDPE

1.1.1 These scoops used for withdrawal and movement of solids, granulates, gels and other sample.

1.1.2 Round bottom for better use.

1.1.3 It has comfortable, stable and robust handle with precise filling edge.

1.1.4 Natural in colour & superior chemical resistance.

1.1.5 Suitable for use with food.

1.2 Design: Round bottom Scoop

1.3 Performance Specifications:

Code	Capacity (ml)	Total Length (mm)	No. Per Case
P11901	100	180	12
P11902	250	240	12
P11903	500	300	6
P11904	1000	360	6
P11905	2500	400	3

1.4 Material:

1.3.1 High Density poly Ethylene (HDPE), Hydrophobic, confirm US FDA 21CFR

1.3.2 Complies with requirement accordance with article of regulation (EC) No. 1907/2006 and its amendments.



Shiv Ganga Industrial Estate,
Plot no.1 Lakeshwari, Bhagwanpur,
Roorkee – 247667 Uttrakhand, India
An ISO 9001 & 13485 Certified
A Division of Abdos

Regd. Office:
13A, Government Place (East)
Kolkata – 700069, India
CIN : U85100WB2008PTC123903
GST : 05AAGCA8605A1Z7

P : +91 (0) 1332 235277
E : labtech@abdosindia.com
W: www.abdoslifesciences.com

Technical Data Sheet

1.3.3 Complies with the requirement of European Pharmacopoeia, 7th edition (2011) & its supplements 7.5 (07/2012) monographs 3.2.2 Plastic Containers and closures for Pharmaceutical use.

1.3.4 The product has passed the USP testing including class VI tests & has been assigned the FDA Drug master file DMF 12123

1.3.5 Materials has successful passed the biological tests according to ISO 10993 external communication devices for indirect blood contact for Prolonged period.

1.3.6 Free from heavy metals according to ELV Direction 2000/53/EC, EU CONEG & it's following amendments. It is safe use in contact with food stuff, Pharmaceuticals & drinking water.

1.5 Color: Natural

1.6 Quality control:

Testing of various characteristics in accordance with the valid Specification


1.7 Other information


1.7.1 Impact resistant and tough.

1.7.2 Superior for its hardness, gloss, toughness, and electrical insulation properties.

1.7.3 **Chemical resistant:** see website.

CE mark	No
Autoclavability	No
Translucence	Yes
Resistant	To be effected by radiations
Sterilization	From – 40°C to + 80°C


Quality Assurance Incharge
Abdos Labtech Pvt. Ltd.



ABDOS

Shiv Ganga Industrial Estate,
Plot no.1 Lakeshwari, Bhagwanpur,
Roorkee – 247667 Uttrakhand, India
An ISO 9001 & 13485 Certified
A Division of Abdos

Regd. Office:
13A, Government Place (East)
Kolkata – 700069, India
CIN : U85100WB2008PTC123903
GST : 05AAGCA8605A1Z7

P : +91 (0) 1332 235277
E : labtech@abdosindia.com
W: www.abdoslifesciences.com