

Technical Data Sheet

Doc No – Ab/TDS/PCR-104

Rev No: 02 Date 11/2019

TECHNICAL DATA SHEET

Abdos Labtech Pvt. Ltd.

Sealing Films, Polyester

Valid for Cat. No. P10718

1. Description/Specification

1.1 Description: Sealing Films, Polyester.

1.1.1 Sealing films are optically clear and are suitable for use in real time PCR.

1.1.2 The film is adhesive coated and sticks to the plate. It can be used for PCR plates.

1.1.3 Sealing films are disposable and for single use only.

1.2 Design:

Sealing Films, Polyester designed to seal PCR plates.

1.3 Performance Specification:

Code	Type	No. per Case
P10718	Clear Seal Film	100 Sheets

1.4 Material: Sealing Films, Polyester, complies 21 CFR

1.5 Color:

Translucent (Natural)

1.6 Quality control:

Testing of various characteristics in accordance with the valid specification

1.7 Other information

1.7.1 DNase, RNase, DNA, Endotoxin and PCR inhibitor.

1.7.2 Recommended for Single use only



Shiv Ganga Industrial Estate,
Plot no.1 Lakeshwari, Bhagwanpur,
Roorkee – 247667 Uttarakhand, India
An ISO 9001 & 13485 Certified
A Division of Abdos

Regd. Office:
13A, Government Place (East)
Kolkata – 700069, India
CIN : U85100WB2008PTC123903
GST : 05AAGCA8605A1Z7

P : +91 (0) 1332 235277
E : labtech@abdosindia.com
W: www.abdoslifesciences.com


Technical Data Sheet

1.8 Basic Features

1.8.1 High dimensional stability, even at high temperature.

1.8.2 Very good thermal stability and resistance to UV, Ozone and weather

1.8.3 Chemical resistance: see website.



Shiv Ganga Industrial Estate,
Plot no.1 Lakeshwari, Bhagwanpur,
Roorkee – 247667 Uttrakhand, India
An ISO 9001 & 13485 Certified
A Division of Abdos

Regd. Office:
13A, Government Place (East)
Kolkata – 700069, India
CIN : U85100WB2008PTC123903
GST : 05AAGCA8605A1Z7

P : +91 (0) 1332 235277
E : labtech@abdosindia.com
W: www.abdoslifesciences.com