

# Technical Data Sheet

Doc No–Ab/TDS/GLP-161

Rev No: 02 Date 11/2019

## TECHNICAL DATA SHEET

**Abdos Labtech Pvt. Ltd.**

**Slide Drain Racks, HIPS**

**Valid for Catalog No. P90112**

### **1. Description/Specification**

#### **1.1 Description:** Slide Drain Racks, HIPS

1.1.1 The most convenient, organized and versatile way of storing microscope slides 75 x 25 mm.

1.1.2 This tray used as drain rack can hold up to 200 microscope slides in 100 individual numbered slots.

1.1.3 All slides are stored upright for easier insertion and removal.

1.1.4 To easily pick up the slide simply tilt the slides forward and backward with one finger.

#### **1.2 Design:**

Rectangular slide Drain Racks

#### **1.3 Performance Specifications:**

<b>Code</b>	<b>Places</b>	<b>L x W x H (mm)</b>	<b>No. Per Case</b>
<b>P90112</b>	100/200	75 x 231 x 25	10

#### **1.4 Material:**

1.4.1 HIPS (high Impact Polystyrene)

1.4.2 Comply with requirement of Directive 2002/72/EC and requirement of FDA Regulation 21 CFR olefin polymers.

1.4.3 Complies with the requirement of Indian Standard IS: 10910 on “ Specification for Polypropylene and its copolymer for its safe use in contact with food stuffs, Pharmaceutical and Drinking water”



**Shiv Ganga Industrial Estate,**  
Plot no.1 Lakeshwari, Bhagwanpur,  
Roorkee – 247667 Uttrakhand, India  
**An ISO 9001 & 13485 Certified**  
A Division of Abdos

**Regd. Office:**  
13A, Government Place (East)  
Kolkata – 700069, India  
**CIN :** U85100WB2008PTC123903  
**GST :** 05AAGCA8605A1Z7

**P :** +91 (0) 1332 235277  
**E :** labtech@abdosindia.com  
**W:** www.abdoslifesciences.com

# Technical Data Sheet

1.4.4 Free from heavy metals according to ELV directive 2000/53/EC and its following amendments.

**1.5 Color:** White

**1.6 Quality control:**

Testing of various characteristics in accordance with the valid Specification

**1.7 Other information**

1.7.1 Impact resistant and tough.


1.7.2 Temperature range from -80°C to 80°C

1.7.3 Superior for its hardness, gloss, toughness, and electrical insulation properties.

**1.8 Basic Features:**

**1.8.1 Autoclavability:** NO, Non autoclavable

**1.8.2 Chemical resistance:** see website



**ABDOS**

**Shiv Ganga Industrial Estate,**  
Plot no.1 Lakeshwari, Bhagwanpur,  
Roorkee – 247667 Uttrakhand, India  
**An ISO 9001 & 13485 Certified**  
A Division of Abdos

**Regd. Office:**  
13A, Government Place (East)  
Kolkata – 700069, India  
**CIN :** U85100WB2008PTC123903  
**GST :** 05AAGCA8605A1Z7

**P :** +91 (0) 1332 235277  
**E :** labtech@abdosindia.com  
**W:** www.abdoslifesciences.com