

Technical Data Sheet

Doc No – Ab/TDS/LH-75

Rev No: 02 Date 11/2019

TECHNICAL DATA SHEET

Abdos Labtech Pvt. Ltd.

Stepmate Stepper

Valid for Catalog No. E11171

1. Description/Specification

1.1 Description: Stepmate Stepper, 0.5-50ml capacity.

1.1.1 Stepmate Stepper is ideal for routine or research applications in diagnostics, molecular biology, environmental analysis and more.

1.1.2 Stepmate Stepper works on positive displacement pressure and offers highest precision dispensing of liquid media with high viscosity, high density and high vapor pressure.

1.1.3 It provides up to 48 dispensing steps and improves work efficiently.

1.1.4 It is lightweight and ergonomic.

1.1.5 Factory fitted calibration

1.2 Design:

Light weight & simple to use operation.

1.3 Material: PP (polypropylene)

1.4 Color: White & Black

1.5 Basic Features

1.5.1 Made up of highly Virgin thermoplastic polypropylene

1.5.2 It can be used with one hand and is maintenance free.

1.5.3 Chemical resistance: see website



Shiv Ganga Industrial Estate,
Plot no.1 Lakeshwari, Bhagwanpur,
Roorkee – 247667 Uttrakhand, India
An ISO 9001 & 13485 Certified
A Division of Abdos

Regd. Office:
13A, Government Place (East)
Kolkata – 700069, India
CIN : U85100WB2008PTC123903
GST : 05AAGCA8605A1Z7

P : +91 (0) 1332 235277
E : labtech@abdosindia.com
W: www.abdoslifesciences.com


Technical Data Sheet


1.6. Suitable Tips & Adapter for Stepmate Stepper

Code	Capacity (ml)	No. Per Case
E11089	0.05	100
E11090	0.5	100
E11091	1.25	100
E11092	2.5	100
E11093	5	100
E11094	12.5	100
E11095	25	25
E11096	50	25
E11080	0.5 sterile	100
E11081	1.25 sterile	100
E11082	2.5 sterile	100
E11083	5 sterile	100
E11084	12.5 sterile	100
E11085	25 sterile	25
E11086	50 sterile	25
E11908	Adapters for 50ml	1
E11909	Adapters for 25ml	1

1.6.1 The first and last dispensing liquid should be discarded for precision.

1.6.2 CE Certified


Quality Assurance Incharge
Abdos Labtech Pvt. Ltd.



ABDOS

Shiv Ganga Industrial Estate,
Plot no.1 Lakeshwari, Bhagwanpur,
Roorkee – 247667 Uttrakhand, India
An ISO 9001 & 13485 Certified
A Division of Abdos

Regd. Office:
13A, Government Place (East)
Kolkata – 700069, India
CIN : U85100WB2008PTC123903
GST : 05AAGCA8605A1Z7

P : +91 (0) 1332 235277
E : labtech@abdosindia.com
W: www.abdoslifesciences.com