

Technical Data Sheet

Doc No – Ab/TDS/LH-77

Rev No: 01 Date 11/2019

TECHNICAL DATA SHEET

Abdos Labtech Pvt. Ltd.

Supreme Plus Bottle Top Dispenser

Valid for Catalog No. E11711, E11712, E11713, E11714, E11715, E11716

1. Description/Specification

1.1 Description: Supreme Plus Bottle Top Dispenser

1.1.1 Supreme Plus is the most advanced next generation highly robust & reliable bottle top dispenser designed through years of research.

1.1.2 Supreme Plus has a smooth and soft plunger for easy working in difficult laboratory conditions.

1.1.3 Fully autoclavable.

1.1.4 Supreme Plus has a recirculation valve for redirecting the liquid into the bottle for bubble free dispensing without any loss of reagent.

1.2 Design:

Light weight & simple to use operation.

1.3 Performance Specifications:

Specifications	
Piston	PTFE & ETP
Cylinder	Borosilicate Glass
Volume Adjustment Knob	PP 180° Rotation
Valve Housing	PFA
Re-Circulation Valve Housing	PFA
Value Assembly	Borosilicate Glass Ball & Seat
Discharge assembly	PTFE
Delivery Tube	FEP
Inlet Tube	FEP
Calibration	Individually Calibrated & Certified
Accuracy & Reproducibility	According to ISO 8655 standards



Shiv Ganga Industrial Estate,
Plot no.1 Lakeshwari, Bhagwanpur,
Roorkee – 247667 Uttrakhand, India
An ISO 9001 & 13485 Certified
A Division of Abdos

Regd. Office:
13A, Government Place (East)
Kolkata – 700069, India
CIN : U85100WB2008PTC123903
GST : 05AAGCA8605A1Z7

P : +91 (0) 1332 235277
E : labtech@abdosindia.com
W: www.abdoslifesciences.com

Technical Data Sheet

Compatibility	Excellent with all reagents except HF
----------------------	---------------------------------------

1.4 Material: PP (polypropylene)

1.5 Color: White & Blue

1.6 Basic Features

1.6.1 Made up of highly Virgin thermoplastic polypropylene


1.6.2 Supreme Plus has a recirculation valve for redirecting the liquid into the bottle for bubble free dispensing without any loss of reagent.


1.6.3 Chemical resistance: see website

1.6.4 CE certified

1.7. Supreme Plus Bottle Top Dispenser

Code	Capacity(ml)	No. Per Case
E11711	0.25-2.5	1
E11712	0.5-5	1
E11713	1-10	1
E11714	2.5-30	1
E11715	5-60	1
E11716	10-100	1


Quality Assurance Incharge
Abdos Labtech Pvt. Ltd.



ABDOS

Shiv Ganga Industrial Estate,
Plot no.1 Lakeshwari, Bhagwanpur,
Roorkee – 247667 Uttrakhand, India
An ISO 9001 & 13485 Certified
A Division of Abdos

Regd. Office:
13A, Government Place (East)
Kolkata – 700069, India
CIN : U85100WB2008PTC123903
GST : 05AAGCA8605A1Z7

P : +91 (0) 1332 235277
E : labtech@abdosindia.com
W: www.abdoslifesciences.com