

Technical Data Sheet

Doc No Ab/TDS/CT-20

Rev No: 02 Date 11/2019

TECHNICAL DATA SHEET

Abdos Labtech Pvt. Ltd.

Swirl Micro Plate Centrifuge

Valid for Catalog No. E11550

1. Description Specifications:

1.1 Description: Swirl Micro Plate Centrifuge

1.1.1 Swirl Micro Plate Centrifuge is astonishingly small footprint (9x10 in.) half the size of conventional centrifuge.

1.1.2 Swirl micro plate centrifuge is easily reachable rotor chamber include two plate carriers that rest at a 75°C angle thus allows plates to be inserted confidently without seal in tapes or caps.

1.1.3 Upon closing the lid of the centrifuge, the rotor accelerates and the centrifugal force “swings” the plates into a horizontal position. Any droplets on the walls of the plate quickly concentrate into the well bottoms.

1.1.4 With a g-force of 600x g, most samples can be spun down in less than 20 sec.

1.1.5 Just open the lid and automatic electric brake brings the rotor to a quick & smooth stop in just 4 sec.

1.2 Design: Swirl Micro Plate Centrifuge

1.3 performance Specifications:

Specifications	E11550
Voltage	AC 220-230V
Frequency	50-60 Hz
Rotor Capacity	2x2 PCR plates, 2xmicroliter plates, 24x0.2ml PCR plates
Max Speed	21000 rpm /400x G
Deceleration	7 seconds



Shiv Ganga Industrial Estate,
Plot no.1 Lakeshwari, Bhagwanpur,
Roorkee – 247667 Uttarakhand, India
An ISO 9001 & 13485 Certified
A Division of Abdos

Regd. Office:
13A, Government Place (East)
Kolkata – 700069, India
CIN : U85100WB2008PTC123903

P : +91 (0) 1332 235277
E : labtech@abdosindia.com
W: www.abdoslifesciences.com

Technical Data Sheet

Dimension	230x260x197mm
Permissible Ambient Temperature	4- 40°C
Weight	4.1 kg
Permissible ambient Humidity	80%

1.4 Material: ABS (Acrylonitrile Butadiene Styrene)


1.5 Color: White

1.6 E11550 – Swirl Micro Plate Centrifuge, CE certified

1.7 Spare Accessories

1.7.1 **E11551** – Adapter, 96x 0.2ml or 12xPCR strips

1.7.2 **E11552** – Rotor with carrier for deepwell plates upto 35mm tall (max speed 1500rpm)


Quality Assurance Incharge
Abdos Labtech Pvt. Ltd.



ABDOS

Shiv Ganga Industrial Estate,
Plot no.1 Lakeshwari, Bhagwanpur,
Roorkee – 247667 Uttrakhand, India
An ISO 9001 & 13485 Certified
A Division of Abdos

Regd. Office:
13A, Government Place (East)
Kolkata – 700069, India
CIN : U85100WB2008PTC123903

P : +91 (0) 1332 235277
E : labtech@abdosindia.com
W: www.abdoslifesciences.com