

Technical Data Sheet

Doc No–Ab/TDS/GLP-169

Rev No: 02 Date 11/2019

TECHNICAL DATA SHEET

Abdos Labtech Pvt. Ltd.

Workstation Storage / Dispenser Bin, Acrylic

Valid for Catalog No. P70101, P70102

1. Description/Specification

1.1 Description: Workstation Storage / Dispenser Bin, Acrylic.

1.1.1 Workstation Storage / Dispenser Bin, is ideal for storing and dispensing small laboratory product such as Pipette tips, Pasture pipette, micro Centrifuge Tubes.

1.1.2 Removable top cover allows easy filling in to the box.

1.1.3 Hinged opening at the front to take out items from the bin.

1.1.4 Holes at the back of the bin help it to be placed or struck on the wall.

1.2 Design:

Workstation Storage / Dispenser Bin

1.3 Performance Specifications:

Code	Size	L x W x H (mm)	No. Per Case
P70101	Small	200 x 140 x 190	1
P70102	Large	200 x 140 x 275	1

1.4 Material:

1.4.1 Poly (methyl methacrylate) (PMMA) Acrylic

1.4.2 Comply with requirement of Directive 2002/72/EC and requirement of FDA Regulation 21 CFR olefin polymers.

1.4.3 Complies with the requirement of Indian Standard IS: 10910 on “ Specification for Polypropylene and its copolymer for its safe use in contact with food stuffs, Pharmaceutical and Drinking water”

1.4.4 Free from heavy metals according to ELV directive 2000/53/EC and its following amendments.



Shiv Ganga Industrial Estate,
Plot no.1 Lakeshwari, Bhagwanpur,
Roorkee – 247667 Uttrakhand, India
An ISO 9001 & 13485 Certified
A Division of Abdos

Regd. Office:
13A, Government Place (East)
Kolkata – 700069, India
CIN : U85100WB2008PTC123903
GST : 05AAGCA8605A1Z7

P : +91 (0) 1332 235277
E : labtech@abdosindia.com
W: www.abdoslifesciences.com

Technical Data Sheet

1.5 Color: White & Yellow

1.6 Quality control:

Testing of various characteristics in accordance with the valid Specification

1.7 Other information

1.7.1 Strong and lightweight material.

1.7.2 Good impact strength, higher than both glass & polystyrene.


1.7.3 Transmits up to 92% of visible light (3 mm thickness)

1.7.4 Better chemical stability.

1.8. Basic Features

1.8.1 Autoclavability: NO, Non autoclavable

1.8.2 Chemical resistance: see website.



ABDOS

Shiv Ganga Industrial Estate,
Plot no.1 Lakeshwari, Bhagwanpur,
Roorkee – 247667 Uttrakhand, India
An ISO 9001 & 13485 Certified
A Division of Abdos

Regd. Office:
13A, Government Place (East)
Kolkata – 700069, India
CIN : U85100WB2008PTC123903
GST : 05AAGCA8605A1Z7

P : +91 (0) 1332 235277
E : labtech@abdosindia.com
W: www.abdoslifesciences.com