



Last Drop™ Low Retention Pipette Tips

Creating For Life Sciences

www.abdoslifesciences.com

Last Drop™ Low Retention Pipette Tips Features




- Retention level is 10 times lesser than standard pipette tips.
- Highly homogeneous surface prevents non-specific binding.
- DNase, RNase, Human DNA, ATP, PCR Inhibitor and Pyrogen free.
- Tips are also free from Bacteriostatic & Fungistatic property that gives tips a unique feature to work in microbiological lab.
- Optimal surface for enhanced dispensing. Minimal sample loss of highly viscous liquids or samples containing detergents.
- The use of the latest injection moulding process produce high precision, thin wall tips with high clarity for better visibility of sample.
- Reduce sample retention for precision and to save on costly reagent.

ABDOS ENSURES LOWEST SAMPLE RETENTION WITH ITS LAST DROP PIPETTE TIPS

- Abdos last drop low retention pipette tips are manufactured from ultra pure USP class VI certified virgin polypropylene to ensure the highest quality & confirming to US FDA 21 CFR. Free from heavy metals, natural rubber & inhibitory plasticizers.
- The universal tip collar fits all popular pipettes.
- Diamond polished moulds for an extremely hydrophobic & smooth surface to ensure lowest retention in its class.
- The use of the latest injection moulding process produce high precision, thin walled tips with high clarity for better visibility of sample & highly homogeneous surface free of defect sites, preventing non-specific binding.
- Sterile tips has been validated according to ENISO 11137-2; SAL 10⁻⁶

Abdos Last Drop pipette tips clearly reduce sample loss. As seen in the graph, it clearly indicates that Abdos Tips retained least amount of residue and saves on costly reagents.



Code

Product Description

No.Per Case

PI0040	ABDOS Last Drop™ Low Retention Pipette Tips, Sterile Racked, Graduated, 0.1 – 10 µl	10 Racks x 96 pcs (960 pcs)
PI0041	ABDOS Last Drop™ Low Retention Pipette Tips, Sterile Racked, Extended Graduated, 0.1 – 10 µl	10 Racks x 96 pcs (960 pcs)
PI0042	ABDOS Last Drop™ Low Retention Pipette Tips, Sterile Racked, Bevelled Graduated, 1 – 200 µl	10 Racks x 96 pcs (960 pcs)
PI0043	ABDOS Last Drop™ Low Retention Pipette Tips, Sterile Racked, Ultra Point Tip Orifice, 1 – 200 µl	10 Racks x 96 pcs (960 pcs)
PI0044	ABDOS Last Drop™ Low Retention Pipette Tips, Sterile Racked, Graduated, 1 – 300 µl	10 Racks x 96 pcs (960 pcs)
PI0045	ABDOS Last Drop™ Low Retention Pipette Tips, Sterile Racked, 100 – 1000 µl	10 Racks x 96 pcs (960 pcs)
PI0046	ABDOS Last Drop™ Low Retention Pipette Tips, Sterile Racked, Graduated, 100 – 1000 µl	10 Racks x 96 pcs (960 pcs)
PI0047	ABDOS Last Drop™ Low Retention Pipette Tips, Sterile Racked, Graduated, 100 – 1250 µl	10 Racks x 96 pcs (960 pcs)
PI0020	ABDOS Last Drop Low Retention Pipette Tips, Racked, Graduated, 0.1 – 10 µl	10 Racks x 96 pcs (960 pcs)
PI0021	ABDOS Last Drop Low Retention Pipette Tips, Racked, Extended Graduated, 0.1 – 10 µl	10 Racks x 96 pcs (960 pcs)
PI0022	ABDOS Last Drop Low Retention Pipette Tips, Racked, Bevelled Graduated, 1 – 200 µl	10 Racks x 96 pcs (960 pcs)
PI0023	ABDOS Last Drop Low Retention Pipette Tips, Racked, Ultra Point Tip Orifice, 1 – 200 µl	10 Racks x 96 pcs (960 pcs)
PI0024	ABDOS Last Drop Low Retention Pipette Tips, Racked, Graduated, 1 – 300 µl	10 Racks x 96 pcs (960 pcs)
PI0025	ABDOS Last Drop Low Retention Pipette Tips, Racked, 100 – 1000 µl	10 Racks x 96 pcs (960 pcs)
PI0026	ABDOS Last Drop Low Retention Pipette Tips, Racked, Graduated, 100 – 1000 µl XL	10 Racks x 96 pcs (960 pcs)
PI0027	ABDOS Last Drop Low Retention Pipette Tips, Racked, Graduated, 100 – 1250 µl	10 Racks x 96 pcs (960 pcs)
PI0060	ABDOS Last Drop Low Retention Pipette Tips, Refill Pack, Graduated, 0.1 – 10 µl	12 Racks x 96 pcs (1152 pcs)
PI0061	ABDOS Last Drop Low Retention Pipette Tips, Refill Pack, Bevelled Graduated, 1 – 200 µl	10 Racks x 96 pcs (960 pcs)
PI0062	ABDOS Last Drop Low Retention Pipette Tips, Refill Pack, Ultra Point Tip Orifice, 1 – 200 µl	10 Racks x 96 pcs (960 pcs)
PI0063	ABDOS Last Drop Low Retention Pipette Tips, Refill Pack, Graduated, 1 – 300 µl	10 Racks x 96 pcs (960 pcs)
PI0064	ABDOS Last Drop Low Retention Pipette Tips, Refill Pack, Graduated, 100 – 1000 µl XL	6 Racks x 96 pcs (576 pcs)
PI0065	ABDOS Last Drop Low Retention Pipette Tips, Refill Pack, Graduated, 100 – 1250 µl	6 Racks x 96 pcs (576 pcs)



Last Drop Low Retention
Pipette Tips, Bulk in
Resealable Bags



Last Drop Low Retention
Pipette Tips, Racked



Last Drop Low Retention
Pipette Tips, Sterile Racked



Last Drop Low Retention
Pipette Tips, Refill Pack

- For additional product or technical information, please visit www.abdoslifesciences.com. For distributor information contact labtech@abdosindia.com.