

Technical Data Sheet

Doc No – Ab/TDS/LH-75

Rev No : 05 Date 10/2021-22

TECHNICAL DATA SHEET

Abdos Labtech Pvt. Ltd.

Rave Pipette Controller, PP

Valid for Catalog No. E11113

1. Description/Specification

1.1 Description:

Rave Pipette controller ideal for Repetitive Pipetting performance

Maximizing Pipetting comfort

1.2 Design:

Robust, Light weight & simple to use operation.

1.3 Performance Specifications:

1.3.1 Can use volumetric or serological pipettes ranging from 0.1 ml, to 100 ml.

1.3.2 It contains integral & replacement 0.3 µm Filter.

1.3.3 Autoclavable silicon pipette holder

1.4 Material : PP (polypropylene), confirm US FDA 21CFR

1.5 Color : Grey

1.6 Basic Fetures

1.6.1 Made up of highly Virgin thermoplastic polypropylene

1.6.2 Autoclavability : Yes (Only Silicon Pipette Holder – Don't autoclave whole instrument)

1.6.3 CE certified

1.6.4 Chemical resistance: See website (www.abdoslifesciences.com)



Shiv Ganga Industrial Estate,
Plot no.1 Lakeshwari, Bhagwanpur,
Roorkee – 247667 Uttrakhand, India

An ISO 9001 ,13485, 14001 & CE Certified

A Division of Abdos

Regd. Office:
13A, Government Place (East)
Kolkata – 700069, India

GST : 05AAGCA8605A1Z7

P : +91 (0) 1332 235277
E : labtech@abdosindia.com
W: www.abdoslifesciences.com
CIN : U85100WB2008PTC123903

Technical Data Sheet


1.7 Dimension:

Code	Dimension
E11113	: Length- 210 mm, ϕ -60mm
E11119	: (Spare) 3.0 μ m Filter

1.8 Quality control:

Product are manufactured, inspected and found to be in compliance with product & quality specification requirements as documented in ISO 9001:2015 & 13485:2016 Quality Management System. Testing of various characteristics in accordance with the valid specification, inspection reviewed & signed for product release.

Warranty: 12 months from the date of Installation or 15 months from the date of Invoicing


Quality Assurance Incharge
Abdos Labtech Pvt. Ltd.



ABDOS

Shiv Ganga Industrial Estate,
Plot no.1 Lakeshwari, Bhagwanpur,
Roorkee – 247667 Uttrakhand, India

An ISO 9001 ,13485, 14001 & CE Certified
A Division of Abdos

Regd. Office:
13A, Government Place (East)
Kolkata – 700069, India

GST : 05AAGCA8605A1Z7

P : +91 (0) 1332 235277

E : labtech@abdosindia.com

W: www.abdoslifesciences.com

CIN : U85100WB2008PTC123903