

Technical Data Sheet

Doc No–Ab/TDS/GLP-176

Rev No: 02 Date 01/2021

TECHNICAL DATA SHEET

Abdos Labtech Pvt. Ltd.

Animal Cage PP/Stainless Steel

Valid for Catalog No. P20501, P20502, P20503, P20504

1. Description/Specification

1.1 Description: Animal Cage PP/Stainless Steel

1.1.1 These provide excellent shelter for animal species in the laboratory.

1.1.2 Metallic top provides good ventilation & provides for feeding the animals species.

1.2 Design:

Animal Cage PP/Stainless Steel

1.3 Performance Specifications:

Code	Places	L x W x H (mm)	No. Per Case
P20501	Single Grill	290 x 220 x 140	1
P20502	Single Grill	430x270x150	1
P20503	Twin Grill	290x220x140	1
P20504	Twin Grill	430x270x150	1

1.4 Material:

1.4.1 Polypropylene

1.4.2 Complies with requirement of Directive 2002/72/EC and requirement of FDA Regulation 21 CFR olefin polymers.

1.4.3 Complies with the requirement of Indian Standard IS: 10910 on “ Specification for Polypropylene and its copolymer for its safe use in contact with food stuffs, Pharmaceutical and Drinking water”



Shiv Ganga Industrial Estate,
Plot no.1 Lakeshwari, Bhagwanpur,
Roorkee – 247667 Uttarakhand, India

An ISO 9001 ,13485, 14001 & CE Certified

A Division of Abdos

Regd. Office:
13A, Government Place (East)
Kolkata – 700069, India

GST : 05AAGCA8605A1Z7

P : +91 (0) 1332 235277

E : labtech@abdosindia.com

W: www.abdoslifesciences.com

CIN : U85100WB2008PTC123903

Technical Data Sheet

1.5 Color: White

1.6 Quality control:

Product are manufactured, inspected and found to be in compliance with product & quality specification requirements as documented in ISO 9001:2015 & 13485:2016 Quality Management System. Testing of various characteristics in accordance with the valid specification, inspection reviewed & signed for product release amendments.

1.7 Other information

1.7.1 Strong and lightweight material.

1.7.2 Good impact strength


1.7.3 Better chemical resistance.

1.8. Basic Features

1.8.1 Autoclavability: Yes

1.8.2 Chemical resistance: see website (www.abdoslifesciences.com)

1.8.3 Manufactured in Aseptic Environment



Shiv Ganga Industrial Estate,
Plot no.1 Lakeshwari, Bhagwanpur,
Roorkee – 247667 Uttrakhand, India

An ISO 9001 ,13485, 14001 & CE Certified

A Division of Abdos

Regd. Office:

13A, Government Place (East)
Kolkata – 700069, India

GST : 05AAGCA8605A1Z7

P : +91 (0) 1332 235277

E : labtech@abdosindia.com

W: www.abdoslifesciences.com

CIN : U85100WB2008PTC123903