

Technical Data Sheet

Doc No–Ab/TDS/GLP-175

Rev No: 03 Date 01/2021

TECHNICAL DATA SHEET

Abdos Labtech Pvt. Ltd.

Drying Rack

Valid for Catalog No. P20302

1. Description/Specification

1.1 Description: Drying Rack

1.1.1 These Drying Racks are used for draining & drying test tubes, beakers, flasks and other lab ware.

1.1.2 There are 55 pegs in the rack, 50 smaller pegs for products with small diameter & 5 large bracket designed to hold Erlenmeyer flasks, beakers etc.

1.1.3 The pegs are plastic coated for better protection.

1.1.4 The rack can be placed on the rack, on the table or near the sink.

1.1.5 The rack comes with a tray at the bottom on order to avoid drops on the surface or table.

1.2 Design:

Drying Rack with 55 pegs

1.3 Performance Specifications:

Code	Places	L x W x H (mm)	No. Per Case
P20302	55	360 x 160 x 680	1

1.4 Material:

Epoxy Coated Steel

1.5 Color: White



Shiv Ganga Industrial Estate,
Plot no.1 Lakeshwari, Bhagwanpur,
Roorkee – 247667 Uttrakhand, India
An ISO 9001 ,13485, 14001 & CE Certified
A Division of Abdos

Regd. Office:
13A, Government Place (East)
Kolkata – 700069, India
GST : 05AAGCA8605A1Z7

P : +91 (0) 1332 235277
E : labtech@abdosindia.com
W: www.abdoslifesciences.com
CIN : U85100WB2008PTC123903

Technical Data Sheet

1.6 Quality control:

Product are manufactured, inspected and found to be in compliance with product & quality specification requirements as documented in ISO 9001:2015 & 13485:2016 Quality Management System. Testing of various characteristics in accordance with the valid specification, inspection reviewed & signed for product release amendments.

1.7 Other information

1.7.1 Strong and lightweight material.

1.7.2 Good impact strength


1.7.3 Better chemical resistance.

1.8. Basic Features

1.8.1 Autoclavability: autoclavable, but these racks are used for drying purpose.

1.8.2 Chemical resistance: see website (www.abdoslifesciences.com)

1.8.3 Manufactured in Aseptic Environment



Shiv Ganga Industrial Estate,
Plot no.1 Lakeshwari, Bhagwanpur,
Roorkee – 247667 Uttarakhand, India

An ISO 9001 ,13485, 14001 & CE Certified

A Division of Abdos

Regd. Office:

13A, Government Place (East)
Kolkata – 700069, India

GST : 05AAGCA8605A1Z7

P : +91 (0) 1332 235277

E : labtech@abdosindia.com

W: www.abdoslifesciences.com

CIN : U85100WB2008PTC123903