

# Technical Data Sheet

Doc No– Ab/TDS/LH-72

Rev No: 03 Date 10/2021-22

## **TECHNICAL DATA SHEET**

**Abdos Labtech Pvt. Ltd.**

**Fast Release Pipette Pump**

**Valid for Catalog No. E11219, E11220, E11221**

### **1. Description/Specification**

#### **1.1 Description:** Fast Release Pipette Pump

1.1.1 Fast Release Pipette Pump ensures quick & accurate operation.

1.1.2 These light weight pipette pump have a fast release trigger that dispenses the entire pipette content.

1.1.3 The thumb wheel easily rotates for precision aspirating or dispensing, down to a volume of one drop at a time.

1.1.4 The thumb wheel can be used to blow out the contents of the pipette.

1.1.5 Resistant to acids, alkalis & easily disassembly for cleaning.

#### **1.2 Design:**

Fast Release Pipette Pump

#### **1.3 Performance Specifications:**

| <b>Code</b>   | <b>Capacity (ml)</b> | <b>No. per Case</b> |
|---------------|----------------------|---------------------|
| <b>E11219</b> | 2                    | 4                   |
| <b>E11220</b> | 10                   | 4                   |
| <b>E11221</b> | 25                   | 4                   |



**Shiv Ganga Industrial Estate,**  
Plot no.1 Lakeshwari, Bhagwanpur,  
Roorkee – 247667 Uttrakhand, India  
**An ISO 9001 ,13485, 14001 & CE Certified**  
A Division of Abdos

**Regd. Office:**  
13A, Government Place (East)  
Kolkata – 700069, India  
**GST : 05AAGCA8605A1Z7**

**P : +91 (0) 1332 235277**  
**E : labtech@abdosindia.com**  
**W: www.abdoslifesciences.com**  
**CIN : U85100WB2008PTC123903**

# Technical Data Sheet

## 1.4 Quality control:


Product are manufactured, inspected and found to be in compliance with product & quality specification requirements as documented in ISO 9001:2015 & 13485:2016 Quality Management System. Testing of various characteristics in accordance with the valid specification, inspection reviewed & signed for product release.

## 1.5. Basic Features:

**1.5.1 Autoclavability:** NO, Non autoclavable

**1.5.2 Chemical resistance:** see website ([www.abdoslifesciences.com](http://www.abdoslifesciences.com))

**1.5.3 CE certified**



**Shiv Ganga Industrial Estate,**  
Plot no.1 Lakeshwari, Bhagwanpur,  
Roorkee – 247667 Uttrakhand, India  
**An ISO 9001 ,13485, 14001 & CE Certified**  
A Division of Abdos

**Regd. Office:**  
13A, Government Place (East)  
Kolkata – 700069, India  
**GST : 05AAGCA8605A1Z7**

**P : +91 (0) 1332 235277**  
**E : labtech@abdosindia.com**  
**W: www.abdoslifesciences.com**  
**CIN : U85100WB2008PTC123903**