

# Technical Data Sheet

Doc No – Ab/TDS/UTL-216

Rev No : 07 Date 02/2022-23

## TECHNICAL DATA SHEET

**Abdos Labtech Pvt. Ltd.**

**Parafilm M Dispenser, Acrylic**

**Valid for Catalog No. P70401**

### **1. Description/Specification**

#### **1.1 Description:**

Parafilm M dispenser – Dispensing 2 “and 4” Rolls  
Ideal for storing, dispensing and cutting Parafilm® M tapes and labels.

**1.2 Design:** Rectangular Box size with easy open lid

#### **1.3 Performance Specifications:**

**1.3.1** Dispenser is designed to store, dispense and cut Parafilm®.

**1.3.2** The Parafilm® dispenser has unique design for the cutting mechanism providing maximum safety.

#### **1.4 Material:**

1.4.1 Acrylic, Hydrophobic, confirm US FDA 21CFR

1.4.2 Complies with requirement accordance with article of regulation (EC) No. 1907/2006 and its amendments.

1.4.3 Complies with the requirement of European Pharmacopoeia, 7<sup>th</sup> edition (2011) & its supplements 7.5 (07/2012) monographs 3.2.2 Plastic Containers and closures for Pharmaceutical use.

1.4.4 The product has passed the USP testing including class VI tests & has been assigned the FDA Drug master file DMF 12123

1.4.5 Materials has successful passed the biological tests according to ISO 10993 external communication devices for indirect blood contact for prolonged period.

1.4.6 Free from Heavy metals, RoHS Free according to meets relevant requirements of EU regulations 1935/2004/CE, 2002/95/EC as amended, UE 10/2011 as amended. EU CONEG & it's following amendments. It is safe use in contact with food stuff, Pharmaceuticals & drinking water

1.4.7 Product is manufactured without the use of Raw Materials of animal origin. Therefore BSE/TSE issue does not concern this product to the best of our Knowledge. No materials containing latex or natural rubber are used in the Manufacturing, handling & packaging process for this product.

**1.5 Color:** Light grey

#### **1.6 Quality control:**

Product are manufactured, inspected and found to be in compliance with product & quality specification requirements as documented in ISO 9001:2015 & 13485:2016 Quality Management System. Testing of various characteristics in accordance with the valid specification, inspection reviewed & signed for product release



**Abdos Labtech Private Limited**  
(A Division of Abdos Group)  
(An ISO 9001,13485,14001,18385,GMP & CE certified company)  
Shiv Ganga Industrial Estate  
Plot No-1, Lakheswari, Bhagwanpur,  
Roorkee-247667, (UK) INDIA

**Registered office**  
13A Government  
Place East, 2nd Floor Kolkata  
West Bengal – 700069  
CIN: U85100WB2008PTC123903

**Corporate Office**  
319, 3rd Floor, DLF Tower A, DDA Dist. Centre  
Jasola, New Delhi - 110025

**Ph. No: +91- (0) 1332235277**  
**Email:labtech@abdosindia.com**  
**Web: [www.abdoslifesciences.com](http://www.abdoslifesciences.com)**  
**GST: 05AAGCA8605A1Z7**

# Technical Data Sheet

## 1.7. Basic Features:

1.7.1 Made up of Highly Virgin Material.

1.7.2 The safety knife having unique design for cutting the film.

1.7.3 Dimensions: Length x width x Height (mm)  
180 x 120 x 160mm

1.7.4 Capacity : 2 & 4 inch

1.7.5 Chemical resistance: See website ([www.abdoslifesciences.com](http://www.abdoslifesciences.com))

1.7.6 Manufactured in Aseptic Environment

## Instruction for Use

Storage temp. Max 30°C

Before using the product check crackness, pouch punctured etc.

Check all the information printed on pouches or boxes.

Always wear gloves while handling specimen.

## WARNING

**DO NOT** store strong oxidizing agents.

**DO NOT** place any plastic labware in a flame.

**DO NOT** mix any chemicals that may result in a thermic reaction, which can cause product failure

Consult local fire codes prior to storage of flammable liquids.



**Abdos Labtech Pvt. Ltd.**  
Plot No. 1 Lakeshwari Bhagwanpur Roorkee 247667 (UK) India  
Email : [labtech@abdosindia.com](mailto:labtech@abdosindia.com), [www.abdoslifesciences.com](http://www.abdoslifesciences.com)  
Ph: +(91) 11 41081897



**Abdos Labtech Private Limited**  
(A Division of Abdos Group)  
(An ISO 9001,13485,14001,18385,GMP & CE certified company)  
Shiv Ganga Industrial Estate  
Plot No-1, Lakshewari, Bhagwanpur,  
Roorkee-247667, (UK) INDIA

**Registered office**  
13A Government  
Place East, 2nd Floor Kolkata  
West Bengal – 700069  
CIN: U85100WB2008PTC123903  
**Corporate Office**  
319, 3rd Floor, DLF Tower A, DDA Dist. Centre  
Jasola, New Delhi - 110025

Ph. No: +91- (0) 1332235277  
Email: [labtech@abdosindia.com](mailto:labtech@abdosindia.com)  
Web: [www.abdoslifesciences.com](http://www.abdoslifesciences.com)  
GST: 05AAGCA8605A1Z7