

# Technical Data Sheet

Doc No– Ab/TDS/LH-79

Rev No: 03 Date 10/2021-22

## **TECHNICAL DATA SHEET**

**Abdos Labtech Pvt. Ltd.**

**Pipette Storage Rack with Magnet, ABS**

**Valid for Catalog No. P70202**

### **1. Description**

1.1.1. Pipette Storage Rack with Magnet, ABS

1.1.2 These racks help to protect & organize glass and plastic pipettes.

1.1.3 Angled four shelved compartments keep 0.1 to 50 ml pipettes or cans within easy reach while preventing damage and clutter.

1.1.4 These have cut-out viewing window and come with two removable attachments that hold pipette bulbs, fillers and pumps.

1.1.5 These unique innovative magnetized racks have strong magnets on both sides and mounts securely to most metal surfaces.

**1.2 Design:** Pipette Storage Rack with Magnet, ABS

### **1.3 Performance Specifications:**

<b>Code</b>	<b>Description</b>	<b>L x W x H (mm)</b>	<b>No. per Case</b>
<b>P70202</b>	Storage 1-50ml pipette or cans	291x88x401	1

### **1.4 Material:**

1.4.1 Acrylonitrile Butadiene Styrene (ABS), confirm US FDA 21CFR

1.4.2 Complies with requirement accordance with article of regulation (EC) No. 1907/2006 and its amendments.



**Shiv Ganga Industrial Estate,**  
Plot no.1 Lakeshwari, Bhagwanpur,  
Roorkee – 247667 Uttarakhand, India  
**An ISO 9001 ,13485, 14001 & CE Certified**  
A Division of Abdos

**Regd. Office:**  
13A, Government Place (East)  
Kolkata – 700069, India  
**GST : 05AAGCA8605A1Z7**

**P : +91 (0) 1332 235277**  
**E : labtech@abdosindia.com**  
**W: www.abdoslifesciences.com**  
**CIN : U85100WB2008PTC123903**

# Technical Data Sheet

1.4.3 Complies with the requirement of European Pharmacopoeia, 7<sup>th</sup> edition (2011) & its supplements 7.5 (07/2012) monographs 3.2.2 Plastic Containers and closures for Pharmaceutical use.

1.4.4 The product has passed the USP testing including class VI tests & has been assigned the FDA Drug master file DMF 12123

1.4.5 Materials has successful passed the biological tests according to ISO 10993 external communication devices for indirect blood contact for Prolonged period.

1.4.6 Free from Heavy metals, RoHS Free according to meets relevant requirements of EU regulations 1935/2004/CE, 2002/95/EC as amended, UE 10/2011 as amended. EU CONEG & it's following amendments. It is safe use in contact with food stuff, Pharmaceuticals & drinking water

1.4.7 Product is manufactured without the use of Raw Materials of animal origin Therefore BSE/TSE issue does not concern this product to the best of our Knowledge. No materials containing latex or natural rubber are used in the Manufacturing, handling & packaging process for this product.

## 1.5 Color:

Blue

## 1.6 Quality control:

Product are manufactured, inspected and found to be in compliance with product & quality specification requirements as documented in ISO 9001:2015 & 13485:2016 Quality Management System. Testing of various characteristics in accordance with the valid specification, inspection reviewed & signed for product release.

## 1.6. Basic Features

**1.6.1 Autoclavability:** NO, Non autoclavable

**1.6.2 Chemical resistance:** see website ([www.abdoslifesciences.com](http://www.abdoslifesciences.com))



**Shiv Ganga Industrial Estate,**

Plot no.1 Lakeshwari, Bhagwanpur,  
Roorkee – 247667 Uttarakhand, India

**An ISO 9001, 13485, 14001 & CE Certified**

A Division of Abdos

**Regd. Office:**

13A, Government Place (East)  
Kolkata – 700069, India

**GST : 05AAGCA8605A1Z7**


**P : +91 (0) 1332 235277**

**E : labtech@abdosindia.com**

**W: www.abdoslifesciences.com**

**CIN : U85100WB2008PTC123903**

# Technical Data Sheet

  
Quality Assurance Incharge  
Abdos Labtech Pvt. Ltd.



  
**ABDOS**

**Shiv Ganga Industrial Estate,**  
Plot no.1 Lakeshwari, Bhagwanpur,  
Roorkee – 247667 Uttrakhand, India  
**An ISO 9001 ,13485, 14001 & CE Certified**  
A Division of Abdos

**Regd. Office:**  
13A, Government Place (East)  
Kolkata – 700069, India  
**GST : 05AAGCA8605A1Z7**

**P : +91 (0) 1332 235277**  
**E : labtech@abdosindia.com**  
**W: www.abdoslifesciences.com**  
**CIN : U85100WB2008PTC123903**