

Technical Data Sheet

Doc No – Ab/TDS/LH-76

Rev No : 05 Date 10/2021-22

TECHNICAL DATA SHEET

Abdos Labtech Pvt. Ltd.

Rave Pipette Filler, PP

Valid for Catalog No. E11114, E11120

1. Description/Specification

1.1 Description: Rave Pipette Filler, PP

1.1.1 Rave Plus Pipette Filler is powerful & can fill a 50 ml pipette in less than 10 seconds motor and pump operate quietly with very low vibration.

1.1.2 Users can use volumetric and serological Pipette ranging from 1 ml to 100 ml.

1.1.3 Rave Pipette Controller is supplied with integral and replaceable 0.45 µm filler.

1.1.4 E11114 – with LCD screen, Rave Plus Pipette Filler, PP

1.1.4 E11120 – NO LCD screen, Rave ECO Pipette Filler, PP

1.2 Design:

Robust, Lightweight & simple to use operation.

1.3 Performance Specifications:

Battery	Lithium-Ion
Battery Service Life	More than 8 hours of intermittence use.
Charging Time	Approx. 2 hour
Pipette Types	Glass & Plastic pipette (0.1-100ml)
Aspirate Speeds	8
Dispense Speed	8 Gravity Dispense
Filter	0.45 µm Hydrophobic

1.4 Material: PP (polypropylene)

1.5 Color: Grey & White

1.6 Basic Features

1.6.1 Made up of highly Virgin thermoplastic polypropylene

1.6.2 Autoclavability : Yes (Only Silicon Pipette Holder – Don't autoclave whole instrument)



Shiv Ganga Industrial Estate,
Plot no.1 Lakeshwari, Bhagwanpur,
Roorkee – 247667 Uttarakhand, India

An ISO 9001 ,13485, 14001 & CE Certified

A Division of Abdos

Regd. Office:
13A, Government Place (East)
Kolkata – 700069, India

GST : 05AAGCA8605A1Z7

P : +91 (0) 1332 235277

E : labtech@abdosindia.com

W: www.abdoslifesciences.com

CIN : U85100WB2008PTC123903

Technical Data Sheet

1.6.3 CE certified

1.6.4 **Chemical resistance:** See website (www.abdoslifesciences.com)


1.7. Spare Filters

Code	Product	No. Per Case
E11116	Spare 0.45 µm Filter	1
E11117	Spare Lithium-ion	1
E11118	Spare Adaptor	1

1.8 Quality control:

Product are manufactured, inspected and found to be in compliance with product & quality specification requirements as documented in ISO 9001:2015 & 13485:2016 Quality Management System. Testing of various characteristics in accordance with the valid specification, inspection reviewed & signed for product release.

Warranty: 12 months from the date of Installation or 15 months from the date of Invoicing


Quality Assurance Incharge
Abdos Labtech Pvt. Ltd.



ABDOS

Shiv Ganga Industrial Estate,
Plot no.1 Lakeshwari, Bhagwanpur,
Roorkee – 247667 Uttrakhand, India

An ISO 9001 ,13485, 14001 & CE Certified

A Division of Abdos

Regd. Office:

13A, Government Place (East)
Kolkata – 700069, India

GST : 05AAGCA8605A1Z7

P : +91 (0) 1332 235277

E : labtech@abdosindia.com

W: www.abdoslifesciences.com

CIN : U85100WB2008PTC123903