

Technical Data Sheet

Doc No– Ab/TDS/LH-74

Rev No: 04 Date 10/2021-22

TECHNICAL DATA SHEET

Abdos Labtech Pvt. Ltd.

Rubber Pipette Pump

Valid for Catalog No. E11217

1. Description/Specification

1.1 Description: Rubber Pipette Pump

1.1.1 Rubber Pipette Pump ensures quick & accurate operation.

1.1.2 Pump opening is standard to fit almost all glass & plastic pipette.

1.1.3 One size can be used for liquid up to 100 ml.

1.1.4 The pump has a heavy wall and easy to grip body.

1.2 Design: Rubber Pipette Pump

1.3 Performance Specifications:

Code	Capacity (ml)	No. per Case
E11217	Up to 100ml	4

1.4 Quality control:

Product are manufactured, inspected and found to be in compliance with product & quality specification requirements as documented in ISO 9001:2015 & 13485:2016 Quality Management System. Testing of various characteristics in accordance with the valid specification, inspection reviewed & signed for product release

1.5. Basic Features:

1.5.1 Autoclavability: NO, Non autoclavable

1.5.2 Chemical resistance: see website (www.abdoslifesciences.com)

1.5.3 Colour- Available in various colour




Shiv Ganga Industrial Estate,
Plot no.1 Lakeshwari, Bhagwanpur,
Roorkee – 247667 Uttrakhand, India
An ISO 9001 ,13485, 14001 & CE Certified
A Division of Abdos

Regd. Office:
13A, Government Place (East)
Kolkata – 700069, India
GST : 05AAGCA8605A1Z7

P : +91 (0) 1332 235277
E : labtech@abdosindia.com
W: www.abdoslifesciences.com
CIN : U85100WB2008PTC123903

Technical Data Sheet


Quality Assurance Incharge
Abdos Labtech Pvt. Ltd.




ABDOS

Shiv Ganga Industrial Estate,
Plot no.1 Lakeshwari, Bhagwanpur,
Roorkee – 247667 Uttrakhand, India
An ISO 9001 ,13485, 14001 & CE Certified
A Division of Abdos

Regd. Office:
13A, Government Place (East)
Kolkata – 700069, India
GST : 05AAGCA8605A1Z7

P : +91 (0) 1332 235277
E : labtech@abdosindia.com
W: www.abdoslifesciences.com
CIN : U85100WB2008PTC123903