

Technical Data Sheet

Doc No– Ab/TDS/Eqpt -240

Rev No : 04 Date 02/2022-23

TECHNICAL DATA SHEET

Abdos Labtech Pvt. Ltd.

SWIRLEX – Microplate Mixer

Valid for Catalog No. E11200

1. Description/Specification

1.1 Description: SWIRLEX - Microplate Mixer

1.1.1 Swirlex Microplate Mixer is used for holding single or double micro plates with an SBS footprint.

1.1.2 Reliability, long continuous service and outstanding performance are its hallmark.

1.1.3 This is widely used in life science including microbiology, cell & molecular biology etc.

1.1.4 The high speed and small diameter combine to offer a true vortexing action in each well of the microplate. This action ensures both complete and uniform mixing.

1.2 Design: SWIRLEX - Microplate Mixer

1.3 Performance Specifications:

Specifications	E11200
Voltage	220-240V
Frequency	50-60Hz
Speed Range	Single Microplate :0-1500 rpm Double Microplate :0-1000 rpm
Dimension	260x150x80 mm
Weight	2.6 kg
Permissible ambient Temperature	5 - 40°C
Permissible ambient	80%
Orbital diameter	4.5mm
Shaking Movement	Orbital
Motor Type	Brushless DC Motor
Max. shaking weight	0.5kg
Protection class according	IP21



Abdos Labtech Private Limited
(A Division of Abdos Group)
(An ISO 9001,13485,14001,18385,GMP & CE certified company)
Shiv Ganga Industrial Estate
Plot No-1, Lakheswari, Bhagwanpur,
Roorkee-247667, (UK) INDIA

Registered office
13A Government
Place East, 2nd Floor Kolkata
West Bengal – 700069
CIN: U85100WB2008PTC123903
Corporate Office
319, 3rd Floor, DLF Tower A, DDA Dist. Centre
Jasola, New Delhi - 110025

Ph. No: +91- (0) 1332235277
Email:labtech@abdosindia.com
Web: www.abdoslifesciences.com
GST: 05AAGCA8605A1Z7

Technical Data Sheet

Code	Product
E11200	SWIRLEX - Microplate Mixer including Microplate Clamp(E11201) as standard

Spare Accessories:

Code	Product
E11201	Microplate Clamp
E11202	Double Microplate Clamp

1.4 Material: Body – ABS Material

1.5 Color: Grey

1.6. Basic features:

1.6.1 **Material:** Made up of highly pure Virgin ABS (Acrylonitrile Butadiene Styrene)

1.6.2 CE certified

1.6.3 **Chemical resistance:** see website (www.abdoslifesciences.com)

Warranty: 12 months from the date of Installation or 15 months from the date of Invoicing

Instruction for Use

Storage temp. Max 30°C

Before using the product check crackness, etc.

Check all the information printed on pouches or boxes.

WARNING

DO NOT store strong oxidizing agents.

DO NOT place any plastic labware in a flame.

DO NOT mix any chemicals that may result in a thermic reaction, which can cause product failure

Consult local fire codes prior to storage of flammable liquids.



Abdos Labtech Private Limited
(A Division of Abdos Group)
(An ISO 9001,13485,14001,18385,GMP & CE certified company)
Shiv Ganga Industrial Estate
Plot No-1, Lakheswari, Bhagwanpur,
Roorkee-247667, (UK) INDIA

Registered office
13A Government
Place East, 2nd Floor Kolkata
West Bengal – 700069
CIN: U85100WB2008PTC123903
Corporate Office
319, 3rd Floor, DLF Tower A, DDA Dist. Centre
Jasola, New Delhi - 110025

Ph. No: +91- (0) 1332235277
Email:labtech@abdosindia.com
Web: www.abdoslifesciences.com
GST: 05AAGCA8605A1Z7

Technical Data Sheet



Abdos Labtech Pvt. Ltd.
Plot No. 1 Lakshwari Bhagwanpur Roorkee 247667 (UK) India
Email : labtech@abdosindia.com, www.abdoslifesciences.com
Ph: +(91) 11 41082897



Am
Quality Assurance Incharge
Abdos Labtech Pvt. Ltd.



Abdos Labtech Private Limited
(A Division of Abdos Group)
(An ISO 9001,13485,14001,18385,GMP & CE certified company)
Shiv Ganga Industrial Estate
Plot No-1, Lakshwari, Bhagwanpur,
Roorkee-247667, (UK) INDIA

Registered office
13A Government
Place East, 2nd Floor Kolkata
West Bengal – 700069
CIN: U85100WB2008PTC123903
Corporate Office
319, 3rd Floor, DLF Tower A, DDA Dist. Centre
Jasola, New Delhi - 110025

Ph. No: +91- (0) 1332235277
Email:labtech@abdosindia.com
Web: www.abdoslifesciences.com
GST: 05AAGCA8605A1Z7