

Technical Data Sheet

Doc No–Ab/TDS/UTL-213

Rev No: 05 Date 02/2022-23

TECHNICAL DATA SHEET

Abdos Labtech Pvt. Ltd.

Wypall X70 Wipers

Valid for Catalog No. U20407, U20413, U20414, U20415

1. Description/Specification

1.1 Description: Wypall X70 Wipers

1.1.1 An ideal product to replace cloth towels or rags.

1.1.2 They can exceptionally absorb oil, grease, solvent & water.

1.1.3 They are designed to perform the wiping function that helps to improve productivity and performance.

1.1.4 They are made from advanced HYDROKNIT material that delivers products that are absorbent, higher in strength and have more bulk.

1.1.5 They are bio-degradable hence reduce disposables coast and are environmental friendly.

1.1.6 Available in 100 meters rood format in white and blue color with perforation making it easy to use.

1.2 Design: Rolled wipers

1.3 Performance Specifications:

Code	Size LxW (cm)	Colour	No. Per Case
U20407	24.7 x 20	White	1roll-100mtr
U20413	24.7 x 20	Blue	1roll-100mtr
U20414	26 x 26	White	50 Sheets
U20415	25.4 x 38.1	Blue	50 Sheets

1.4 Color: White

1.5 Quality control:

Product are manufactured, inspected and found to be in compliance with product & quality specification requirements as documented in ISO 9001:2015 & 13485:2016 Quality Management System. Testing of various characteristics in accordance with the valid specification, inspection reviewed & signed for product release.



Abdos Labtech Private Limited
(A Division of Abdos Group)
(An ISO 9001,13485,14001,18385,GMP & CE certified company)
Shiv Ganga Industrial Estate
Plot No-1, Lakheswari, Bhagwanpur,
Roorkee-247667, (UK) INDIA

Registered office
13A Government
Place East, 2nd Floor Kolkata
West Bengal – 700069
CIN: U85100WB2008PTC123903
Corporate Office
319, 3rd Floor, DLF Tower A, DDA Dist. Centre
Jasola, New Delhi - 110025

Ph. No: +91- (0) 1332235277
Email:labtech@abdosindia.com
Web: www.abdoslivesciences.com
GST: 05AAGCA8605A1Z7

Technical Data Sheet

1.6 Other information

1.6.1 These wipers are used for cleaning & wiping delicate surfaces.

1.6.2 These wipers are soft and wipes easily for liquid & dust on the surfaces.

1.6.3 These Wipers are suitable for cleaning & wiping delicate surfaces.

1.6.4 **Chemical resistance:** See website (www.abdoslifesciences.com)

1.6.5 Manufactured in Aseptic Environment

Instruction for Use

Storage temp. Max 30°C

Before using the product check crackness etc.

Check all the information printed on pouches or boxes.

Always wear gloves while handling specimen.


WARNING

DO NOT store strong oxidizing agents.

DO NOT place any plastic labware in a flame.

DO NOT mix any chemicals that may result in a thermic reaction, which can cause product failure

Consult local fire codes prior to storage of flammable liquids.



Abdos Labtech Pvt. Ltd.

Plot No. 1 Lakeshwari Bhagwanpur Roorkee 247667 (UK) India
Email : labtech@abdosindia.com, www.abdoslifesciences.com
Ph: +(91) 11 41081897

