

ABDOS TIPS & LAST DROP PIPETTE TIPS

COMPATIBILITY CHART

Pipettor Brand / Make			Abdos Premium (Single Channel)					Abdos Premium Plus (Single Channel)					Abdos Fixed (Single Channel)					Abdos Premium (Multi Channel)			Abdos Premium Plus (Multi Channel)							
Abdos Code for Pipette Tips	Abdos Code for Last Drop Pipette Tips	Volume	0.5 - 10 µl	2-20 µl	10-100 µl	20-200µl	100-1000 µl	0.5 - 10 µl	2-20 µl	10-100 µl	20-200 µl	100-1000 µl	5 µl	10 µl	20 µl	25 µl	50 µl	100 µl	200 µl	250 µl	500 µl	1000 µl	0.5 - 10 µl	5 - 50 µl	50-300 µl	0.5 - 10 µl	5 - 50 µl	50-300 µl
			P10115	P10001	0.2 - 10 µl	✓					✓				✓	✓										✓		
P10137	P10002	0.2 - 10 µl	✓					✓				✓	✓										✓			✓		
P10129	-	2 - 200 µl		✓	✓	✓			✓	✓	✓			✓	✓	✓	✓	✓	✓	✓	✓			✓	✓		✓	✓
P10130	P10003	2 - 200 µl		✓	✓	✓			✓	✓	✓			✓	✓	✓	✓	✓	✓	✓	✓			✓	✓		✓	✓
P10101	-	2 - 200 µl		✓	✓	✓			✓	✓	✓			✓	✓	✓	✓	✓	✓	✓	✓			✓	✓		✓	✓
P10105	P10004	2 - 200 µl		✓	✓	✓			✓	✓	✓			✓	✓	✓	✓	✓	✓	✓	✓			✓	✓		✓	✓
P10138	P10005	2 - 300 µl		✓	✓	✓			✓	✓	✓			✓	✓	✓	✓	✓	✓	✓	✓			✓	✓		✓	✓
P10102	-	100 - 1000 µl					✓														✓	✓						
P10106	P10006	100 - 1000 µl					✓														✓	✓						

Pipettor Brand/ Make			Biohit eLine (Single channel)				Biohit eLine (multi channel)				Biohit m-Line® (Single channel)				Biohit m-Line® (multi channel)				Biohit Proline (Single channel)				Biohit Proline (Multi channel)				Biohit ePet (Single channel)	
Abdos Code for Pipette Tips	Abdos Code for Last Drop Pipette Tips	Volume	0.2-10µl	5-120µl	10-300µl	50-1000µl	0.2-10µl	5-120µl	10-300µl	0.5-10µl	2-20µl	10-100µl	20-200µl	100-1000µl	0.5-10µl	5-100µl	30-300µl	0.2-10µl	5-100µl	10-250µl	10-500µl	50-1000µl	0.2-10µl	5-100µl	25-250µl	0.2-10µl	5-100µl	50-1200µl
			P10115	P10001	0.2 - 10 µl	✓				✓		✓					✓			✓					✓			✓
P10137	P10002	0.2 - 10 µl	✓				✓		✓					✓			✓					✓			✓			
P10129	-	2 - 200 µl		✓							✓	✓			✓				✓					✓			✓	
P10130	P10003	2 - 200 µl		✓							✓	✓			✓				✓					✓			✓	
P10101	-	2 - 200 µl		✓							✓	✓			✓				✓					✓			✓	
P10105	P10004	2 - 200 µl		✓							✓	✓			✓				✓					✓			✓	
P10138	P10005	2 - 300 µl		✓	✓			✓			✓	✓			✓	✓		✓	✓	✓	✓			✓	✓		✓	
P10102	-	100 - 1000 µl				✓							✓									✓						
P10106	P10006	100 - 1000 µl				✓							✓									✓						

Pipettor Brand/ Make			Biohit ePet (multi channel)				Biohit Proline electronics (Single channel)				Biohit Proline electronics (Multi channel)				Biohit Proline (Single channel)				Biohit Proline (Multi channel)				Biohit Proline Plus (Single channel)					
Abdos Code for Pipette Tips	Abdos Code for Last Drop Pipette Tips	Volume	0.2-10µl	5-100µl	50-1200µl	25-250 µl	0.2-10µl	5-100µl	10-250µl	50-1000µl	50-1200µl	0.2-10µl	5-100µl	10-250µl	50-1200µl	0.2-10µl	5-50µl	10-100µl	50-200µl	10-250µl	200-1000µl	0.5-10µl	5-50µl	50-250µl	50-300µl	0.5-10µl	20-200µl	100-1000 µl
			P10115	P10001	0.2 - 10 µl	✓				✓				✓				✓						✓				✓
P10137	P10002	0.2 - 10 µl	✓				✓				✓				✓						✓				✓			
P10129	-	2 - 200 µl		✓				✓				✓				✓	✓	✓				✓				✓		
P10130	P10003	2 - 200 µl		✓				✓				✓				✓	✓	✓				✓				✓		
P10101	-	2 - 200 µl		✓				✓				✓				✓	✓	✓				✓				✓		
P10105	P10004	2 - 200 µl		✓				✓				✓				✓	✓	✓				✓				✓		
P10138	P10005	2 - 300 µl		✓	✓		✓	✓			✓	✓			✓	✓	✓	✓	✓			✓	✓	✓		✓		
P10102	-	100 - 1000 µl							✓											✓								✓
P10106	P10006	100 - 1000 µl							✓											✓								✓

Pipettor Brand/ Make			Brand Transferpette												Brand Electronics												
			Single Channel						Multi Channel						Single Channel			Multi Channel									
Abdos Code for Pipette Tips	Abdos Code for Last Drop Pipette Tips	Volume	0.1-1µl	0.5-10µl	2-20µl 'ULTRA MICRO'	2-20µl	5-50µl	10-100µl	20-200µl	25-250µl	100-1000µl	0.5-10µl	2-20µl	2.5-25µl	5-50µl	10-100µl	20-200µl	0.5-10µl	10-100µl	20-200µl	30-300µl	100-1000µl	0.5-10µl	5-100µl	10-200µl	15-300µl	
P10115	P10001	0.2 - 10 µl		√								√						√					√				
P10137	P10002	0.2 - 10 µl		√								√						√					√				
P10129	-	2 - 200 µl					√	√	√					√	√	√			√	√					√	√	
P10130	P10003	2 - 200 µl					√	√	√					√	√	√			√	√					√	√	
P10101	-	2 - 200 µl					√	√	√					√	√	√			√	√					√	√	
P10105	P10004	2 - 200 µl					√	√	√					√	√	√			√	√					√	√	
P10138	P10005	2 - 300 µl					√	√	√	√				√	√	√	√		√	√	√				√	√	√
P10102	-	100 - 1000 µl									√												√				
P10106	P10006	100 - 1000 µl									√												√				

Pipettor Brand/ Make			Brand Transferpette S						Capp Aero						Gilson Pipetman				Gilson Pipetman Ultra									
			Single Channel			Multi Channel			Single Channel / Multi Channel						Single Channel				Multi Channel									
Abdos Code for Pipette Tips	Abdos Code for Last Drop Pipette Tips	Volume	0.1-1µl	0.5-10µl	2-20µl 'ULTRA MICRO'	10-100µl	20-200µl	100-1000µl	0.5-10µl	10-100µl	20-200µl	30-300µl	0.2-2µl	0.5-10µl	5-50µl	10-100µl	25-200µl	50-300µl	100-1000µl	0.2-2µl	1-10µl	2-20µl	20-100µl	50-200µl	200-1000µl	1-20µl	20-300µl	
P10115	P10001	0.2 - 10 µl		√				√					√	√						√	√							
P10137	P10002	0.2 - 10 µl		√				√					√	√						√	√							
P10129	-	2 - 200 µl				√	√		√	√					√	√							√	√				
P10130	P10003	2 - 200 µl				√	√		√	√					√	√							√	√				
P10101	-	2 - 200 µl				√	√		√	√					√	√							√	√				
P10105	P10004	2 - 200 µl				√	√		√	√					√	√							√	√				
P10138	P10005	2 - 300 µl				√	√		√	√	√				√	√	√					√	√				√	
P10102	-	100 - 1000 µl						√											√						√			
P10106	P10006	100 - 1000 µl						√											√						√			

Pipettor Brand / Make			Gilson Pipetman Ultra					Gilson Concept			Gilson Concept			Gilson Neo				
			Single Channel					Single Channel			Multi Channel			Single Channel				
Abdos Code for Pipette Tips	Abdos Code for Last Drop Pipette Tips	Volume	0.2-2µl	1-10µl	2-20µl	20-100µl	20-200µl	200-1000µl	0.5-10µl	5-100µl	20-300µl	1-10µl	10-100µl	30-300µl	0.2-2µl	1-10µl	20-100µl	100-1000µl
P10115	P10001	0.2 - 10 µl		√					√			√				√		
P10137	P10002	0.2 - 10 µl		√					√			√				√		
P10129	-	2 - 200 µl				√	√			√			√				√	
P10130	P10003	2 - 200 µl				√	√			√			√				√	
P10101	-	2 - 200 µl				√	√			√			√				√	
P10105	P10004	2 - 200 µl				√	√			√			√				√	
P10138	P10005	2 - 300 µl				√	√			√	√		√	√			√	
P10102	-	100 - 1000 µl						√										√
P10106	P10006	100 - 1000 µl						√										√

Pipettor Brand/ Make			Vistalab Ovation QS					Vistalab Ovation M					Vistalab Ovation E				Eppendorf Reference				"Eppendorf Electronic"								
			Single Channel					Single Channel					Multi Channel				Single Channel												
Abdos Code for Pipette Tips	Abdos Code for Last Drop Pipette Tips	Volume	0.2-2µL	1-10µL	2-20µL	10-100µL	20-200µL	100-1000µL	0.2-2µL	1-10µL	2-20µL	10-100µL	20-200µL	100-1000µL	0.5-20µL	5-250µL	2-125µL	25-1250µL	0.1-10µL	"0.1-10µL ULTRA MICRO"	2-20µL	10-100µL	20-200µL	100-1000µL	0.5-10µL	5-100µL	2-200µL	50-1000µL	
P10115	P10001	0.2 - 10 µl	√	√	√			√	√	√				√				√		√									
P10137	P10002	0.2 - 10 µl	√	√	√			√	√	√				√				√		√									
P10129	-	2 - 200 µl				√	√					√	√		√	√						√	√			√	√		
P10101	-	2 - 200 µl				√	√					√	√		√	√						√	√			√	√		
P10105	P10004	2 - 200 µl				√	√					√	√		√	√						√	√			√	√		
P10138	P10005	2 - 300 µl				√	√					√	√		√	√						√	√			√	√		
P10102	-	100 - 1000 µl					√						√				√							√				√	
P10106	P10006	100 - 1000 µl					√						√				√							√				√	

Pipettor Brand/ Make			Eppendorf Research					"Eppendorf "					"Thermo Finnpiquette"				"Thermo Finnpiquette"		"Matrix Impact/Impact2"										
			Single Channel					Multi Channel					Single Channel				Multi Channel		s /mc										
Abdos Code for Pipette Tips	Abdos Code for Last Drop Pipette Tips	Volume	0.1-2.5µl	0.5-10µl	2-20µl	10-100µl	20-200µl	100-1000µl	0.5-10µl	5-100µl	50-100µl	0-300µl	30-300µl	50-1200µl	0.2-10µl	0.5-10µl	5-40µl	40-200µl	100-1000µl	200-1000µl	0.1-10µl	10-100µl	30-300µl	2-125µl	2-20µl	15-1250µl			
P10115	P10001	0.2 - 10 µl		√	√				√						√	√					√								
P10137	P10002	0.2 - 10 µl		√	√				√						√	√					√								
P10129	-	2 - 200 µl				√	√			√	√							√				√			√				
P10101	-	2 - 200 µl				√	√			√	√							√				√			√				
P10105	P10004	2 - 200 µl				√	√			√	√							√				√			√				
P10138	P10005	2 - 300 µl				√	√			√	√	√	√					√				√	√		√				
P10102	-	100 - 1000 µl					√						√*						√	√									
P10106	P10006	100 - 1000 µl					√						√*						√	√									

Pipettor Brand/ Make			Labnet Labpette					Labnet Labpette				Labnet BioPette E				"Nichiryo / Oxford Benchmate"													
			Single Channel					Multi Channel				Single Channel		Multi Channel		Single Channel & Multi Channel													
Abdos Code for Pipette Tips	Abdos Code for Last Drop Pipette Tips	Volume	0.1-2µl	0.5-10µl	2-20µl	10-100µl	20-200µl	100-1000µl	1-10µl	5-50µl	20-200µl	50-300µl	0.5-10µl	2-20µl	10-200µl	100-1000µl	0.5-10µl	2-20µl	10-200µl	100-1000µl	0.5-10µl	2-20µl	10-50µl	10-100µl	20-200µl	100-1000µl			
P10115	P10001	0.2 - 10 µl	√	√					√				√			√					√								
P10137	P10002	0.2 - 10 µl	√	√					√				√			√					√								
P10129	-	2 - 200 µl				√	√							√					√					√	√				
P10101	-	2 - 200 µl				√	√							√					√					√	√				
P10105	P10004	2 - 200 µl				√	√							√					√					√	√				
P10138	P10005	2 - 300 µl				√	√							√					√					√	√				
P10102	-	100 - 1000 µl					√							√					√								√		
P10106	P10006	100 - 1000 µl					√							√					√								√		

Pipettor Brand/ Make			Socorex Calibra /Accura												Axygen Axypet												
			Single Channel						Multi Channel						Single Channel				Multi Channel								
Abdos Code for Pipette Tips	Abdos Code for Last Drop Pipette Tips	Volume	0.1-20µl	0.2-20µl	0.5-10µl	1-10µl	2-20µl	10-100µl	20-200µl	100-1000µl	1-10µl	5-50µl	10-100µl	20-200µl	50-200µl	0.1-2µl	0.5-10µl	2-20µl	5-50µl	10-100µl	20-200µl	50-250µl	100-1000µl	1-10µl	5-50µl	20-200µl	50-300µl
P10115	P10001	0.2 - 10 µl			√	√					√					√								√			
P10137	P10002	0.2 - 10 µl			√	√					√					√								√			
P10129	-	2 - 200 µl						√	√		√	√	√	√				√	√	√					√	√	√
P10130	P10003	2 - 200 µl						√	√		√	√	√	√				√	√	√					√	√	√
P10101	-	2 - 200 µl						√	√		√	√	√	√				√	√	√					√	√	√
P10105	P10004	2 - 200 µl						√	√		√	√	√	√				√	√	√					√	√	√
P10138	P10005	2 - 300 µl						√	√		√	√	√	√				√	√	√	√				√	√	√
P10102	-	100 - 1000 µl								√													√				
P10106	P10006	100 - 1000 µl								√													√				

Pipettor Brand/ Make		BRAND		Eppendorf		Biohit Proline plus		Gilson		Thermo		Abdos	
Abdos Code for Pipette Tips	Volume	1-5 ml	1-10 ml	1-5 ml	1-10 ml	1-5 ml	1-10 ml	1-5 ml	1-10 ml	1-5 ml	1-10 ml	1-5 ml	1-10 ml
P10144	1 - 5 ml			√		√		√				√	
P10145	1 - 10 ml		√		√		√		√				√

Pipettor Brand/ Make			Abdos premium	Abdos premium Plus	ABDOS D. Pipette (Electronic)	Brand Transferpette S	Gilson Pipetman	Eppendorf Research
			Single channel			Single channel		
Abdos Code for Pipette Tips	Abdos Code for Last Drop Pipette Tips	Volume	100-1000µl	100-1000µl	100-1000µl	100-1000µl	100-1000µl	100-1000µl
P10174	P10026	100 - 1000 µl XL	√	√	√	√	√	√
P10175	P10027	100 - 1250 µl	√	√	√	√	√	√

Pipettor Brand/ Make			THERMO (Electronic) FinnpiPETTE	Thermo FinnpiPETTE F3	Nichiryo NICHIPET EX	SOCOREX ACURA 825	BIOHIT Proline	ErgoOne Star Lab
			Multi channel	Single channel				
Abdos Code for Pipette Tips	Abdos Code for Last Drop Pipette Tips	Volume	100-1200µl	100-1000µl	100-1000µl	100-1000µl	100-1000µl	100-1000µl
P10174	P10026	100 - 1000 µl XL	√	√	√	√	√	√
P10175	P10027	100 - 1250 µl	√	√	√	√	√	√

ABDOS FILTER TIPS & LAST DROP FILTER PIPETTE TIPS

COMPATIBILITY CHART

Pipettor Brand/ Make			Abdos Premium (Single Channel)					Abdos Premium Plus (Single Channel)					Abdos Fixed (Single Channel)					Abdos Premium (Multi Channel)			Abdos Premium Plus (Multi Channel)							
Abdos Code for Filter Pipette Tips	Abdos Code for Last Drop Filter Pipette Tips	Volume	0.5 - 10 µl	2-20 µl	10-100 µl	20-200µl	100-1000 µl	0.5 - 10 µl	2-20 µl	10-100 µl	20-200 µl	100-1000 µl	5 µl	10 µl	20 µl	25 µl	50 µl	100 µl	200 µl	250 µl	500 µl	1000 µl	0.5 - 10 µl	5 - 50 µl	50-300 µl	0.5 - 10 µl	5 - 50 µl	50-300 µl
P11001	P11050	0.2 - 10 µl	√					√				√	√										√			√		
P11002	P11051	0.2 - 10 µl	√					√				√	√										√			√		
P11003	P11052	2 - 20 µl		√	√	√			√	√	√				√	√	√	√	√	√			√	√		√	√	√
P11026	P11066	2 - 50 µl		√	√	√			√	√	√				√	√	√	√	√	√			√	√		√	√	√
P11004	P11053	2 - 100 µl		√	√	√			√	√	√				√	√	√	√	√	√			√	√		√	√	√
P11005	P11054	2 - 200 µl		√	√	√			√	√	√				√	√	√	√	√	√			√	√		√	√	√
P11006	P11055	2 - 300 µl		√	√	√			√	√	√				√	√	√	√	√	√			√	√		√	√	√
P11007	P11056	100 - 1000 µl					√					√									√	√						

Pipettor Brand / Make			Biohit eLine (Single Channel)				Biohit eLine (Multi Channel)				Biohit m- Line® (Single channel)				Biohit m-Line® (Multi Channel)				Biohit Proline (Single Channel)				Biohit Proline (Multi Channel)				Biohit ePet (Single Channel)	
Abdos Code for Filter Pipette Tips	Abdos Code for Last Drop Filter Pipette Tips	Volume	0.2-10µl	5-120µl	10-300µl	50-1000µl	0.2-10µl	5-120µl	10-300µl	0.5-10µl	2-20µl	10-100µl	20-200µl	100-1000µl	0.5-10µl	5-100µl	30-300µl	0.2-10µl	5-100µl	10-250µl	10-500µl	50-1000µl	0.2-10µl	5-100µl	25-250µl	0.2-10µl	5-100µl	50-1200µl
P11001	P11050	0.2 - 10 µl	√				√		√					√			√						√			√		
P11002	P11051	0.2 - 10 µl	√				√		√					√			√						√			√		
P11003	P11052	2 - 20 µl		√							√	√			√			√						√			√	
P11026	P11066	2 - 50 µl		√							√	√			√			√						√			√	
P11004	P11053	2 - 100 µl		√							√	√			√			√						√			√	
P11005	P11054	2 - 200 µl		√							√	√			√			√						√			√	
P11006	P11055	2 - 300 µl		√	√			√			√	√			√	√		√	√				√	√		√	√	√
P11007	P11056	100 - 1000 µl				√							√									√						

Pipettor Brand/ Make			Biohit ePet Multi Channel				Biohit Proline electronics (Single Channel)				Biohit Proline Electronics (Multi Channel)				Biohit Proline (Single Channel)				Biohit Proline (Multi Channel)				Biohit Proline Plus (Single Channel)					
Abdos Code for Filter Pipette Tips	Abdos Code for Last Drop Filter Pipette Tips	Volume	0.2-10µl	5-100µl	50-1200µl	25-250 µl	0.2-10µl	5-100µl	10-250µl	50-1000µl	50-1200µl	0.2-10µl	5-100µl	10-250µl	50-1200µl	0.2-10µl	5-50µl	10-100µl	50-200µl	10-250µl	200-1000µl	0.5-10µl	5-50µl	50-250µl	50-300µl	0.5-10µl	20-200µl	100-1000 µl
P11001	P11050	0.2 - 10 µl	√				√				√				√							√				√		
P11002	P11051	0.2 - 10 µl	√				√				√				√							√				√		
P11003	P11052	2 - 20 µl		√				√				√				√	√	√					√				√	
P11026	P11066	2 - 50 µl		√				√				√				√	√	√					√				√	
P11004	P11053	2 - 100 µl		√				√				√				√	√	√					√				√	
P11005	P11054	2 - 200 µl		√				√				√				√	√	√					√				√	
P11006	P11055	2 - 300 µl		√	√		√	√			√	√			√	√	√	√	√				√	√	√	√	√	√
P11007	P11056	100 - 1000 µl							√											√								√

Pipettor Brand/ Make			Brand Transferpette												Brand Electronics												
			Single Channel						Multi Channel						Single Channel			Multi Channel									
Abdos Code for Filter Pipette Tips	Abdos Code for Last Drop Filter Pipette Tips	Volume	0.1-1µl	0.5-10µl	2-20µl 'ULTRA MICRO'	2-20µl	5-50µl	10-100µl	20-200µl	25-250µl	100-1000µl	0.5-10µl	2-20µl	2.5-25µl	5-50µl	10-100µl	20-200µl	30-300µl	0.5-10µl	10-100µl	20-200µl	30-300µl	100-1000µl	0.5-10µl	5-100µl	10-200µl	15-300µl
P11001	P11050	0.2 - 10 µl		√								√							√					√			
P11002	P11051	0.2 - 10 µl		√								√							√					√			
P11003	P11052	2 - 20 µl					√	√	√						√	√	√			√	√					√	√
P11026	P11066	2 - 50 µl					√	√	√						√	√	√			√	√					√	√
P11004	P11053	2 - 100 µl					√	√	√						√	√	√			√	√					√	√
P11005	P11054	2 - 200 µl					√	√	√						√	√	√			√	√					√	√
P11006	P11055	2 - 300 µl					√	√	√	√				√	√	√	√			√	√	√			√	√	√
P11007	P11056	100 - 1000 µl								√													√				

Pipettor Brand/ Make			Brand Transferpette S						Capp Aero						Gilson Pipetman			Gilson Pipetman Ultra										
			Single Channel			Multi Channel			Single Channel /Multi Channel						Single Channel			Multi Channel										
Abdos Code for Filter Pipette Tips	Abdos Code for Last Drop Filter Pipette Tips	Volume	0.1-1µl	0.5-10µl	2-20µl 'ULTRA MICRO'	10-100µl	20-200µl	100-1000µl	0.5-10µl	10-100µl	20-200µl	30-300µl	0.2-2µl	0.5-10µl	5-50µl	10-100µl	25-200µl	50-300µl	100-1000µl	0.2-2µl	1-10µl	2-20µl	20-100µl	50-200µl	200-1000µl	1-20µl	20-300µl	
P11001	P11050	0.2 - 10 µl		√					√				√	√						√	√							
P11002	P11051	0.2 - 10 µl		√					√				√	√						√	√							
P11003	P11052	2 - 20 µl				√	√			√	√					√	√						√	√				
P11026	P11066	2 - 50 µl				√	√			√	√					√	√						√	√				
P11004	P11053	2 - 100 µl				√	√			√	√					√	√						√	√				
P11005	P11054	2 - 200 µl				√	√			√	√					√	√						√	√				
P11006	P11055	2 - 300 µl				√	√			√	√	√				√	√	√					√	√			√	
P11007	P11056	100 - 1000 µl						√											√					√				

Pipettor Brand/ Make			Gilson Pipetman Ultra					Gilson Concept			Gilson Concept			Gilson Neo				
			Single Channel					Single Channel			Multi Channel			Single Channel				
Abdos Code for Filter Pipette Tips	Abdos Code for Last Drop Filter Pipette Tips	Volume	0.2-2µl	1-10µl	2-20µl	20-100µl	20-200µl	200-1000µl	0.5-10µl	5-100µl	20-300µl	1-10µl	10-100µl	30-300µl	0.2-2µl	1-10µl	20-100µl	100-1000µl
P11001	P11050	0.2 - 10 µl		√					√			√				√		
P11002	P11051	0.2 - 10 µl		√					√			√				√		
P11003	P11052	2 - 20 µl				√	√			√			√				√	
P11026	P11066	2 - 50 µl				√	√			√			√				√	
P11004	P11053	2 - 100 µl				√	√			√			√				√	
P11005	P11054	2 - 200 µl				√	√			√			√				√	
P11006	P11055	2 - 300 µl				√	√			√	√		√	√			√	
P11007	P11056	100 - 1000 µl						√										√

Pipettor Brand/ Make			Vistalab Ovation QS					Vistalab Ovation M					Vistalab Ovation E				Eppendorf Reference				Eppendorf Electronic								
			Single Channel					Single Channel					Multi Channel				Single Channel												
Abdos Code for Filter Pipette Tips	Abdos Code for Last Drop Filter Pipette Tips	Volume	0.2-2µL	1-10µL	2-20µL	10-100µL	20-200µL	100-1000µL	0.2-2µL	1-10µL	2-20µL	10-100µL	20-200µL	100-1000µL	0.5-20µL	5-250µL	2-125µL	25-1250µL	0.1-10µL	'0.1-10µL 'ULTRA MICRO'	2-20µL	10-100µL	20-200µL	100-1000µL	0.5-10µL	5-100µL	2-200µL	50-1000µL	
P11001	P11050	0.2 - 10 µl	✓	✓	✓				✓	✓	✓			✓				✓		✓									
P11002	P11051	0.2 - 10 µl	✓	✓	✓				✓	✓	✓			✓				✓		✓									
P11003	P11052	2 - 20 µl				✓	✓				✓	✓			✓	✓						✓	✓			✓	✓		
P11026	P11066	2 - 50 µl				✓	✓				✓	✓			✓	✓						✓	✓			✓	✓		
P11004	P11053	2 - 100 µl				✓	✓				✓	✓			✓	✓						✓	✓			✓	✓		
P11005	P11054	2 - 200 µl				✓	✓				✓	✓			✓	✓						✓	✓			✓	✓		
P11006	P11055	2 - 300 µl				✓	✓				✓	✓			✓	✓						✓	✓			✓	✓		
P11007	P11056	100 - 1000 µl						✓					✓				✓							✓				✓	

Pipettor Brand/ Make			Eppendorf Research					Eppendorf					Thermo Finnpipette					Thermo Finnpipette		Matrix Impact/Impact2								
			Single Channel					Multi Channel					Single Channel					Multi Channel		s/mc								
Abdos Code for Filter Pipette Tips	Abdos Code for Last Drop Filter Pipette Tips	Volume	0.1-2.5µl	0.5-10µl	2-20µl	10-100µl	20-200µl	100-1000µl	0.5-10µl	5-100µl	50-100µl	0-300µl	30-300µl	50-1200µl	0.2-10µl	0.5-10µl	5-40µl	40-200µl	100-1000µl	200-1000µl	0.1-10µl	10-100µl	30-300µl	2-125µl	2-20µl	15-1250µl		
P11001	P11050	0.2 - 10 µl		✓	✓				✓					✓	✓						✓							
P11002	P11051	0.2 - 10 µl		✓	✓				✓					✓	✓						✓							
P11003	P11052	2 - 20 µl				✓	✓			✓	✓						✓					✓		✓				
P11026	P11066	2 - 50 µl				✓	✓			✓	✓						✓					✓		✓				
P11004	P11053	2 - 100 µl				✓	✓			✓	✓						✓					✓		✓				
P11005	P11054	2 - 200 µl				✓	✓			✓	✓						✓					✓		✓				
P11006	P11055	2 - 300 µl				✓	✓			✓	✓	✓	✓				✓					✓	✓	✓				
P11007	P11056	100 - 1000 µl						✓					✓*						✓	✓								

Pipettor Brand/ Make			Labnet Labpette					Labnet Labpette				Labnet BioPette E					Nichiryo /Oxford Benchmate										
			Single Channel					Multi Channel				Single Channel			Multi Channel		Single Channel & Multi Channel										
Abdos Code for Filter Pipette Tips	Abdos Code for Last Drop Filter Pipette Tips	Volume	0.1-2µl	0.5-10µl	2-20µl	10-100µl	20-200µl	100-1000µl	1-10µl	5-50µl	20-200µl	50-300µl	0.5-10µl	2-20µl	10-200µl	100-1000µl	0.5-10µl	2-20µl	10-200µl	100-1000µl	0.5-10µl	2-20µl	10-50µl	10-100µl	20-200µl	100-1000µl	
P11001	P11050	0.2 - 10 µl	✓	✓				✓				✓				✓					✓						
P11002	P11051	0.2 - 10 µl	✓	✓				✓				✓				✓					✓						
P11003	P11052	2 - 20 µl				✓	✓			✓				✓										✓	✓		
P11026	P11066	2 - 50 µl				✓	✓			✓				✓									✓	✓			
P11004	P11053	2 - 100 µl				✓	✓			✓				✓									✓	✓			
P11005	P11054	2 - 200 µl				✓	✓			✓				✓									✓	✓			
P11006	P11055	2 - 300 µl				✓	✓			✓	✓			✓									✓	✓			
P11007	P11056	100 - 1000 µl						✓						✓						✓							✓

Pipettor Brand/ Make			Socorex Calibra/Accura								Axygen Axypet																
			Single Channel				Multi Channel				Single Channel				Multi Channel												
Abdos Code for Filter Pipette Tips	Abdos Code for Last Drop Filter Pipette Tips	Volume	0.1-20µl	0.2-20µl	0.5-10µl	1-10µl	2-20µl	10-100µl	20-200µl	100-1000µl	1-10µl	5-50µl	10-100µl	20-200µl	50-200µl	0.1-2µl	0.5-10µl	2-20µl	5-50µl	10-100µl	20-200µl	50-250µl	100-1000µl	1-10µl	5-50µl	20-200µl	50-300µl
P11001	P11050	0.2 - 10 µl			√	√					√						√							√			
P11002	P11051	0.2 - 10 µl			√	√					√						√							√			
P11003	P11052	2 - 20 µl						√	√		√	√	√	√					√	√	√				√	√	√
P11026	P11066	2 - 50 µl						√	√		√	√	√	√					√	√	√				√	√	√
P11004	P11053	2 - 100 µl						√	√		√	√	√	√					√	√	√				√	√	√
P11005	P11054	2 - 200 µl						√	√		√	√	√	√					√	√	√				√	√	√
P11006	P11055	2 - 300 µl						√	√		√	√	√	√					√	√	√	√			√	√	√
P11007	P11056	100 - 1000 µl								√												√					

Pipettor Brand/ Make			BRAND		Eppendorf		Biohit Proline plus		Gilson		Thermo		Abdos	
Abdos Code for Filter Pipette Tips	Volume		1-5 ml	1-10 ml	1-5 ml	1-10 ml	1-5 ml	1-10 ml	1-5 ml	1-10 ml	1-5 ml	1-10 ml	1-5 ml	1-10 ml
P10144	1 - 5 ml				√		√		√				√	
P10145	1 - 10 ml			√		√		√		√				√

Pipettor Brand/ Make			Abdos premium	Abdos premium Plus	ABDOS D. Pipette (Electronic)	Brand Transferpette S	Gilson Pipetman	Eppendorf Research
Abdos Code for Filter Pipette Tips	Abdos Code for Last Drop Filter Pipette Tips	Volume	Single channel			Single channel		
			100-1000µl	100-1000µl	100-1000µl	100-1000µl	100-1000µl	100-1000µl
P11015	-	100 - 1000 µl XL	√	√	√	√	√	√
P11016	-	100 - 1250 µl	√	√	√	√	√	√

Pipettor Brand/ Make			THERMO (Electronic) Finn pipette	Thermo Finn pipette F3	Nichiryo NICHIPET EX	SOCOREX ACURA 825	BIOHIT Proline	ErgoOne Star Lab
Abdos Code for Filter Pipette Tips	Abdos Code for Last Drop Filter Pipette Tips	Volume	Multi channel	Single channel				
			100-1200µl	100-1000µl	100-1000µl	100-1000µl	100-1000µl	100-1000µl
P11015	P11057	100 - 1000 µl XL	√	√	√	√	√	√
P11016	P11058	100 - 1250 µl	√	√	√	√	√	√

Pipettor Brand/ Make		BRAND		Eppendorf		Biohit Proline plus		Gilson		Thermo		Abdos	
Abdos Code for Filter Pipette Tips	Volume	1-5 ml	1-10 ml	1-5 ml	1-10 ml	1-5 ml	1-10 ml	1-5 ml	1-10 ml	1-5 ml	1-10 ml	1-5 ml	1-10 ml
P10170	1 - 5 ml			√		√		√				√	
P10171	1 - 10 ml		√		√		√		√		√		√

Note: Variants of all tips are compatible with the pipette's brand mentioned in the above list.

Eg . Type : Tips - Bulk in Resealable Bags

Variants : Pipette Tips, (P10115, P10137, P10129, P10130, P10101, P10105, P10138, P10102, P10106, P10174, P10175)
 Last Drop™ Low Retention Pipette Tips (P10001, P10002, P10003, P10004, P10005, P10006)
 Filter Pipette Tips (P11001, P11002, P11003, P11004, P11005, P11006, P11007, P11015, P11016)
 Last Drop™ Low Retention Filter Pipette Tips (P11050, P11051, P11052, P11053, P11054, P11055, P11056)

Followed with other variants & Types: Tips - Individually wrapped, Racked, Racked Sterile & Refill pack.