


# ABDOS Screw Cap Tubes

# Features

- Abdos offers PCR clean screw cap tubes free from DNase/ RNase, pyrogen, endo-toxins & ATP.
- Gamma Irradiated, offers free of RNase/DNase, DNA, Pyrogen, Endotoxin, ATP and PCR Inhibitor.
- 100% virgin polypropylene, resin offers high clarity and ideal for biological & enzyme applications.
- O-ring design, thickness and material selection is optimized to provide protective features for maximum sealing ability, biologic interaction and temperature capability. The black colored O-ring is made from EPDM.
- Automated System Compatible. Whether automate is used for filling, capping or labeling, the product is designed to withstand high torque pressures associated with automation. Ideal for storage, transport of samples or reagents & for other biological applications.
- Available in various volumes ranging from 0.5 , 1.5 to 2 ml with high clarity , natural or in amber colour.
- Tube Color choices are an ultra-high clarity, natural color or amber option which provides ultra violet (UV) protection.
- Standard or Tethered cap, available in 4 colors ( Red, Yellow, Blue & Green).
- Knurl tubes feature the ability to lock tube into specialty racks for single-hand operation.
- Skirted or Self Standing tubes can stand upright on their own, no rack required.
- Conical tubes requires a rack to support them in the upright position.
- RCF upto 20,000 x g.
- Temperature Range -80°C to +95°C. Short-term storage in -80°C, -20°C, dry ice and wet ice. Short-term storage in vapor phase LN<sub>2</sub>. Suitable for a quick dip flash freeze in LN<sub>2</sub>.
- Water Bath compatible up to +95°C.

## Technical Drawings



# Ordering Details

Code	Product	Outer Diameter, mm	Height, mm	Style	Cap Color	Assembled	Sterile	No./ Case
			<b>0.5 ml</b>					
P10223	0.5 ml Screw Cap Tube	10.15	43.85	Self Standing	Natural	No	No	1000
P10224	0.5 ml Screw Cap Tube	10.15	43.85	Self Standing	Natural	No	Yes	1000
P10224A	0.5 ml Screw Cap Tube AMBER	10.15	43.85	Self Standing	Natural	No	Yes	1000
P10225 (R,Y,B,G)	0.5 ml Screw Cap Tube	10.15	43.85	Self Standing	Red, Yellow, Bue, Green	No	No	1000
P10226 (R,Y,B,G)	0.5 ml Screw Cap Tube	10.15	43.85	Self Standing	Red, Yellow, Bue, Green	No	Yes	1000
			<b>1.5 ml</b>					
P10227	1.5 ml Screw Cap Tube	10.15	43.95	Self Standing, Skirted base	Natural	No	No	1000
P10228	1.5 ml Screw Cap Tube	10.15	43.95	Self Standing, Skirted base	Natural	No	Yes	1000
P10228A	1.5 ml Screw Cap Tube AMBER	10.15	43.95	Self Standing, Skirted base	Natural	No	Yes	1000
P10229 (R,Y,B,G)	1.5 ml Screw Cap Tube	10.15	43.95	Self Standing, Skirted base	Red, Yellow, Bue, Green	No	No	1000
P10230 (R,Y,B,G)	1.5 ml Screw Cap Tube	10.15	43.95	Self Standing, Skirted base	Red, Yellow, Bue, Green	No	Yes	1000
			<b>1.5 ml</b>					
P10231	1.5 ml Screw Cap Tube	10.15	43.05	Conical	Natural	No	No	1000
P10232	1.5 ml Screw Cap Tube	10.15	43.05	Conical	Natural	No	Yes	1000
P10232A	1.5 ml Screw Cap Tube AMBER	10.15	43.05	Conical	Natural	No	Yes	1000
P10233 (R,Y,B,G)	1.5 ml Screw Cap Tube	10.15	43.05	Conical	Red, Yellow, Bue, Green	No	No	1000
P10234 (R,Y,B,G)	1.5 ml Screw Cap Tube	10.15	43.05	Conical	Red, Yellow, Bue, Green	No	Yes	1000
			<b>2.0 ml</b>					
P10235	2.0 ml Screw Cap Tube	10.15	45.05	Self Standing, Skirted base	Natural	No	No	1000
P10236	2.0 ml Screw Cap Tube	10.15	45.05	Self Standing, Skirted base	Natural	No	Yes	1000
P10236A	2.0 ml Screw Cap Tube AMBER	10.15	45.05	Self Standing, Skirted base	Natural	No	Yes	1000
P10237 (R,Y,B,G)	2.0 ml Screw Cap Tube	10.15	45.05	Self Standing, Skirted base	Red, Yellow, Bue, Green	No	No	1000
P10238 (R,Y,B,G)	2.0 ml Screw Cap Tube	10.15	45.05	Self Standing, Skirted base	Red, Yellow, Bue, Green	No	Yes	1000
			<b>2.0 ml</b>					
P10239	2.0 ml Screw Cap Tube	10.15	44.30	Conical	Natural	No	No	1000
P10240	2.0 ml Screw Cap Tube	10.15	44.30	Conical	Natural	No	Yes	1000
P10240A	2.0 ml Screw Cap Tube AMBER	10.15	44.30	Conical	Natural	No	Yes	1000
P10241 (R,Y,B,G)	2.0 ml Screw Cap Tube	10.15	44.30	Conical	Red, Yellow, Bue, Green	No	No	1000
P10242 (R,Y,B,G)	2.0 ml Screw Cap Tube	10.15	44.30	Conical	Red, Yellow, Bue, Green	No	Yes	1000