

Technical Data Sheet

Doc No Ab/TDS/ELECTRO-122

Rev No: 02 Date 02/2022-23

TECHNICAL DATA SHEET

Abdos Labtech Pvt. Ltd.

ABDOS Vertical Midi Gel System

Valid for Catalog No. E12276

1. Description Specifications:

1.1 Description: Vertical Midi Gel System Includes – Midi Vertical Assembly, 2 sets of Glass Plates (Rectangular & Notched), 13 wells 2 Combs each of 1 mm & 1.5mm thickness, 4 nos. spacers of 1mm & 1.5 mm Thickness, Gel Casting Stand - 1 No, Supporting Screws-2 No's, Safety Lid - 1 No, Connecting Cord - 1 Pair

1.1.1 The system is designed to run One 16x14cm midi hand cast PAGE & SDS PAGE gels.

1.1.2 Ventilation holes in the system prevent the build-up of condensation during electrophoresis.

1.1.3 The system comes with standard casting unit, combs & spacers of 1 & 1.5mm thickness & a set of glass plates.

1.1.4 Low buffer consumption design having small footprint area, occupies minimal lab bench space.

1.1.5 Unique “T” spacers for leak proof casting.

1.2 Design: Gel Electrophoresis system

1.3 Performance Specifications:

Specifications	E12276
Dimensions (LxWxH)	25 x 20 x 17cm
Gel Dimension (WxH)	16x14cm (1 gels)
Buffer volume	500ml
Electrode wires	Brass, Platinum
Standard Combs-Teflon	13 Well-1.0mm and 1.5mm (2 Nos each)
Standard Combs-Teflon	1.0mm and 1.5mm (2Nos each)
Max. Sample Capacity	40
Time Consumed	4 to 5 hrs
Recommended Power Supply	E12377 or E12379



Abdos Labtech Private Limited
(A Division of Abdos Group)
(An ISO 9001,13485,14001,18385,GMP & CE certified company)
Shiv Ganga Industrial Estate
Plot No-1, Lakheswari, Bhagwanpur,
Roorkee-247667, (UK) INDIA

Registered office
13A Government
Place East, 2nd Floor Kolkata
West Bengal – 700069
CIN: U85100WB2008PTC123903
Corporate Office
319, 3rd Floor, DLF Tower A, DDA Dist. Centre
Jasola, New Delhi - 110025

Ph. No: +91- (0) 1332235277
Email: labtech@abdosindia.com
Web: www.abdoslivesciences.com
GST: 05AAGCA8605A1Z7

Technical Data Sheet

1.4 Material: Acrylic

1.5 Color: Natural

1.6 CE Certified

1.7 Quality control:

Product are manufactured, inspected and found to be in compliance with product & quality specification requirements as documented in ISO 9001:2015 & 13485:2016 Quality Management System. Testing of various characteristics in accordance with the valid specification, inspection reviewed & signed for product release

1.8 Spare Accessories

Code	Product	No. Per Case
E12277	Glass Plates(Notched & Rectangular)	1set
E12278	Regular Midi Gel Casting Unit(For 2 Gels)	1
E12279	Gel Developing Tray with Gel Holder 16 x 14 cm	1
E12280	Gradient Maker Total Capacity 25+25ml	1
E12281	Teflon Comb 13 well 0.5 mm	1
E12282	Teflon Comb 13 well 1.0 mm	1
E12283	Teflon Comb 10 well 1.5 mm	1
E12284	Teflon Comb 20 well 0.5 mm	1
E12295	Teflon Comb 20 well 1.0 mm	1
E12296	Teflon Comb 10 well 1.5 mm	1
E12285	Teflon Spacers 0.5 mm	6
E12286	Teflon Spacers 1.0 mm	6
E12287	Teflon Spacers 1.5 mm	6
E12291	ACRYLIC Clamp(1 pcs) and Screws(2 pcs)	1set
E12292	Levelling Screws	3
E12293	Glass Plate Stand	1
E12244	Connecting Cords	2

Warranty: 12 months from the date of Installation or 15 months from the date of Invoicing



Abdos Labtech Private Limited
(A Division of Abdos Group)
(An ISO 9001,13485,14001,18385,GMP & CE certified company)
Shiv Ganga Industrial Estate
Plot No-1, Lakheswari, Bhagwanpur,
Roorkee-247667, (UK) INDIA

Registered office
13A Government
Place East, 2nd Floor Kolkata
West Bengal – 700069
CIN: U85100WB2008PTC123903
Corporate Office
319, 3rd Floor, DLF Tower A, DDA Dist. Centre
Jasola, New Delhi - 110025

Ph. No: +91- (0) 1332235277
Email:labtech@abdosindia.com
Web: www.abdoslivesciences.com
GST: 05AAGCA8605A1Z7

Technical Data Sheet

Instruction for Use

Storage temp. Max 30°C
Before using the product check crackness etc.
Check all the information printed on pouches or boxes.
Always wear gloves while handling specimen.

WARNING

DO NOT store strong oxidizing agents.
DO NOT place any plastic labware in a flame.
DO NOT mix any chemicals that may result in a thermic reaction, which can cause product failure
Consult local fire codes prior to storage of flammable liquids.



Abdos Labtech Pvt. Ltd.

Plot No. 1, Lakeshwari Bhagwanpur Roorkee 247667 (UK) India
Email : labtech@abdosindia.com, www.abdoslivesciences.com
Ph: +{91} 11 41081897



Quality Assurance Incharge
Abdos Labtech Pvt. Ltd.



Abdos Labtech Private Limited
(A Division of Abdos Group)
(An ISO 9001,13485,14001,18385,GMP & CE certified company)
Shiv Ganga Industrial Estate
Plot No-1, Lakshewari, Bhagwanpur,
Roorkee-247667, (UK) INDIA

Registered office
13A Government
Place East, 2nd Floor Kolkata
West Bengal – 700069
CIN: U85100WB2008PTC123903
Corporate Office
319, 3rd Floor, DLF Tower A, DDA Dist. Centre
Jasola, New Delhi - 110025

Ph. No: +91- (0) 1332235277
Email: labtech@abdosindia.com
Web: www.abdoslivesciences.com
GST: 05AAGCA8605A1Z7

Technical Data Sheet



Abdos Labtech Private Limited
(A Division of Abdos Group)
(An ISO 9001,13485,14001,18385,GMP & CE certified company)
Shiv Ganga Industrial Estate
Plot No-1, Lakheswari, Bhagwanpur,
Roorkee-247667, (UK) INDIA

Registered office
13A Government
Place East, 2nd Floor Kolkata
West Bengal – 700069
CIN: U85100WB2008PTC123903
Corporate Office
319, 3rd Floor, DLF Tower A, DDA Dist. Centre
Jasola, New Delhi - 110025

Ph. No: +91- (0) 1332235277
Email:labtech@abdosindia.com
Web: www.abdoslifesciences.com
GST: 05AAGCA8605A1Z7