

Technical Data Sheet

Doc No–Ab/TDS/UTL-208

Rev No: 07 Date 02/2022-23

TECHNICAL DATA SHEET

Abdos Labtech Pvt. Ltd.

Abdos Brand – Derma Gaurd Blue Nitrile Gloves – 9.5 inches length

Valid for Catalog No. U20341, U20342, U20343, U20447, U20448

1. Description/Specification

1.1 Description: Abdos Brand – Derma Gaurd Blue Nitrile Gloves – 9.5 inches length

1.1.1 These gloves are ideal product to protect you exposure to bacteria, viruses and chemicals. It also enhances your ability to offer maximum security against cross contamination.

1.1.2 Its unique blue color provides quick, visual differentiation from natural rubber latex and other synthetic gloves.

1.1.3 Ambidextrous, textured fingertips and beaded cuff enhances wet & dry grip.

1.1.4 Custom design enhances comfort & fit.

1.1.5 These gloves are powder free, minimizing the potential for powder related complications.

1.1.6 Conforms to ASTM D6319 and EN 455 Standards, EU standard directive 93/42/EEC Class I

1.1.7 Derma Gaurd Blue Nitrile Gloves Powder free- sample meets the overall migration requirement under Commission Regulation (EU) No 10/2011- Plastic materials and articles shall not transfer their constituents to foodstuff in quantities exceeding 10mg of total constituents released per dm² of food contact surface (mg/dm²) (overall migration limit)

1.1.8 Palm Thickness :0.08mm (3.15mil) & Finger Thickness : 0.11mm (4.33mil)

1.1.9 Average Length : 240mm

1.2 Design: Hand gloves design

1.3 Performance Specifications:

Code	Size	No. Per Case
U20341	Small	100
U20342	Medium	100
U20343	Large	100
U20447	X Small	100
U20448	X Large	100



Abdos Labtech Private Limited
(A Division of Abdos Group)
(An ISO 9001,13485,14001,18385,GMP & CE certified company)
Shiv Ganga Industrial Estate
Plot No-1, Lakhewari, Bhagwanpur,
Roorkee-247667, (UK) INDIA

Registered office
13A Government
Place East, 2nd Floor Kolkata
West Bengal – 700069
CIN: U85100WB2008PTC123903
Corporate Office
319, 3rd Floor, DLF Tower A, DDA Dist. Centre
Jasola, New Delhi - 110025

Ph. No: +91- (0) 1332235277
Email:labtech@abdosindia.com
Web: www.abdoslivesciences.com
GST: 05AAGCA8605A1Z7

Technical Data Sheet

1.4 Material: Nitrile gloves

1.5 Color: Blue

1.6 Quality control:

Product are manufactured, inspected and found to be in compliance with product & quality specification requirements as documented in ISO 9001:2015 & 13485:2016 Quality Management System. Testing of various characteristics in accordance with the valid specification, inspection reviewed & signed for product release

1.7 Other information

1.7.1 These gloves are used for working on lab, machine and any other work.

1.7.2 These gloves are ideal product to protect you any of microbial contamination.

1.7.3 These Flexible gloves are fit in hand & easy to use on working.

1.7.4 Powder free Chlorinated gloves

1.7.5 Suitable for food handling

1.7.6 **Chemical resistance:** See website (www.abdoslifesciences.com)

1.7.7 Manufactured in Aseptic Environment

Instruction for Use

Product can be used for food contact also.

Storage temp. Max 30°C

Before using the product check crackness, punctured etc.

Check all the information printed on pouches or boxes.

Always wear gloves while handling specimen.

WARNING

DO NOT store strong oxidizing agents.

DO NOT place any plastic labware in a flame.

DO NOT mix any chemicals that may result in a thermic reaction, which can cause product failure

Consult local fire codes prior to storage of flammable liquids.




Abdos Labtech Private Limited
(A Division of Abdos Group)
(An ISO 9001,13485,14001,18385,GMP & CE certified company)
Shiv Ganga Industrial Estate
Plot No-1, Lakheswari, Bhagwanpur,
Roorkee-247667, (UK) INDIA

Registered office
13A Government
Place East, 2nd Floor Kolkata
West Bengal – 700069
CIN: U85100WB2008PTC123903
Corporate Office
319, 3rd Floor, DLF Tower A, DDA Dist. Centre
Jasola, New Delhi - 110025

Ph. No: +91- (0) 1332235277
Email:labtech@abdosindia.com
Web: www.abdoslifesciences.com
GST: 05AAGCA8605A1Z7

Technical Data Sheet



Abdos Labtech Pvt. Ltd.

Plot No. 1 Lakshwari Bhagwanpur Roorkee 247667 (UK) India
Email : labtech@abdosindia.com, www.abdoslifesciences.com
Ph: +(91) 11 41081897



Abdos Labtech Private Limited
(A Division of Abdos Group)
(An ISO 9001,13485,14001,18385,GMP & CE certified company)
Shiv Ganga Industrial Estate
Plot No-1, Lakshwari, Bhagwanpur,
Roorkee-247667, (UK) INDIA

Registered office
13A Government
Place East, 2nd Floor Kolkata
West Bengal – 700069
CIN: U85100WB2008PTC123903
Corporate Office
319, 3rd Floor, DLF Tower A, DDA Dist. Centre
Jasola, New Delhi - 110025

Ph. No: +91- (0) 1332235277
Email:labtech@abdosindia.com
Web: www.abdoslifesciences.com
GST: 05AAGCA8605A1Z7