

Technical Data Sheet

Doc No–Ab/TDS/UTL-209

Rev No: 06 Date 02/2022-23

TECHNICAL DATA SHEET

Abdos Labtech Pvt. Ltd.

Abdos Brand – Derma Gaurd White Latex Gloves – 9.5 inches length

Valid for Catalog No. U20344, U20345, U20346, U20449, U20450

1. Description/Specification

1.1 Description: Abdos Brand – Derma Guard White Latex Gloves – 9.5 inches length

1.1.1 These gloves are ideal product to protect you exposure to bacteria, viruses and chemicals. It also enhances your ability to offer maximum security against cross contamination.

1.1.2 Ambidextrous, textured fingertips and beaded cuff enhances wet & dry grip.

1.1.3 Softness provides superior comfort & natural fit.

1.1.4 These gloves are powder free, minimizing the potential for powder related complications, irritation and dermatitis.

1.1.5 Lower level of extractable proteins and chemical residue, hence lower incidences of allergic reactions.

1.1.6 Conforms to ASTM D3578 and EN 455 Standards, EU standard directive 93/42/EEC Class I

1.1.7 Palm Thickness :0.10mm (3.93mil) & Finger Thickness : 0.10mm (3.93mil)

1.1.8 Average Length : 240mm

1.2 Design: Hand gloves design.

1.3 Performance Specifications:

Code	Size	No. Per Case
U20344	Small	100
U20345	Medium	100
U20346	Large	100
U20449	X Small	100
U20450	X Large	100

1.4 Material: Latex gloves



Abdos Labtech Private Limited
(A Division of Abdos Group)
(An ISO 9001,13485,14001,18385,GMP & CE certified company)
Shiv Ganga Industrial Estate
Plot No-1, Lakheswari, Bhagwanpur,
Roorkee-247667, (UK) INDIA

Registered office
13A Government
Place East, 2nd Floor Kolkata
West Bengal – 700069
CIN: U85100WB2008PTC123903
Corporate Office
319, 3rd Floor, DLF Tower A, DDA Dist. Centre
Jasola, New Delhi - 110025

Ph. No: +91- (0) 1332235277
Email:labtech@abdosindia.com
Web: www.abdoslivesciences.com
GST: 05AAGCA8605A1Z7

Technical Data Sheet

1.5 Color: White

1.6 Quality control:

Product are manufactured, inspected and found to be in compliance with product & quality specification requirements as documented in ISO 9001:2015 & 13485:2016 Quality Management System. Testing of various characteristics in accordance with the valid specification, inspection reviewed & signed for product release

1.7 Other information

1.7.1 These gloves are used for working on lab, machine and any other work.

1.7.2 These gloves are ideal product to protect you any of microbial contamination.

1.7.3 These Flexible gloves are fit in hand & easy to use on working.

1.7.4 Powder free Chlorinated gloves

1.7.5 **Chemical resistance:** See website (www.abdoslifesciences.com)

1.7.6 Manufactured in Aseptic Environment

Instruction for Use

Product can be used for food contact also.

Storage temp. Max 30°C

Before using the product check crackness, punctured etc.

Check all the information printed on pouches or boxes.

Always wear gloves while handling specimen.

WARNING

DO NOT store strong oxidizing agents.

DO NOT place any plastic labware in a flame.

DO NOT mix any chemicals that may result in a thermic reaction, which can cause product failure

Consult local fire codes prior to storage of flammable liquids.



Abdos Labtech Pvt. Ltd.

Plot No. 1 Lakshwari Bhagwanpur Roorkee 247667 (UK) India
Email : labtech@abdosindia.com, www.abdoslifesciences.com
Ph: +(91) 11 41081897




Abdos Labtech Private Limited
(A Division of Abdos Group)
(An ISO 9001,13485,14001,18385,GMP & CE certified company)
Shiv Ganga Industrial Estate
Plot No-1, Lakshwari, Bhagwanpur,
Roorkee-247667, (UK) INDIA

Registered office
13A Government
Place East, 2nd Floor Kolkata
West Bengal – 700069
CIN: U85100WB2008PTC123903

Corporate Office
319, 3rd Floor, DLF Tower A, DDA Dist. Centre
Jasola, New Delhi - 110025

Ph. No: +91- (0) 1332235277
Email:labtech@abdosindia.com
Web: www.abdoslifesciences.com
GST: 05AAGCA8605A1Z7

Technical Data Sheet


Quality Assurance Incharge
Abdos Labtech Pvt. Ltd.



Abdos Labtech Private Limited
(A Division of Abdos Group)
(An ISO 9001,13485,14001,18385,GMP & CE certified company)
Shiv Ganga Industrial Estate
Plot No-1, Lakheswari, Bhagwanpur,
Roorkee-247667, (UK) INDIA

Registered office
13A Government
Place East, 2nd Floor Kolkata
West Bengal – 700069
CIN: U85100WB2008PTC123903

Corporate Office
319, 3rd Floor, DLF Tower A, DDA Dist. Centre
Jasola, New Delhi - 110025

Ph. No: +91- (0) 1332235277
Email:labtech@abdosindia.com
Web: www.abdoslivesciences.com
GST: 05AAGCA8605A1Z7