



ABDOS
LIFE SCIENCES
Advancing Science with Quality

Deep Well Plates

**IDEAL FOR SAMPLE
STORAGE, SAMPLE
TREATMENT &
AUTOMATED OPERATION**



Features

- Compatible with many applications such as HTS assays, sample storage, nucleic acid extraction, sample preparation, elution, and liquid handling stations
- Compatible with common Automated Nucleic Acid Extractors and Liquid Handling systems
- Multichannel pipette and automatic equipment compatible
- Standard SBS footprint
- 96 well plates available in square and round wells with conical and round bottoms
- Alphanumeric marking for easy identification
- Low extractables and non-leachable

Abdos Deep Well Plates

All Abdos Deep Well Plates conform to SBS International Standard and are compatible with most of the pipetting workstations and extraction systems available. Manufactured in a state-of-the-art precision mould facility from USP class VI medical grade virgin polypropylene conforming USFDA 21CFR, free from natural rubber and heavy metals. The range includes high and low profile Deep Well Plates in different formats - Comb, Silicone Mats, and Reservoirs.

Compatibility Chart

[illegible]

Quality Compliance




Ordering Information

Code	Description	Nos/Case
P10802	2.2 ml, 96 Square Well, Round Conical Bottom ("V") Deep Well Plate	50
P10808	96 Magnetic Tip Comb Suitable for P10802	100
P10809	0.5 ml, 96 Square Well Conical ("V") Bottom (200 ul Elution Plate) Deep Well Plate	50
P10810	2.0 ml, 96 Round Well, Round ("U") Bottom Deep Well Plate	50
P10811	1.0 ml, 96 Round Well, Round ("U") Bottom Deep Well Plate	50
P10803	1.2 ml, 96 Square Well, Conical ("V") Bottom, Deep Well Plate	100
P10801	12 Channel Troughs Reservoir, Each Well 15 ml	25
P10807	Adhesive Sealing Film, Polyester	500
P10804	Silicone Mat for 96 Square Well Plate, Pierceable	50
P10805	Silicone Mat for 2.0 ml 96 Round Well Plate, Pierceable	50
P10806	Silicone Mat for 96 Round Well Plate 1.0 ml, Pierceable	50


Other Items of Interest from PCR & Deep Well Plate Range



AmpQlear
PCR Products



PCR Plate Accessories



PCR Plates

What Abdos Life Sciences stands for?



Precision



Reliability



Responsiveness



Innovation



Sustainability



Relationships

Get in touch with us:

☎ 011-41081897

✉ labtech@abdosindia.com

🌐 www.abdoslifesciences.com