

# Abdos PET PETG Media Bottles

## Ergonomic SQUARE BOTTLES

Narrow mouth to provide better control when pouring

Square shape, easy-to-hold, and saves shelf counter space compared with round bottles of the same volume

Lightweight impact-resistant material that can withstand -70°C ~ +50°C



### Features:

Available in two polymer variants: Polyethylene Terephthalate Glycol co-polymer (PETG) and Polyethylene Terephthalate (PET)

Moulded embossed graduations for easy visual observation of the solution volume

Highly transparent and optically clear with smooth inner & outer walls

Abdos PET & PETG Media Bottles have excellent gas barrier properties, making them best-suited for storing and sampling active pharmaceutical ingredients and bulk intermediates.

These media bottles are also ideal for preparing and storing buffers, serum, culture fluids, or long-term storage of pH-sensitive liquids.

## Wide Volume Options Available



**60 ml**  
(PETG)



**125 ml**  
(PET/PETG)



**250 ml**  
(PET/PETG)




**500 ml**  
(PET/PETG)



**1000 ml**  
(PET/PETG)

## Certifications and Compliances

**SAL 10<sup>-6</sup>**



ATP  
FREE



CYTOXIN

FREE



ENDOTOXIN

FREE



DNase  
RNase  
PYROGEN

FREE



ISO  
13485



ISO  
9001



TSE &  
BSE

FREE

# Ordering Information


PETG Media Bottles		
Code	Description	Nos/Case
P21170	60ml, Non-Sterile	192
P21171	125ml, Non-Sterile	96
P21172	250ml, Non-Sterile	60
P21173	500ml, Non-Sterile	48
P21174	1000ml, Non-Sterile	24
P21180	60ml, Sterile	192
P21181	125ml, Sterile	96
P21182	250ml, Sterile	60
P21183	500ml, Sterile	48
P21184	1000ml, Sterile	24

PET Media Bottles		
Code	Description	Nos/Case
P21176	125ml, Non-Sterile	96
P21177	250ml, Non-Sterile	224
P21178	500ml, Non-Sterile	120
P21179	1000ml, Non-Sterile	40
P21186	125ml, Sterile	96
P21187	250ml, Sterile	224
P21188	500ml, Sterile	120
P21189	1000ml, Sterile	40


# Other Items of Interest in Cell Culture Range



BIOFILL Bottle  
Top Filtration System



Cell Culture Plates



Cell Culture Dishes



Cell Culture Flasks

## What Abdos Life Sciences stands for?



Precision



Reliability



Responsiveness



Innovation



Sustainability



Relationships

Get in touch with us:

☎ 011-41081897

✉ [labtech@abdosindia.com](mailto:labtech@abdosindia.com)

🌐 [www.abdoslifesciences.com](http://www.abdoslifesciences.com)