



CENTRIFUGE TUBES

**BEST-IN-CLASS QUALITY
BREAKAGE-RESISTANT
DESIGNS**

**Revolutionise Your
Research with ABDOS
Centrifuge Tubes**

ABDOS Centrifuge Tubes are ergonomically designed with easy-to-grip caps, manufactured from high-quality Polypropylene ensuring durability and resistance to breakage



Features:

Sterile tubes have been validated according to **EN ISO 11137-2: SAL 10-6**

All new **"MAX FREEZE"** indicator for maximum liquid volume for samples that will be frozen

Highly Transparent - Allowing researchers to easily monitor and observe the contents inside the tube

Downward Scale Graduation for precise measurements



Even Wall Thickness offers maximum centrifugal forces, up to 18,000 RCF.

Precise Lid Sealing & Perfect Grip on the screw prevents any contamination or loss of samples during centrifugation

Accurate Graduation Marks enable precise measurement of sample volumes & separation of specific layers during centrifugation.

Suitable for a temperature range between **-80°C (deep freeze)** to **+121 °C (autoclavable)**

Large White Frosted Portion for easy writing **Hydrodynamic Conical Bottom** for perfect pelleting of sediments

Resins are Certified from USP Class VI
Medical grade 100 % virgin Polypropylene, conforming to US FDA 21 CFR assuring the highest quality in its class

Compatibility with Various Centrifuge Models and Rotors - Allowing flexibility in laboratory experimental setups

THE CAPS


Flat Thread Confirms Vacuum Tight Closure – leak-proof and prevents override of cap



Unique Plug Seal – ensuring perfect sealing & make it leak-proof










Advantages of Abdos Centrifuge Tubes with Even Wall Thickness

Even wall thickness easily distributes the G Force within the tube and thus **Abdos Centrifuge Tubes** can take greater RCF



Abdos Centrifuge Tubes in the Leading Journals



| | | |
|---|--|---|
|  |  |  |
|  |  |  |
|  |  |  |

Quality Compliances



Ordering Information

| Code | Product Description | Pack Size |
|----------------|---|-----------------------------------|
| P10430 | Centrifuge Tube, Conical, PP, 15 mL, Bulk, DNase/RNase free - 500 un. – ABDOS | 500 |
| P10401 | Centrifuge Tube, Conical, PP, 15 mL, DNase/RNase free - 20 x 25 un. – ABDOS | 500 |
| P10525 | Last Drop Low Retention Centrifuge Tube, Conical, PP, 15 mL, DNase/RNase free - 20 x 25 un. - ABDOS | 500 |
| P10441 | Centrifuge Tube, Conical, PP, 15 mL, Red Cap, DNase/RNase free - 20 x 25 un. - ABDOS | 500 |
| P10402 | Centrifuge Tube, Conical, PP, 15 mL, Sterile, DNase/RNase free - 20 x 25 un. – ABDOS | 500 |
| P10443 | Centrifuge Tube, Conical, PP, 50 mL, Sterile, Red Cap, DNase/RNase free - 20 x 25 un. - ABDOS | 500 |
| P10404 | Centrifuge Tube, Conical, PP, 50 mL, Sterile, DNase/RNase free - 20 x 25 un. – ABDOS | 500 |
| P10442 | Centrifuge Tube, Conical, PP, 15 mL, Sterile, Red Cap, DNase/RNase free - 20 x 25 un. - ABDOS | 500 |
| P10431 | Centrifuge Tube, Conical, PP, 50 mL, Bulk, DNase/RNase free - 500 un. – ABDOS | 500 |
| P10403 | Centrifuge Tube, Conical, PP, 50 mL, DNase/RNase free - 20 x 25 un. – ABDOS | 500 |
| P10444 | Centrifuge Tube, Conical, PP, 50 mL, Sterile, Red Cap, DNase/RNase free - 20 x 25 un. – ABDOS | 500 |
| P10405 | Centrifuge Tube, Conical, PP, 15 mL, Amber, DNase/RNase free - 20 x 25 un. – ABDOS | 500 |
| P10406 | Centrifuge Tube, Conical, PP, 15 mL, Sterile, Amber, DNase/RNase free - 20 x 25 un. – ABDOS | 500 |
| P10407 | Centrifuge Tube, Conical, PP, 50 mL, Amber, DNase/RNase free - 20 x 25 un. – ABDOS | 500 |
| P10408 | Centrifuge Tube, Conical, PP, 50 mL, Sterile, Amber, DNase/RNase free - 20 x 25 un. – ABDOS | 500 |
| P10405T | Centrifuge Tube, Conical, PP, 15 mL, Transparent Amber, DNase/RNase free - 20 x 25 un. - ABDOS | 500 |
| P10406T | Centrifuge Tube, Conical, PP, 15 mL, Sterile, Transparent Amber, DNase/RNase free - 20 x 25 un. - ABDOS | 500 |
| P10407T | Centrifuge Tube, Conical, PP, 50 mL, Transparent Amber, DNase/RNase free - 20 x 25 un. - ABDOS | 500 |
| P10408T | Centrifuge Tube, Conical, PP, 50 mL, Sterile, Transparent Amber, DNase/RNase free - 20 x 25 un. - ABDOS | 500 |
| P10419 | Centrifuge Tube, Conical, PS/HDPE, 15 mL, Sterile, Cell Culture, DNase/RNase free - 20 x 25 un. - ABDOS | 500 |
| P10421 | Centrifuge Tube, Conical, Racked, PS/HDPE, 15 mL, Sterile, Cell Culture, DNase/RNase free - 20 x 25 un. - ABDOS | 500 (25 PCS X 20 RACKS) |
| P10432 | Centrifuge Tube, Conical, PP, 50 mL, Self-Standing, Bulk, DNase/RNase free - 500 un. - ABDOS | 500 |
| P10423 | Centrifuge Tube, Conical, PP, 50 mL, Self-Standing, DNase/RNase free - 20 x 25 un. – ABDOS | 500 |
| P10449 | Centrifuge Tube, Conical, PP, 50 mL, Self-Standing, Red Cap, DNase/RNase free - 20 x 25 un. - ABDOS | 500 |
| P10424 | Centrifuge Tube, Conical, PP, 50 mL, Self-Standing, Sterile, DNase/RNase free - 20 x 25 un. - ABDOS | 500 |
| P10450 | Centrifuge Tube, Conical, PP, 50 mL, Self-Standing, Sterile, Red Cap, DNase/RNase free - 20 x 25 un. - ABDOS | 500 |

What Abdos Life Sciences stands for?



Precision



Reliability



Responsiveness



Innovation



Sustainability



Relationships

Get in touch with us:

011-41081897



labtech@abdosindia.com



www.abdoslifesciences.com