

Centrifuge Ware



- ◆ Abdos centrifuge tubes are manufactured from ultra-pure USP class VI certified virgin polypropylene to ensure the highest quality & confirming to US FDA 21 CFR. Free from heavy metals, natural rubber & inhibitory plasticizers.
- ◆ Robust polymer selection process ensures Abdos tube do not leech & migrate unwanted substances into your sample.
- ◆ Unique printing technology makes the surface chemical resistant & scratch proof.
- ◆ Perfect grip on screw cap provides secure control on the cap in wet conditions.
- ◆ Even wall thickness to withstand maximum G-force, 15000 RCF. Suitable for temperatures of -80° to +121° C.
- ◆ Sterile tubes have been validated according to ENISO 11137-2; SAL 10⁻⁶
- ◆ The larger white frosted portion for easy writing.
- ◆ Highly transparent, 5% less haze percentage than a competitor.
- ◆ Hydrodynamic conical bottom for perfect pelleting of sediments.

THE CAP

- ◆ Specially design flat thread confirm perfect vacuum tight closure make leak proof & prevent from override of cap.
- ◆ Unique plug seal ensure perfect sealing & make it leak proof.



Creating For Life Sciences

www.abdoslifesciences.com

Experience the Consistency & Purity with Abdos Centrifuge Tubes

Code	Product Description	Pack Size	No. Per Case
P10430	Abdos Centrifuge Tubes Conical Bottom, PP /HOPE, 15 ml Bulk, Blue Cap	25 Pcs / Bag	500
P10401	Abdos Centrifuge Tubes Conical Bottom, PP /HOPE, 15 ml, Blue Cap	25 Pcs / Bag	500
P10441	Abdos Centrifuge Tubes Conical Bottom, PP /HOPE, 15 ml, Red Cap	25 Pcs / Bag	500
P10402	Abdos Centrifuge Tubes Conical Bottom, PP /HOPE, 15 ml, Sterile, Blue Cap	25 Pcs / Bag	500
P10442	Abdos Centrifuge Tubes Conical Bottom, PP /HOPE, 15 ml, Sterile, Red Cap	25 Pcs / Bag	500
P10431	Abdos Centrifuge Tubes Conical Bottom, PP /HOPE, 50 ml, Bulk, Blue cap	25 Pcs / Bag	500
P10403	Abdos Centrifuge Tubes Conical Bottom, PP /HOPE, 50 ml, Blue Cap	25 Pcs / Bag	500
P10443	Abdos Centrifuge Tubes Conical Bottom, PP /HOPE, 50 ml, Red Cap	25 Pcs / Bag	500
P10404	Abdos Centrifuge Tubes Conical Bottom, PP /HOPE, 50 ml, Sterile, Blue Cap	25 Pcs / Bag	500
P10444	Abdos Centrifuge Tubes Conical Bottom, PP /HOPE, 50 ml, Sterile, Red Cap	25 Pcs / Bag	500
P10405	Abdos Centrifuge Tubes Conical Bottom, PP/ Amber, 15 ml, Amber, Blue Cap	25 Pcs / Bag	500
P10406	Abdos Centrifuge Tubes Conical Bottom, PP/ Amber, 15 ML, Sterile Amber, Blue Cap	25 Pcs / Bag	500
P10407	Abdos Centrifuge Tubes Conical Bottom, PP/ Amber, 50 ml, Amber, Blue Cap	25 Pcs / Bag	500
P10408	Abdos Centrifuge Tubes Conical Bottom, PP/ Amber, 50 ml, Sterile Amber, Blue Cap	25 Pcs / Bag	500
P10432	Abdos Self Standing Centrifuge Tubes, PP/HOPE, 50 ml, Bulk, Blue Cap	25 Pcs / Bag	500
P10423	Abdos Self Standing Centrifuge Tubes, PP/HOPE, 50 ml, Blue Cap	25 Pcs / Bag	500
P10449	Abdos Self Standing Centrifuge Tubes, PP/HOPE, 50 ml, Red Cap	25 Pcs / Bag	500
P10424	Abdos Self Standing Centrifuge Tubes, PP/HOPE, 50 ml, Sterile, Blue Cap	25 Pcs / Bag	500
P10450	Abdos Self Standing Centrifuge Tubes, PP/HOPE, 50 ml, Sterile, Red Cap	25 Pcs / Bag	500
P10419	Abdos Centrifuge Tubes Conical Bottom, PS /HOPE, 15 ml, Sterile, Blue Cap	25 Pcs / Bag	500
P10421	Abdos Centrifuge Tubes Conical Bottom, PS /HOPE, 15 ml, Racked Sterile, Blue Cap	25 Pcs / Bag	750
P10425	Abdos Individually Wrapped Centrifuge Tubes, PP /HOPE, 15 ml, Sterile, Blue Cap	1 Ps / Bag	500
P10451	Abdos Individually Wrapped Centrifuge Tubes, PP /HOPE, 15 ml, Sterile, Red Cap	1 Ps / Bag	500
P10426	Abdos Individually Wrapped Centrifuge Tubes, PP /HOPE, 50 ml, Sterile, Blue Cap	1 Ps / Bag	500
P10452	Abdos Individually Wrapped Centrifuge Tubes, PP /HOPE, 50 ml, Sterile, Red Cap	1 Ps / Bag	500
P10427	Abdos Individually Wrapped Centrifuge Tubes, PP /HOPE, 50 ml, Sterile, Self Standing, Blue Cap	1 Ps / Bag	450
P10453	Abdos Individually Wrapped Centrifuge Tubes, PP /HOPE, 50 ml, Sterile, Self Standing, Red Cap	1 Ps / Bag	450



Choose from Abdos range of Centrifuge Boxes & Racks. For more details please visit us at www.abdoslifesciences.com/centrifuge-ware.html



Centrifuge Tube Box



Click Together Conical Tube Racks



Foldable Space Saving Rack

Abdos centrifuge tubes are available in environment friendly paper rack. For more information please log on to www.abdoslifesciences.com/racked-centrifuge-tube-conical-bottom-paper.html

For additional product or technical information please log on to www.abdoslifesciences.com. For distributor information contact labtech@abdosindia.com

Abdos Labtech Pvt Ltd

318 - 320, DLF Tower A, Jasola, New Delhi - 110025, India.
T: +91-11-44800980-986 E: labtech@abdosindia.com