

CENTRIFUGE TUBES


BEST-IN-CLASS QUALITY
BREAKAGE-RESISTANT
DESIGNS

**Revolutionise Your
Research with ABDOS
Centrifuge Tubes**

ABDOS Centrifuge Tubes are ergonomically designed with easy-to-grip caps, manufactured from high-quality Polypropylene ensuring durability and resistance to breakage



Features:



Resins are Certified from USP Class VI
Medical grade 100 % virgin Polypropylene, conforming to US FDA 21 CFR assuring the highest quality in its class

Compatibility with Various Centrifuge Models and Rotors -
Allowing flexibility in laboratory experimental setups

THE CAPS


Flat Thread Confirms Vacuum Tight Closure – leak-proof and prevents override of cap



Unique Plug Seal – ensuring perfect sealing & make it leak-proof

Advantages of Abdos Centrifuge Tubes with Even Wall Thickness

Even wall thickness easily distributes the G Force within the tube and thus **Abdos Centrifuge Tubes can take greater RCF**



Abdos Centrifuge Tubes in the Leading Journals



ROYAL SOCIETY
OF CHEMISTRY




Springer
the language of science



bio-protocol
Bio protocol



Quality Compliances



Ordering Information

Code	Product Description	Pack Size
P10430	Centrifuge Tube, Conical, PP, 15 mL, Bulk, DNase/RNase free - 500 un. - ABDOS	500
P10401	Centrifuge Tube, Conical, PP, 15 mL, DNase/RNase free - 20 x 25 un. - ABDOS	500
P10525	Last Drop Low Retention Centrifuge Tube, Conical, PP, 15 mL, DNase/RNase free - 20 x 25 un. - ABDOS	500
P10441	Centrifuge Tube, Conical, PP, 15 mL, Red Cap, DNase/RNase free - 20 x 25 un. - ABDOS	500
P10402	Centrifuge Tube, Conical, PP, 15 mL, Sterile, DNase/RNase free - 20 x 25 un. - ABDOS	500
P10443	Centrifuge Tube, Conical, PP, 50 mL, Sterile, Red Cap, DNase/RNase free - 20 x 25 un. - ABDOS	500
P10404	Centrifuge Tube, Conical, PP, 50 mL, Sterile, DNase/RNase free - 20 x 25 un. - ABDOS	500
P10442	Centrifuge Tube, Conical, PP, 15 mL, Sterile, Red Cap, DNase/RNase free - 20 x 25 un. - ABDOS	500
P10431	Centrifuge Tube, Conical, PP, 50 mL, Bulk, DNase/RNase free - 500 un. - ABDOS	500
P10403	Centrifuge Tube, Conical, PP, 50 mL, DNase/RNase free - 20 x 25 un. - ABDOS	500
P10444	Centrifuge Tube, Conical, PP, 50 mL, Sterile, Red Cap, DNase/RNase free - 20 x 25 un. - ABDOS	500
P10405	Centrifuge Tube, Conical, PP, 15 mL, Amber, DNase/RNase free - 20 x 25 un. - ABDOS	500
P10406	Centrifuge Tube, Conical, PP, 15 mL, Sterile, Amber, DNase/RNase free - 20 x 25 un. - ABDOS	500
P10407	Centrifuge Tube, Conical, PP, 50 mL, Amber, DNase/RNase free - 20 x 25 un. - ABDOS	500
P10408	Centrifuge Tube, Conical, PP, 50 mL, Sterile, Amber, DNase/RNase free - 20 x 25 un. - ABDOS	500
P10405T	Centrifuge Tube, Conical, PP, 15 mL, Transparent Amber, DNase/RNase free - 20 x 25 un. - ABDOS	500
P10406T	Centrifuge Tube, Conical, PP, 15 mL, Sterile, Transparent Amber, DNase/RNase free - 20 x 25 un. - ABDOS	500
P10407T	Centrifuge Tube, Conical, PP, 50 mL, Transparent Amber, DNase/RNase free - 20 x 25 un. - ABDOS	500
P10408T	Centrifuge Tube, Conical, PP, 50 mL, Sterile, Transparent Amber, DNase/RNase free - 20 x 25 un. - ABDOS	500
P10419	Centrifuge Tube, Conical, PS/HDPE, 15 mL, Sterile, Cell Culture, DNase/RNase free - 20 x 25 un. - ABDOS	500
P10421	Centrifuge Tube, Conical, Racked, PS/HDPE, 15 mL, Sterile, Cell Culture, DNase/RNase free - 20 x 25 un. - ABDOS	500 (25 PCS X 20 RACKS)
P10432	Centrifuge Tube, Conical, PP, 50 mL, Self-Standing, Bulk, DNase/RNase free - 500 un. - ABDOS	500
P10423	Centrifuge Tube, Conical, PP, 50 mL, Self-Standing, DNase/RNase free - 20 x 25 un. - ABDOS	500
P10449	Centrifuge Tube, Conical, PP, 50 mL, Self-Standing, Red Cap, DNase/RNase free - 20 x 25 un. - ABDOS	500
P10424	Centrifuge Tube, Conical, PP, 50 mL, Self-Standing, Sterile, DNase/RNase free - 20 x 25 un. - ABDOS	500
P10450	Centrifuge Tube, Conical, PP, 50 mL, Self-Standing, Sterile, Red Cap, DNase/RNase free - 20 x 25 un. - ABDOS	500

What Abdos Life Sciences stands for?



Precision



Reliability



Responsiveness



Innovation



Sustainability



Relationships

Get in touch with us:

011-41081897



labtech@abdosindia.com



www.abdoslifesciences.com