

# Technical Data Sheet

Doc No – Ab/TDS/MRB-130

Rev No: 03 Date 01/2021

## TECHNICAL DATA SHEET

**Abdos Labtech Pvt. Ltd.**

**FROSTY MATE™ Freezer Storage Racks, PC**

**Valid For Cat. No. P11104, P11105**

### **1. Description/Specification**

#### **1.1 Description:** Freezer Storage Racks, PC – FROSTY MATE™

1.1.1 These racks have consistent thermal protection to maintain temperature at low temperature.

1.1.2 These racks are suitable for storage of PCR reagents, cells, enzymes or other temperature sensitive reagents.

1.1.3 These racks are available in two temperature range:

(A): -20°C unit will maintain a temperature between -10°C and -20°C for up to 3 hours.

(B): 0°C unit will maintain a temperature of 0°C for 5 hours.

1.1.4 These racks are filled with an inert, non-toxic gel.

#### **1.2 Design:**

Yellow color lid permit viewing indexing makes easy to identification.

#### **1.3 Performance Specifications:**

1.3.1 **P11104:** -20°C Freezer Storage racks (240x120x120mm)

**Capacity:** 0.5, 1.5-2 Micro tubes & 1-2 Cryo tubes

**Places:** 32(8x4 Array) with 12 adaptors

1.3.2 **P11105:** 0°C Freezer Storage racks (240x120x120mm)

**Capacity:** 0.5, 1.5-2 Micro tubes & 1-2 Cryo tubes

**Places:** 32(8x4 Array) with 12 adaptors

#### **1.4 Material:**

1.4.1 Polycarbonate (PPC, Hydrophobic, confirm US FDA 21CFR)



**Shiv Ganga Industrial Estate,**  
Plot no.1 Lakeshwari, Bhagwanpur,  
Roorkee – 247667 Uttarakhand, India  
**An ISO 9001 ,13485, 14001 & CE Certified**  
A Division of Abdos

**Regd. Office:**  
13A, Government Place (East)  
Kolkata – 700069, India  
**GST : 05AAGCA8605A1Z7**

**P :** +91 (0) 1332 235277  
**E :** labtech@abdosindia.com  
**W:** www.abdoslifesciences.com  
**CIN :** U85100WB2008PTC123903

# Technical Data Sheet

1.4.2 Complies with requirement accordance with article of regulation (EC) No. 1907/2006 and its amendments.

1.4.3 Complies with the requirement of European Pharmacopoeia, 7<sup>th</sup> edition (2011) & its supplements 7.5 (07/2012) monographs 3.2.2 Plastic Containers and closures for Pharmaceutical use

1.4.4 The product has passed the USP testing including class VI tests & has been assigned the FDA Drug master file DMF 12123

1.4.5 Materials has successful passed the biological tests according to ISO 10993 external communication devices for indirect blood contact for Prolonged period.

1.4.6 Free from Heavy metals, RoHS Free according to meets relevant requirements of EU regulations 1935/2004/CE, 2002/95/EC as amended, UE 10/2011 as amended. EU CONEG & it's following amendments. It is safe use in contact with food stuff, Pharmaceuticals & drinking water

1.4.7 Product is manufactured without the use of Raw Materials of animal origin. Therefore BSE/TSE issue does not concern this product to the best of our Knowledge. No materials containing latex or natural rubber are used in the Manufacturing, handling & packaging process for this product

**1.5 Color:** Yellow

**1.6 Indexing:** Molded indexing makes identification faster & easier.

**1.7 Adaptor** for using 0.5ml Micro tube.

**1.8 Quality control:**

Product are manufactured, inspected and found to be in compliance with product & quality specification requirements as documented in ISO 9001:2015 & 13485:2016 Quality Management System. Testing of various characteristics in accordance with the valid specification, inspection reviewed & signed for product release.

**1.9. Basic Features**

**1.9.1 Temperature range:** 0°C to -20°C

**1.9.2 Autoclavability:** NA

**1.9.3 Chemical resistance:** See website ([www.abdoslifesciences.com](http://www.abdoslifesciences.com))

**1.9.4** Manufactured in Aseptic Environment

**2.0 Other Information:**

Reusable




**Shiv Ganga Industrial Estate,**  
Plot no.1 Lakeshwari, Bhagwanpur,  
Roorkee – 247667 Uttarakhand, India  
**An ISO 9001 ,13485, 14001 & CE Certified**  
A Division of Abdos

**Regd. Office:**  
13A, Government Place (East)  
Kolkata – 700069, India  
**GST : 05AAGCA8605A1Z7**

**P :** +91 (0) 1332 235277  
**E :** labtech@abdosindia.com  
**W:** www.abdoslifesciences.com  
**CIN :** U85100WB2008PTC123903

# Technical Data Sheet

**NOTE:-** FROSTY MATE™ Freezer Storage Rack should be kept 24 hours in the deep freeze before the use to get improved result.





**Shiv Ganga Industrial Estate,**  
Plot no.1 Lakeshwari, Bhagwanpur,  
Roorkee – 247667 Uttrakhand, India  
**An ISO 9001 ,13485, 14001 & CE Certified**  
A Division of Abdos


**Regd. Office:**  
13A, Government Place (East)  
Kolkata – 700069, India  
**GST : 05AAGCA8605A1Z7**

**P : +91 (0) 1332 235277**  
**E : labtech@abdosindia.com**  
**W: www.abdoslifesciences.com**  
**CIN : U85100WB2008PTC123903**

# Technical Data Sheet



  
Quality Assurance Incharge  
Abdos Labtech Pvt. Ltd.



  
**ABDOS**

**Shiv Ganga Industrial Estate,**  
Plot no.1 Lakeshwari, Bhagwanpur,  
Roorkee – 247667 Uttarakhand, India

**An ISO 9001 ,13485, 14001 & CE Certified**

A Division of Abdos

**Regd. Office:**

13A, Government Place (East)  
Kolkata – 700069, India

**GST : 05AAGCA8605A1Z7**

**P : +91 (0) 1332 235277**

**E : labtech@abdosindia.com**

**W: www.abdoslifesciences.com**

**CIN : U85100WB2008PTC123903**