

**ABDOS**  
LIFE SCIENCE

**10**  
YEARS

**ABDOS**  
LIFE SCIENCE



# Last Drop™ Filter Pipette Tips

Creating For Life Sciences

[www.abdoslifesciences.com](http://www.abdoslifesciences.com)



- Minimum retention with viscous liquid thus reducing sample retention and save on costly reagents.
- Superior non-leaching performance in its class assures contamination free results.
- Optimal surface for enhanced dispensing and minimal sample loss of highly viscous liquids or samples containing detergents.
- Precise graduation marks for accurate pipetting.
- Diamond polished moulds creates a hydrophobic surface ensuring the lowest retention in its class.
- Highly homogeneous surface prevents non-specific binding.
- Suited for biological application.
- Abdos Filter is 100% pure porous PE. It makes Recovery of Samples from Wetted, Filter without Contamination.
- Non Cytotoxic (ISO 10993-5 MRM, Minimal Essential Media Elution, GLP Compliant)  
99.9% Bacterial Aerosol Filtration Efficiency (BFE), (ASTM F21012)
- Non Haemolytic (Modified ASTM F 756-08 Haemolysis Direct Contact GLP with Human Blood)
- Inert, Hydrophobic Filter does not trap or block the liquid.

# ABDOS ENSURES LOWEST SAMPLE RETENTION WITH ITS LAST DROP FILTER PIPETTE TIPS



- Precise safety gap between filter and liquid, filter and pipette shaft to avoid contamination & worry free pipetting.
- Abdos super hydrophobic material ensures retention level 10 times lesser with viscous samples than standard pipette tips to save an costly reagents.
- Abdos Filter are 100% pure hydrophobic polyethylene (PE) makes recovery of samples from wetted, filtered without contamination.
- 5 Stage quality control procedures ensures 100% quality & performance.
- Super hydrophobic filter that does not absorb samples, thus save your precious sample loss in pipetting error.
- Sterile tips has been validated according to ENISO 11137-2; 2012 for SAL  $10^{-6}$ .
- Abdos filter pipette tips are manufactured from ultra pure USP class VI certified virgin polypropylene to ensure the highest quality & confirming to US FDA 21 CFR. Free from heavy metals, natural rubber & inhibitory plasticizers.
- Proprietary polyethylene super filter provides aerosol barrier and aqueous liquid by-pass prevention in over-pipetting hereby preventing cross contamination.
- Best suited for biological application.



Filters used in Abdos Last Drop filter pipette tips are non cytotoxic (ISO 10993-5 MRM, minimal essential media elution, GLP compliant). Abdos filters provide 99.9% bacterial aerosol filtration efficiency (BFE), (ASTM F21012), are non hemolytic (modified ASTM F 756-08 hemolytic direct contact GLP with human blood), inert and hydrophobic filter does not trap or block the liquid.

Code	Product Description	No.Per Case
P11061	ABDOS Filter Pipette Tips, Racked Sterile, Graduated, 0.1 – 10 µl	10 Racks x 96 pcs (960 pcs)
P11063	ABDOS Filter Pipette Tips, Racked Sterile, Extended Graduated, 0.1 – 10 µl	10 Racks x 96 pcs (960 pcs)
P11065	ABDOS Filter Pipette Tips, Racked Sterile, Graduated, 1 – 20 µl	10 Racks x 96 pcs (960 pcs)
P11067	ABDOS Filter Pipette Tips, Racked Sterile, Graduated, 1 – 50 µl	10 Racks x 96 pcs (960 pcs)
P11069	ABDOS Filter Pipette Tips, Racked Sterile, Graduated, 1 – 100 µl	10 Racks x 96 pcs (960 pcs)
P11073	ABDOS Filter Pipette Tips, Racked Sterile, Graduated, 1 – 200 µl	10 Racks x 96 pcs (960 pcs)
P11075	ABDOS Filter Pipette Tips, Racked Sterile, Graduated, 1 – 300 µl	10 Racks x 96 pcs (960 pcs)
P11077	ABDOS Filter Pipette Tips, Racked Sterile, 100 – 1000 µl	10 Racks x 96 pcs (960 pcs)
P11079	ABDOS Filter Pipette Tips, Racked Sterile, Graduated, 100 – 1000 µl XL	10 Racks x 96 pcs (960 pcs)
P11081	ABDOS Filter Pipette Tips, Racked Sterile, Graduated, 100 – 1250 µl	10 Racks x 96 pcs (960 pcs)
P11060	ABDOS Filter Pipette Tips, Racked, Graduated, 0.1 – 10 µl	10 Racks x 96 pcs (960 pcs)
P11062	ABDOS Filter Pipette Tips, Racked Extended Graduated, 0.1 – 10 µl	10 Racks x 96 pcs (960 pcs)
P11064	ABDOS Filter Pipette Tips, Racked, Bevelled Graduated, 1 – 20 µl	10 Racks x 96 pcs (960 pcs)
P11066	ABDOS Filter Pipette Tips, Racked, Bevelled Graduated, 1 – 50 µl	10 Racks x 96 pcs (960 pcs)
P11068	ABDOS Filter Pipette Tips, Racked, Bevelled Graduated, 1 – 100 µl	10 Racks x 96 pcs (960 pcs)
P11072	ABDOS Filter Pipette Tips, Racked, Graduated, 1 – 200 µl	10 Racks x 96 pcs (960 pcs)
P11074	ABDOS Filter Pipette Tips, Racked, 100 – 300 µl	10 Racks x 96 pcs (960 pcs)
P11076	ABDOS Filter Pipette Tips, Racked, 100 – 1000 µl	10 Racks x 96 pcs (960 pcs)
P11078	ABDOS Filter Pipette Tips, Racked, Graduated, 100 – 1000 µl XL	10 Racks x 96 pcs (960 pcs)
P11080	ABDOS Filter Pipette Tips, Racked, Graduated, 100 – 1250 µl	10 Racks x 96 pcs (960 pcs)



**Last Drop™ Low Retention Filter Pipette Tips, Bulk in Resealable Bags**



**Last Drop™ Low Retention Filter Pipette Tips, Racked**



**Last Drop™ Low Retention Filter Pipette Tips, Sterile Racked**



**Last Drop™ Low Retention Filter Pipette Tips, Refill Pack**