

Abdos PET PETG Media Bottles

EASY-TO-HOLD SQUARE BOTTLES

Narrow mouth to provide better control when pouring

Square shape, easy-to-hold, and saves shelf counter space compared with round bottles of the same volume


Lightweight impact-resistant material that can withstand $-70^{\circ}\text{C} \sim +50^{\circ}\text{C}$

Features:

Available in two polymer variants: Polyethylene Terephthalate Glycol co-polymer (PETG) and Polyethylene Terephthalate (PET)

Moulded embossed graduations for easy visual observation of the solution volume

Highly transparent and optically clear with smooth inner & outer walls



Abdos PET & PETG Media Bottles have excellent gas barrier properties, making them best-suited for storing and sampling active pharmaceutical ingredients and bulk intermediates.

These media bottles are also ideal for preparing and storing buffers, serum, culture fluids, or long-term storage of pH-sensitive liquids.

Wide Volume Options Available



60 ml
(PETG)



125 ml
(PET/PETG)



250 ml
(PET/PETG)




500 ml
(PET/PETG)



1000 ml
(PET/PETG)

Certifications and Compliances

SAL 10⁻⁶



FREE



FREE



FREE



FREE



FREE

Ordering Information

PETG Media Bottles		
Code	Description	Nos/Case
P21170	60ml, Non-Sterile	192
P21171	125ml, Non-Sterile	96
P21172	250ml, Non-Sterile	60
P21173	500ml, Non-Sterile	48
P21174	1000ml, Non-Sterile	24
P21180	60ml, Sterile	192
P21181	125ml, Sterile	96
P21182	250ml, Sterile	60
P21183	500ml, Sterile	48
P21184	1000ml, Sterile	24

PET Media Bottles		
Code	Description	Nos/Case
P21176	125ml, Non-Sterile	96
P21177	250ml, Non-Sterile	224
P21178	500ml, Non-Sterile	120
P21179	1000ml, Non-Sterile	40
P21186	125ml, Sterile	96
P21187	250ml, Sterile	224
P21188	500ml, Sterile	120
P21189	1000ml, Sterile	40


Other Items of Interest in Cell Culture Range




BIOFILL Bottle
Top Filtration System



Cell Culture Plates



Cell Culture Dishes



Cell Culture Flasks

What Abdos Life Sciences stands for?



Precision



Reliability



Responsiveness



Innovation



Sustainability



Relationships

Get in touch with us:

☎ 011-41081897

✉ labtech@abdosindia.com

🌐 www.abdoslifesciences.com