

Technical Data Sheet

Doc No Ab/TDS/ELECTRO-124

Rev No: 01 Date 10/2021-22

TECHNICAL DATA SHEET

Abdos Labtech Pvt. Ltd.

Microprocessor Based Programmable Power Supply

Valid for Catalog No. E12377

1. Description Specifications:

1.1 Description: Microprocessor Based Programmable Power Supplies 0-300V;0-400mA

- 1.1.1 Compact and Lightweight
- 1.1.2 LCD Digital Display & Digital Control
- 1.1.3 Runs Four Gels simultaneously
- 1.1.4 Short circuit & overload protection
- 1.1.5 Pause Key function to hold the output power
- 1.1.6 Auto Recovery after Power failure
- 1.1.7 Alarm & Memory option

1.2 Design: Power Supply

1.3 Performance Specifications:

Specifications	E12377
Body	MS with Powder Coated Body
Dimension (LxWxH)	240 x 260 x 120 mm
Power Requirement	220-240 V, 50-60 Hz, AC
Output Voltage	0 to 300 V
Output Current	0 to 400 mA
Timer	0 minutes to 24 Hrs
Parameters	16 programmable
Display	LCD Digital Display
Output Terminal	4 Parallel (4mm Socket)

1.4 Material: MS with Power Coated Body

1.5 CE Certified



Shiv Ganga Industrial Estate,
Plot no.1 Lakeshwari, Bhagwanpur,
Roorkee – 247667 Uttarakhand, India
An ISO 9001 ,13485, 14001 & CE Certified
A Division of Abdos

Regd. Office:
13A, Government Place (East)
Kolkata – 700069, India
GST : 05AAGCA8605A127

P : +91 (0) 1332 235277
E : labtech@abdosindia.com
W: www.abdoslifesciences.com
CIN : U85100WB2008PTC123903

Technical Data Sheet

1.6 Quality control:

Product are manufactured, inspected and found to be in compliance with product & quality specification requirements as documented in ISO 9001:2015 & 13485:2016 Quality Management System. Testing of various characteristics in accordance with the valid specification, inspection reviewed & signed for product release

Warranty: 12 months from the date of Installation or 15 months from the date of Invoicing



Shiv Ganga Industrial Estate,
Plot no.1 Lakeshwari, Bhagwanpur,
Roorkee – 247667 Uttarakhand, India

An ISO 9001 ,13485, 14001 & CE Certified

A Division of Abdos

Regd. Office:

13A, Government Place (East)
Kolkata – 700069, India

GST : 05AAGCA8605A1Z7

P : +91 (0) 1332 235277

E : labtech@abdosindia.com

W: www.abdoslifesciences.com

CIN : U85100WB2008PTC123903