

Technical Data Sheet

Doc No–Ab/TDS/Micro-88

Rev No: 02 Date 04/2023-24

TECHNICAL DATA SHEET

Abdos Labtech Pvt. Ltd.

Petridish Stand, Acrylic

Valid for Catalog No. P10970

SRN (Single Registration Number): IN-MF-000022272

1. Description/Specification

1.1 Description: Petridish Stand, Acrylic

1.1.1 Can hold & secure glass & plastic petridishes.

1.1.2 This stand makes transport safe & easier.

1.1.3 Clear racks give the visibility to check culture during incubation.

1.2 Design: Rectangular stand

1.3 Performance Specifications:

Code	No. of Places	L x W x H (mm)	Compatible	No. Per Case
P10970	28	252 x 105 x 206	90mm Petridish	1

1.4 Material:

1.4.1 Made up of Acrylic, confirm US FDA 21CFR

1.4.2 Complies with the requirement of European Pharmacopoeia, Plastic Containers, and closures for pharmaceutical use.

1.4.3 The product has passed the USP testing including class VI tests, Materials has successful passed the biological tests according to ISO 10993 external communication devices for indirect blood contact for Prolonged period.

1.4.4 Free from Heavy metals, RoHS Free according to meets relevant requirements of EU regulations 1935/2004/CE, 2002/95/EC as amended, UE 10/2011 as amended. EU CONEG & it's following amendments. It is safe use in contact with food stuff, Pharmaceuticals & drinking water.

1.4.5 Product is manufactured without the use of Raw Materials of animal origin. Therefore BSE/TSE issue does not concern this product to the best of our knowledge. No materials containing latex or natural rubber are used in the Manufacturing, handling & packaging process for this product.

1.5 Color: Transparent Natural



Abdos Labtech Private Limited
(A Division of Abdos Group)
(An ISO 9001,13485,14001,18385,GMP & CE certified company)
Shiv Ganga Industrial Estate
Plot No-1, Lakheswari, Bhagwanpur,
Roorkee-247667, (UK) INDIA

Registered office
13A Government
Place East, 2nd Floor Kolkata
West Bengal – 700069
CIN: U85100WB2008PTC123903
Corporate Office
319, 3rd Floor, DLF Tower A, DDA Dist. Centre
Jasola, New Delhi - 110025

Ph. No: +91- (0) 1332235277
Email: labtech@abdosindia.com
Web: www.abdoslivesciences.com
GST: 05AAGCA8605A1Z7

Technical Data Sheet

1.6 Quality control:

Product are manufactured, inspected and found to be in compliance with product & quality specification requirements as documented in ISO 9001:2015 & 13485:2016 Quality Management System. Testing of various characteristics in accordance with the valid specification, inspection reviewed & signed for product release.

1.7 Basic Features

1.7.1 Autoclavability: NO

1.7.2 Chemical resistance: See website (www.abdoslifesciences.com)

1.7.3 Manufactured in Aseptic Environment

Instruction for Use

Product can be used for food contact also.
Storage temp. Max 30°C
Before using the product check crackness, etc.
Check all the information printed on pouches or boxes.
Always wear gloves while handling specimen.

WARNING

DO NOT store strong oxidizing agents.
DO NOT place any plastic labware in a flame.
DO NOT mix any chemicals that may result in a thermic reaction, which can cause product failure
Consult local fire codes prior to storage of flammable liquids.



Quality Manager

Abdos Labtech Private Limited

Employee Code 30110

