

SEROLOGICAL PIPETTES


BEST-IN-CLASS ULTRA-CLEAR AND TRANSPARENT PIPETTES

Features

- ▶ Polystyrene polymer for maximum transparency
- ▶ Easy-to-read high contrast graduations
- ▶ Ascending | descending graduation allows volume to be read at glance
- ▶ Negative graduation offers increased volume capacity thereby providing more flexibility
- ▶ Graduated pipette marked with a pipette colour code for easy volume identification

Volume Capacity (ml)	Pipette Colour Code Marking
1	Yellow
2	Green
5	Blue
10	Orange
25	Red
50	Black

- ▶ PE filters plug for protection against liquid aspiration into the pipetting device while further reducing sample-to-sample contamination risks
- ▶ Volumetric accuracy +/- 2%



Serological Pipettes

Abdos Life Sciences introduces the best-in-class, ultra-clear and transparent serological pipette made from polystyrene with working volume variant, ranging from 1 ml to 50 ml.

Application: Whether you are working in a tissue culture lab or working with a non-toxic solution that requires aliquoting or aspirating | dispensing in a sterile environment, Abdos Serological Pipettes come quite handy.

Quality & Compliance




FREE



FREE



FREE



FREE

STERILE R

Key Points

- Working temperature range -10°C to $+70^{\circ}\text{C}$
- Colour coding on individual packaging
- Gamma irradiated with SAL 10^{-6}

Packaging Format

- Individual poly-paper packaging with peel off function
- Bulk packaging for an ecological and economical solution

Ordering Information


Product Code	Product Description	No. Case
P10601	Serological Pipettes Sterile, PS, Capacity 1 ml	500
P10607	Serological Pipettes, Sterile, PS, 1 ml, Bulk Sleeve of 100	500
P10602	Serological Pipettes Sterile, PS, Capacity 2 ml	400
P10608	Serological Pipettes, Sterile, PS, 2 ml, Bulk Sleeve of 100	500
P10603	Serological Pipettes Sterile, PS, Capacity 5 ml	200
P10609	Serological Pipettes, Sterile, PS, 5 ml, Bulk Sleeve of 25	200
P10604	Serological Pipettes Sterile, PS, Capacity 10 ml	200
P10610	Serological Pipettes, Sterile, PS, 10 ml, Bulk Sleeve of 25	200
P10605	Serological Pipettes Sterile, PS, Capacity 25 ml	200
P10611	Serological Pipettes, Sterile, PS, 25ml, Bulk Sleeve of 25	200
P10606	Serological Pipettes Sterile, PS, Capacity 50 ml	100
P10612	Serological Pipettes, Sterile, PS, 50ml, Bulk Sleeve of 25	100

Full traceability on individual packed serological pipette via lot number and expiry date

Scan for more information



Other products of Interest



Cell culture plate



Cell Culture Flask



Petridish



Erlenmeyer Flask

What Abdos Life Sciences stands for?



Precision



Reliability



Responsiveness



Innovation



Sustainability



Relationships

Get in touch with us:

☎ 011-41081897

✉ labtech@abdosindia.com

🌐 www.abdoslifesciences.com