

# Technical Data Sheet

Doc No Ab/TDS/CT-24

Rev No: 00 Date 02/2022-23

## TECHNICAL DATA SHEET

**Abdos Labtech Pvt. Ltd.**

**Swirl Mini Plate Centrifuge**

**Valid for Catalog No. E11560**

### **1. Description Specifications:**

#### **1.1 Description: Swirl Mini Plate Centrifuge**

1.1.1 Unique design for the horizontal rotor centrifuge. Short-spin and timing of two operations make the tests more convenient.

1.1.2 Braking time from max. speed: less than 30sec.Timer can be set up to 10 minutes.

1.1.3 Streamline design with more compact shape and transparent PC lid.

1.1.4 It can accept below formats:

PCR standard plates with skirted, non-skirted

Elisa Plates

8x0.2ml PCR strip tubes

0.2ml PCR tubes

#### **1.2 Design: Swirl Mini Plate Centrifuge**

#### **1.3 Performance Specifications:**

<b>Specifications</b>	<b>E11560</b>
<b>Voltage</b>	220V/120V
<b>Frequency</b>	50-60 Hz
<b>Rotor Capacity</b>	2x2 PCR plates, 2xMicroliter plates, 24x0.2ml 8 strops PCR tubes 96x0.2ml tubes
<b>Max Speed</b>	2200 rpm /480x G
<b>Deceleration</b>	Less than 30 sec



Abdos Labtech Private Limited  
(A Division of Abdos Group)  
(An ISO 9001,13485,14001,18385,GMP & CE certified company)  
Shiv Ganga Industrial Estate  
Plot No-1, Lakheswari, Bhagwanpur,  
Roorkee-247667, (UK) INDIA

**Registered office**  
13A Government  
Place East, 2nd Floor Kolkata  
West Bengal – 700069  
CIN: U85100WB2008PTC123903  
**Corporate Office**  
319, 3rd Floor, DLF Tower A, DDA Dist. Centre  
Jasola, New Delhi - 110025

Ph. No: +91- (0) 1332235277  
Email:labtech@abdosindia.com  
Web: [www.abdoslivesciences.com](http://www.abdoslivesciences.com)  
GST: 05AAGCA8605A1Z7

# Technical Data Sheet

<b>Dimension</b>	290x360x140mm
<b>Weight</b>	3.6 kg

**1.4 Material:** ABS (Acrylonitrile Butadiene Styrene)

**1.5 Color:** White

**1.6 E11560** – Swirl Mini Plate Centrifuge, CE certified

**1.7 Spare Accessories**

1.7.1 **E11561** – Adaptor for PCR plates with skirt and without skirt

1.7.2 **E11562** – Adaptor for 8x0.2ml strip or 0.2ml PCR tubes

**Warranty:** 12 months from the date of Installation or 15 months from the date of Invoicing.

## Instruction for Use

Storage temp. Max 30°C

Before using the product check .

Check all the information printed on pouches or boxes.

Always wear gloves while handling specimen.

## WARNING

**DO NOT** store strong oxidizing agents.

**DO NOT** place any plastic labware in a flame.

**DO NOT** mix any chemicals that may result in a thermic reaction, which can cause product failure

Consult local fire codes prior to storage of flammable liquids.



**Abdos Labtech Pvt. Ltd.**

Plot No. 1 Lakshwari Bhagwanpur Roorkee 247667 (UK) India  
Email : labtech@abdosindia.com, www.abdoslifesciences.com  
Ph: +(91) 11 41081897



**Abdos Labtech Private Limited**  
(A Division of Abdos Group)  
(An ISO 9001,13485,14001,18385,GMP & CE certified company)  
Shiv Ganga Industrial Estate  
Plot No-1, Lakshwari, Bhagwanpur,  
Roorkee-247667, (UK) INDIA

**Registered office**  
13A Government  
Place East, 2nd Floor Kolkata  
West Bengal – 700069  
CIN: U85100WB2008PTC123903  
**Corporate Office**  
319, 3rd Floor, DLF Tower A, DDA Dist. Centre  
Jasola, New Delhi - 110025

**Ph. No: +91- (0) 1332235277**  
**Email:labtech@abdosindia.com**  
**Web: [www.abdoslifesciences.com](http://www.abdoslifesciences.com)**  
**GST: 05AAGCA8605A1Z7**

# Technical Data Sheet



**Abdos Labtech Private Limited**  
(A Division of Abdos Group)  
(An ISO 9001,13485,14001,18385,GMP & CE certified company)  
Shiv Ganga Industrial Estate  
Plot No-1, Lakheswari, Bhagwanpur,  
Roorkee-247667, (UK) INDIA

**Registered office**  
13A Government  
Place East, 2nd Floor Kolkata  
West Bengal – 700069  
CIN: U85100WB2008PTC123903  
**Corporate Office**  
319, 3rd Floor, DLF Tower A, DDA Dist. Centre  
Jasola, New Delhi - 110025

**Ph. No: +91- (0) 1332235277**  
**Email: labtech@abdosindia.com**  
**Web: [www.abdoslivesciences.com](http://www.abdoslivesciences.com)**  
**GST: 05AAGCA8605A1Z7**