

Technical Data Sheet

Doc No Ab/TDS/CT-24

Rev No: 01 Date 04/2023-24

TECHNICAL DATA SHEET

Abdos Labtech Pvt. Ltd.

Swirl Mini Plate Centrifuge

Valid for Catalog No. E11560

SRN (Single Registration Number): IN-MF-000022272

1. Description Specifications:

1.1 Description: Swirl Mini Plate Centrifuge

1.1.1 Unique design for the horizontal rotor centrifuge. Short-spin and timing of two operations make the tests more convenient.

1.1.2 Braking time from max. speed: less than 30sec.Timer can be set up to 10 minutes.

1.1.3 Streamline design with more compact shape and transparent PC lid.

1.1.4 It can accept below formats:

PCR standard plates with skirted, non-skirted

Elisa Plates

8x0.2ml PCR strip tubes

0.2ml PCR tubes

1.2 Design: Swirl Mini Plate Centrifuge

1.3 Performance Specifications:

Specifications	E11560
Voltage	220V/120V
Frequency	50-60 Hz
Rotor Capacity	2x2 PCR plates, 2xMicroliter plates, 24x0.2ml 8 strops PCR tubes 96x0.2ml tubes
Max Speed	2200 rpm /480x G



Abdos Labtech Private Limited
(A Division of Abdos Group)
(An ISO 9001,13485,14001,18385,GMP & CE certified company)
Shiv Ganga Industrial Estate
Plot No-1, Lakheswari, Bhagwanpur,
Roorkee-247667, (UK) INDIA

Registered office
13A Government
Place East, 2nd Floor Kolkata
West Bengal – 700069
CIN: U85100WB2008PTC123903
Corporate Office
319, 3rd Floor, DLF Tower A, DDA Dist. Centre
Jasola, New Delhi - 110025

Ph. No: +91- (0) 1332235277
Email:labtech@abdosindia.com
Web: www.abdoslivesciences.com
GST: 05AAGCA8605A1Z7

Technical Data Sheet

Deceleration	Less than 30 sec
Dimension	290x360x140mm
Weight	3.6 kg

1.4 Material: ABS (Acrylonitrile Butadiene Styrene)

1.5 Color: White

1.6 E11560 – Swirl Mini Plate Centrifuge, CE certified

1.7 Spare Accessories

1.7.1 **E11561** – Adaptor for PCR plates with skirt and without skirt

1.7.2 **E11562** – Adaptor for 8x0.2ml strip or 0.2ml PCR tubes

Warranty: 12 months from the date of Installation or 15 months from the date of Invoicing.

Instruction for Use

Storage temp. Max 30°C

Before using the product check .

Check all the information printed on pouches or boxes.

Always wear gloves while handling specimen.

WARNING

DO NOT store strong oxidizing agents.

DO NOT place any plastic labware in a flame.

DO NOT mix any chemicals that may result in a thermic reaction, which can cause product failure

Consult local fire codes prior to storage of flammable liquids.



Abdos Labtech Pvt. Ltd.

Plot No. 1 Lakeshwari Bhagwanpur Roorkee 247667 (UK) India
Email : labtech@abdosindia.com, www.abdoslifesciences.com
Ph: +(91) 11 41081897



Abdos Labtech Private Limited
(A Division of Abdos Group)
(An ISO 9001,13485,14001,18385,GMP & CE certified company)
Shiv Ganga Industrial Estate
Plot No-1, Lakshewari, Bhagwanpur,
Roorkee-247667, (UK) INDIA

Registered office
13A Government
Place East, 2nd Floor Kolkata
West Bengal – 700069
CIN: U85100WB2008PTC123903

Corporate Office
319, 3rd Floor, DLF Tower A, DDA Dist. Centre
Jasola, New Delhi - 110025

Ph. No: +91- (0) 1332235277
Email:labtech@abdosindia.com
Web: www.abdoslifesciences.com
GST: 05AAGCA8605A1Z7

Technical Data Sheet

Quality Manager

Abdos Labtech Private Limited

Employee Code 30110



Abdos Labtech Private Limited
(A Division of Abdos Group)
(An ISO 9001,13485,14001,18385,GMP & CE certified company)
Shiv Ganga Industrial Estate
Plot No-1, Lakheswari, Bhagwanpur,
Roorkee-247667, (UK) INDIA

Registered office
13A Government
Place East, 2nd Floor Kolkata
West Bengal – 700069
CIN: U85100WB2008PTC123903
Corporate Office
319, 3rd Floor, DLF Tower A, DDA Dist. Centre
Jasola, New Delhi - 110025

Ph. No: +91- (0) 1332235277
Email: labtech@abdosindia.com
Web: www.abdoslifesciences.com
GST: 05AAGCA8605A1Z7