

TECHNICAL DATA SHEET

Doc No – Ab/TDS/CT-23
Rev No: 02 Date- 05/2024-25

Abdos Labtech Pvt. Ltd.
Swirl Mini Plate Centrifuge

Valid for Catalog No. E11560

SRN (Single Registration Number): IN-MF-000022272

1. Description Specifications:

1.1 Description: Swirl Mini Plate Centrifuge

1.1.1 Unique design for the horizontal rotor centrifuge. Short-spin and timing of two operations make the tests more convenient.

1.1.2 Braking time from max. speed: less than 30sec. Timer can be set up to 10 minutes.

1.1.3 Streamline design with more compact shape and transparent PC lid.

1.1.4 It can accept below formats:

**PCR standard plates with skirted, non-skirted Elisa Plates 8x0.2ml PCR strip tubes
0.2ml PCR tubes**

1.2 Design: Swirl Mini Plate Centrifuge

1.3 Performance Specifications:

Specifications	E11560
Voltage	220V/120V
Frequency	50-60 Hz
Rotor Capacity	2x2 PCR plates, 2xMicroliter plates, 24x0.2ml 8 strops PCR tubes 96x0.2ml tubes

TECHNICAL DATA SHEET

Doc No – Ab/TDS/CT-23
Rev No: 02 Date- 05/2024-25

Abdos Labtech Pvt. Ltd. Swirl Mini Plate Centrifuge

Max Speed	2200 rpm /480x G
Deceleration	Less than 30 sec
Dimension	290x360x140mm
Weight	3.6 kg

1.4 Material: ABS (Acrylonitrile Butadiene Styrene)

1.5 Color: White

1.6 E11560 – Swirl Mini Plate Centrifuge, CE certified

1.7 Spare Accessories

1.7.1 **E11561** – Adaptor for PCR plates with skirt and without skirt

1.7.1 **E11562** – Adaptor for 8x0.2ml strip or 0.2ml PCR tubes

1.8 Intended use: Widely used in scientific studies of molecular biology, clinical chemistry, biochemistry and any type of scientific laboratory.

Warranty: 12 months from the date of Installation or 15 months from the date of Invoicing.

Instruction for Use

Storage temp. Max 30°C

Before using the product check .

Check all the information printed on pouches or boxes.

Always wear gloves while handling specimen.


TECHNICAL DATA SHEET

Doc No – Ab/TDS/CT-23
Rev No: 02 Date- 05/2024-25

Abdos Labtech Pvt. Ltd. Swirl Mini Plate Centrifuge

WARNING

DO NOT store strong oxidizing agents.
DO NOT place any plastic labware in a flame.
DO NOT mix any chemicals that may result in a thermic reaction, which can cause product failure
Consult local fire codes prior to storage of flammable liquids.



Quality Manager
Abdos Labtech Private Limited
Employee Code 30110



TECHNICAL DATA SHEET

Doc No – Ab/TDS/CT-23
Rev No: 02 Date- 05/2024-25

Abdos Labtech Pvt. Ltd.
Swirl Mini Plate Centrifuge

