

TECHNICAL DATA SHEET

Doc No – Ab/TDS/Eqpt-233
Rev No: 09 Date- 05/2024-25

Abdos Labtech Pvt. Ltd. SWIRLEX – VOTREX SHAKER

Valid for Catalog No. E11190

SRN (Single Registration Number): IN-MF-000022272

1. Description/Specification

1.1 Description:

- 1.1.1 Swirlex-Vortex is suitable for touch or continuous or fixed operation.
- 1.1.2 Used for gentle mixing to vigorous re-suspension of biological cells and canonical liquid components in tubes using eccentric mechanism.
- 1.1.3 These vortex shakers contain specially designed silicone feet and steel base excellent for damping vibration.
- 1.1.4 Swirlex vortex Shaker is eccentric with oil less ball bearing. A wide range of accessories are also available.

1.2 Design: SWIRLEX – VOTREX SHAKER

1.3 Performance Specifications:

1.3.1 Speed Range	: 0-3000 rpm.
1.3.2 Voltage	: 220 ~ 240 AC
1.3.3 Frequency	: 50-60 Hz
1.3.4 Power	: 60 W
1.3.5 Shaking movement	: Orbital
1.3.6 Outer Dia	: 4 mm
1.3.7 Speed display	: Analog Type
1.3.8 Run type	: Both Mode – continuous, touch type

TECHNICAL DATA SHEET

Doc No – Ab/TDS/Eqpt-233
Rev No: 09 Date- 05/2024-25

Abdos Labtech Pvt. Ltd. SWIRLEX – VOTREX SHAKER

1.3.8 Permissible Ambient Temp. : 5-40 ⁰ C
1.3.9 Permissible Relative humidity: 80 %
1.4.0 Weight : 3.5 kg

1.4 Material: Body – ABS

1.5 Color : Grey

1.6. Basic features:

1.6.1 Material: Made up of highly pure Virgin ABS (Acrylonitrile Butadine Styrene)

1.6.2 Dimension: 127 X 130 X 160 mm

1.6.3 It contains special designed silicon feet base excellent for damping vibration.

1.6.4 Swirlex is eccentric with oil less ball bearings.

1.6.5 Built from sturdy cast casing.

1.6.6 CE certified

1.6.7 Chemical resistance: see website (www.abdoslifesciences.com)

1.6.8 Intended use: Widely used in scientific studies of molecular biology, clinical chemistry, biochemistry and any type of scientific laboratory.

1.7 Spare Accessories

E11191 Tube holding rod.
E11192 Universal 100 mm plate with 4 vacuum suction feet
E11198 Platform pad for tubes and vessel Dia 99 mm
E11197 Tube adaptor for 48 holes test tubes Dia 6 mm

TECHNICAL DATA SHEET

Doc No – Ab/TDS/Eqpt-233
Rev No: 09 Date- 05/2024-25

Abdos Labtech Pvt. Ltd. SWIRLEX – VOTREX SHAKER

E11193 Tube adaptor for 15 holes test tubes Dia 10 mm
E11194 Tube adopter for 12 holes test tubes Dia 12 mm
E11195 Tube adopter for 8 holes test tube Dia 16 mm
E11196 Tube adaptor for holes test tube Dia 20 mm

*** E11192 is required to run Abdos E11198, E11197, E11193, E11194, E11195 & E11196 adapters.**

Warranty: 12 months from the date of Installation or 15 months from the date of Invoicing

Instruction for Use

Storage temp. Max 30°C
Before using the product check crackness, etc.
Check all the information printed on pouches or boxes.

WARNING

DO NOT store strong oxidizing agents.
DO NOT place any plastic labware in a flame.
DO NOT mix any chemicals that may result in a thermic reaction, which can cause product failure
Consult local fire codes prior to storage of flammable liquids.

Manufactured & Marketed By :
Abdos Labtech Private Limited
Plot No. 1 Lakshewari Bhagwanpur Roorkee 247667 (UK) India
For Complaints/Feedback Contact Customer Care at above address or
Ph.: +91-011-41081897 or Email us at labtech@abdosindia.com
Lic. No.: MFG/IVD/2024/000115 MADE IN INDIA



Scan to visit website

TECHNICAL DATA SHEET

Doc No – Ab/TDS/Eqpt-233
Rev No: 09 Date- 05/2024-25

Abdos Labtech Pvt. Ltd.
SWIRLEX – VOTREX SHAKER

Quality Manager
Abdos Labtech Private Limited
Employee Code 30110

